



# 2020 Liberty Mutual Insurance Industrial Hygiene Lab sampling guide

AIHA-LAP, LLC Laboratory ID # LAP-100045

**Contact us for assistance:**

[lmihlaboratory@libertymutual.com](mailto:lmihlaboratory@libertymutual.com)

or

1-800-230-6263



Liberty Mutual Insurance Industrial Hygiene Laboratory (LMIHL) provides analytical services to assist customers with safety issues related to occupational disease exposure in their workplace.

**Comprehensive service includes:**

- Method development or validation
- Technical consultation
- Training
- Sample Analyses:
  - Aldehydes
  - Organic and Inorganic acids
  - Amines
  - Pharmaceuticals
  - Asbestos fiber counting and bulk identification
  - Silica
  - Medical gas analysis
  - Common and specialty metals
  - Grade D breathing air analysis
  - Common and specialty organic solvents
  - Scans – metals, solvent, acids, aldehydes, isocyanates, etc.
  - Gravimetrics for total, respirable or inhalable dust

Liberty Mutual SafetyNet™ is your source for chain of custody and media order forms as well as the Fee Schedule. You may also contact the lab.

**Telephone assistance:** 1-800-230-6263

Telephone assistance is available to all our customers from 7:30AM – 6PM, Eastern Time. We encourage our customers to call us with their technical questions relating to proper sampling and media selection.

Listed below is our support team.

- Ethel T. Patricio, MS, Laboratory Director, ext. 27352 or 508-544-5352
- Eva Longo, MS, Prod/QC Chemist Consultant, ext 27338 or 508-544-5338
- Laura Melton, MS, Prod/QC Chemist Consultant, ext 27348 or 508-544-5348
- Srma Gulavita, MS, Prod/QC Chemist Consultant, ext. 27380 or 508-544-5380

**Analysis turnaround time**

The Laboratory's standard turnaround time is 5 business days or better from the receipt of samples. Samples received after 3PM are considered next day's samples. Results are provided by email only. Rush analysis for routine samples can be arranged with prior notification. Verbal results are also available after data has been approved for release by the QC coordinator.

**The following surcharges apply for RUSH requests:**

<b>One business day</b>	200% surcharge
<b>Two business days</b>	150% surcharge
<b>Three business days</b>	75% surcharge

**Note:** Certain Specialty Analysis cannot be done on a rush basis. Contact the laboratory at least 24 hours before shipping your samples for RUSH analysis.

**Discounts for multiple analytes on the same media**

The price listed for each individual analyte is the price for a single analyte per sample. Analysis for multiple compatible analytes collected on the same media is available at a discounted fee. The higher priced analyte is billed at the published fee with additional analytes discounted as follows:

Organics on OVM/sorbent tubes	\$32
Metals by ICP	\$26
IC	\$44
HPLC/Specialty	see Fee Schedule

## Payment terms

The total cost of the services provided by the LMIH Laboratory will be based on the quoted rate. Prices are subject to change without notice. Payment term is net 30 days if you are paying by check.

## Blanks

Blanks are required by all analytical methods and are good insurance for dealing with contamination. Contamination can occur during handling, storing or shipping samples. Therefore, as standard good IH practice, field blanks are recommended. The recommended number of field blanks per sample set is 10% of the total number of samples, or a minimum of 1 blank per set of 10 samples. Since blanks are analyzed as any other sample, they are priced according to the analysis requested.

## Media

Common media such as charcoal tubes, pre-weighed filters, etc. are provided at no extra cost as long as the samples and unused media are returned to the laboratory for analysis within 30 days. Unreturned media after 30 days will be charged at cost. Specialty media that cannot be reused, such as isocyanate filters, aldehyde badges and tubes, ozone filters, etc. will be invoiced at cost. Media are shipped UPS ground at no charge to the client. Rush and international shipments will be billed to the client at cost.

- Request by email using the media order form.
- Return shipping is the responsibility of the customer.

**Please note that there will be media charge for specialty filters (e.g. hydrogen peroxide) and sampling devices such as PPI.**

## Sample collection and sample submission

### Sampling collection supplies

We provide free collection supplies (wipes, templates, gloves, plastic bags, coolers, etc.) and ship them to you via UPS ground. Clients requesting Rush delivery either within or outside the U.S. will be billed at cost. It is the client's responsibility to ensure proper sampling, handling, packing (returning the cooler with frozen packs) and return shipping to the laboratory.

Chain-of-custody forms will be provided with every media and equipment order; please complete the form and submit with your samples. The chain-of-custody forms are also available on Liberty Mutual SafetyNet™.

### Sample minimum

When Preparing samples for analysis, please note that some analyses require a three (3) sample minimum fee to cover the cost of method set-up for non-routine analysis.

## Special handling and shipping

Most NIOSH and OSHA methods indicate that solvent samples are stable at room temperature. LM IH Laboratory recommends shipping samples cold overnight especially during the summer months. Please refer to the list of Analytes for additional instructions. These guidelines must be followed in order to ensure the integrity and validity of the samples.

## Bulk sample submissions

Certain analysis requires the submission of bulks. Bulks should always be packaged and shipped separately from the samples. Safety Data Sheets (SDS) must accompany the bulk samples.

## Analytical results

LMIHL reserves the right to determine appropriate format in which the analytical results are reported. All results are provided for the exclusive use of our client. LMIHL accepts no responsibility or liability for the client's use of the analytical results.

LMIHL may release verbal or email results in ahead of the written report. These results are tentative and subject to subsequent confirmation or modification during LMIHL peer review process.

LMIHL requires precise and complete instruction before it releases any reports. Any later request will require written permission from submitter.

## Free pump loan program

We loan pre-calibrated sampling pumps and accessories to customers for a period of two (2) weeks at no cost if the samples are returned to the lab for analysis. All equipment is loaned on a first-come, first-serve basis. As a result, there may be instances when the equipment is not available at the time of your request. We suggest submitting your request for equipment at least one week prior to your scheduled sampling activity.

At the end of the two week period, all equipment must be returned and collected samples sent for analysis to the LMIH Lab. Equipment not returned after the two week period will incur an equipment rental charge.

Contact your local Loss Control Representative or the LMIH lab for information/questions ([lmihlaboratory@libertymutual.com](mailto:lmihlaboratory@libertymutual.com) or 800-230-6263 ext. 27371) regarding your sampling needs.

## Scans

The scans performed by the laboratory are listed below. We can also customize scans.

Aldehyde scan (Glutaraldehyde is collected and analyzed separately)			
Acetaldehyde	Benzaldehyde	Formaldehyde	Valeraldehyde
Acrolein	n-Butyraldehyde	Propionaldehyde	
Aliphatic amine scan (NIOSH 2010)			
(Please call Lab for other aliphatic amines)			
Ethylamine	Diethylamine	Triethylamine	
Aromatic amine scan			
(Please call Lab for other aromatic amines)			
Aniline	Methyl aniline	o-Toluidine	
Anesthetic gases scan (Nitrous oxide is collected and analyzed separately)			
Desflurane (Suprane)	Halothane (Fluothane)	Sevoflurane (Sevofrane)	
Enflurane (Ethrane)	Isoflurane (Forane)		
Inorganic acid scan			
Bromide	Fluoride	Phosphate	
Chloride	Nitrate	Sulfate	
Isocyanate scan			
Hexamethylene diisocyanate (HDI)		Isophorone diisocyanate (IPDI)	
4,4-Methylene bisphenyl isocyanate (MDI)		2,4-Toluene diisocyanate (2,4-TDI)	
2,6-Toluene diisocyanate (2,6-TDI)			
Metals – Liberty's 20-Metal scan			
Aluminum	Cadmium	Lead	Tin
Antimony	Chromium	Magnesium	Titanium
Arsenic	Cobalt	Manganese	Thallium
Beryllium	Copper	Nickel	Vanadium pentoxide as V
Calcium	Iron	Selenium	Zinc

Metals – 13-Metal scan (welding fume scan)			
Antimony	Cobalt	Manganese	Vanadium pentoxide as V
Beryllium	Copper	Molybdenum	Zinc
Cadmium	Iron	Nickel	
Chromium	Lead		

Organic acid scan			
Acetic acid	Butyric acid	Formic acid	Propionic acid

Organic solvents (GC/MS) scan		
(Qualitative analysis. For quantitative analysis please call lab.)		
Acetone	Dimethyl formamide (DMF)	Methyl amyl ketone (MAK)
Acetonitrile	1,4-Dioxane	Methyl cyclopentane
Acrylonitrile	Epichlorohydrin	Methyl ethyl ketone (MEK)
n-Amyl acetate	Ethanol	Methyl isoamyl ketone (MIAK)
Benzene	Ethyl acetate	Methyl isobutyl ketone (MIBK)
Sec-Butanol	Propylene glycol methyl ether acetate (PGMEA)	Methyl propyl ketone (MPK)
n-Butyl acetate		Methyl t-butyl ether (MTBE)
s-Butyl acetate	Propylene glycol monomethyl ether (PGME)	a-Methyl styrene
n-Butyl acrylate		Methylene chloride
n-Butyl alcohol	Styrene	Methyl methacrylate
Butyl cellosolve (2-Butoxyethanol)	Tetrahydrofuran	Nitromethane
Butyl cellosolve acetate	Toluene	3-Pentanone (Diethyl ketone)
Carbon tetrachloride	1,1,1-Trichloroethane	Perchloroethylene
Cellosolve (2-Ethoxyethanol)	Propylene glycol butyl ether (PGBE)	a-Pinene
Cellosolve acetate	Ethyl lactate	n-Propanol
Chlorobenzene	Ethyl benzene	2-Propoxyethanol
Chloroform	1-Ethoxy-2-propyl acetate	n-Propyl acetate
o-Chlorotoluene	Heptane	Ethyl acrylate
Cumene	Hexane	1,1,2-Trichloroethane
Cyclohexanone	Isobutanol	Trichloroethylene
Diacetone alcohol	Isobutyl acetate	1,2,4-Trimethylbenzene
1,1-Dichloroethane	Isopropanol	1,3,5-Trimethylbenzene
1,2-Dichloroethane	Isopropyl acetate	Vinyl acetate
1,2-Dichloroethylene	d-Limonene	Xylene
Diisobutyl ketone	Methyl acetate	

PNAs (NIOSH 5506)			
Acenaphthene	Benzo[b]fluoranthene	Chrysene	Indeno[1,2,3-cd]pyrene
Acenaphthylene	Benzo[k]fluoranthene	Dibenz[a,h]anthracene	Naphthalene
Anthracene	Benzo[ghi]perylene	Fluoranthene	Phenanthrene
Benz[a]anthracene	Benzo[a]pyrene	Fluorene	Pyrene

PNAs (OSHA 58)		
Anthracene	Chrysene	Pyrene
Benzo[a]pyrene	Phenanthrene	

## Sampling guide analyte descriptions and abbreviations

The following information is included in the Sample Guide's alphabetical listing of analytes:

- Analyte** Analytes are listed by their common name in alphabetical order in this sampling guide. The synonyms are listed within parenthesis. TVOC stands for Total Volatile Organic Compounds. VM&P Naphtha stands for Varnish Makers and Painters Naphtha.
- CAS#** Chemical Abstract Service number for the compound to be sampled.

## Analytical method

This specifies the preferred analytical method used by Liberty Mutual for the analysis of the compound.

- ASTM** American Society for Testing and Materials
- Ciba-Geigy** Ciba-Geigy In-House Method for 1,3,5-triglycidyl isocyanurate
- DuPont** DuPont In-House Method for Perfluorooctanoic Acid
- HSE MDHS** UK Health & Safety Executive Methods for the Determination of Hazardous Substances
- LMI** Liberty Mutual Insurance "in-house" analytical methods
- NIOSH** NIOSH Manual of Analytical Methods
- OSHA** OSHA Manual of Analytical Methods
- YAMATE LEVEL II** Airborne Asbestos by Transmission Electron Microscopy

## Analytical technique

- AA** Atomic absorption spectrophotometry
- AA-CV** Cold vapor atomic absorption spectrophotometry
- EGA-TDA** Evolved gas analysis-thermo dilatometric analyzer
- FTIR** Fourier Transform Infrared spectrophotometry
- GC-DID** Gas chromatography with pulsed discharge ionization detector
- GC-ECD** Gas chromatography with electron capture detector
- GC-FID** Gas chromatography with flame ionization detector
- GC-MS** Gas chromatography-mass spectrometry
- GC-NPD** Gas chromatography with nitrogen phosphorus detector
- GC-TCD** Gas chromatography with thermal conductivity detector
- GC-XSD** Gas chromatography with halogen specific detector
- GFAA** Graphite furnace atomic absorption spectrophotometry
- GRAV** Gravimetric analysis
- HPLC** High performance liquid chromatography
- IC** Ion chromatography
- ICP** Inductively coupled plasma spectroscopy
- ICP-MS** Inductively coupled plasma spectroscopy – mass spectrometry

<b>ISE</b>	Ion selective electrode
<b>LC-MS</b>	High performance liquid chromatography – mass spectrometry
<b>PCM</b>	Phase contrast microscopy
<b>PLM</b>	Polarized light microscopy
<b>TEM</b>	Transmission Electron Microscopy; Energy dispersive X-ray (EDX) analyzer
<b>UV/VIS</b>	UV/VIS spectrophotometry
<b>XRD</b>	X-ray diffractometry

## Sampling media

The recommended sampling media for each of the methods in this sampling guide are:

<b>AgMF</b>	25mm 0.45um Silver membrane filter (SKC 225-1802)
<b>Anasorb708</b>	SKC Anasorb 708 sorbent tube (SKC 226-30-08)
<b>Anasorb747</b>	SKC Anasorb 747 sorbent tube (SKC 226-81A)
<b>Anasorb747/Anasorb747</b>	SKC Anasorb 747 sorbent tube in series (SKC 226-82)
<b>Anasorb747, Treated</b>	SKC Anasorb 747 treated with tert-Butyl catechol (SKC 575-006)
<b>AT Monitor</b>	Assay Technology monitor for aldehydes (N517AT)
<b>AT N2O Monitor</b>	Assay Technology monitor for nitrous oxide (X575AT)
<b>Bulk</b>	Bulk sample tube
<b>Cellulose Nitrate,Na2CO3</b>	Cellulose nitrate filter treated with sodium carbonate (SKC 225-9031)
<b>Carulite</b>	Sorbent tube for mercury, replaces Hopcalite (SKC 226-17-1A/3A)
<b>CS106</b>	Chromosorb 106 sorbent tube (SKC 226-111A)
<b>CT</b>	Charcoal tube (SKC 226-01, 226-09)
<b>CT-CT</b>	2 Charcoal tubes in series
<b>CT, KOH</b>	Potassium hydroxide treated Anasorb CSC coconut charcoal tube (SKC 226-67)
<b>Cylinder</b>	300cc Aluminum cylinder
<b>di H2O</b>	Deionized water
<b>GFF</b>	Glass fiber filter
<b>GFF, 1-2PP</b>	1-(2-Pyridyl)piperazine treated glass fiber filter
<b>GFF, Acid</b>	Sulfuric acid treated glass fiber filter
<b>GFF, HBr</b>	Hydrogen bromide treated glass fiber filter
<b>GFF-Florisil</b>	Millipore Swinnex 13 with glass fiber filter (SX0001300/-01/ AP2001300) -Florisil tube (SKC 226-39) in series
<b>GFF, IOM</b>	IOM sampler with glass fiber filter
<b>GFF-SGT</b>	Millipore Swinnex 13 with glass fiber filter (SX0001300/-01 / AP2001300) -silica gel tube (SKC 226-10) in series
<b>GFF, NaNO2</b>	Sodium nitrite treated glass fiber filter
<b>GFF,Vamine</b>	Glass fiber filter coated with 10 mg of veratrylamine
<b>GFF- PE,HgCl2</b>	Glass fiber filter -polyester filter treated with mercuric chloride in series (SKC 225-9018)
<b>Ghost wipe</b>	SKC wipe for surface lead or other metals (SKC 225-2414)
<b>IABC</b>	Impregnated activated beaded carbon (SKC 226-80)



<b>Impinger 4</b>	Midget fritted glass bubbler containing 0.02% potassium iodide in sodium carbonate/sodium bicarbonate buffer
<b>IOM</b>	Personal inhalable sampler developed at the Institute of Occupational Medicine (IOM) in Scotland
<b>MCE</b>	Mixed cellulose ester filter membrane ((Zefon 728 MCE))
<b>MCE2</b>	Mixed cellulose ester filter membrane, 0.8µm, 25mm (Zefon 528 MCE) Mixed cellulose nitrate (SKC 225-9032)
<b>MCE, 25mm 0.8um</b>	Zefon PCM air sampling cassettes with conductive cowl (ZEFON Z008BA)
<b>MCE, 25mm 0.45um</b>	Zefon TEM air sampling cassettes with conductive cowl (ZEFON Z045BA)
<b>MCE,carbonate</b>	Mixed cellulose ester filter treated with 20:1 sodium carbonate: glycerol solution
<b>MCE-CS102</b>	Mixed cellulose ester filter-Chromosorb 102 tube (SKC 226-104)
<b>MCE, SGT**</b>	Mixed cellulose ester filter — SKC specially cleaned silica gel tube (226-10-03)
<b>ORBO 34</b>	Supelco specially treated charcoal tube for H2S (20211)
<b>ORBO 77</b>	Supelco carbon beads treated with sulfuric acid for ammonia(20036)
<b>ORBO 78</b>	Supelco Carboxen-564 carbon molecular sieve treated with hydrogen bromide (20355)
<b>ORBO 92</b>	Supelco Carboxen-564 carbon molecular sieve (20362)
<b>OVM</b>	3M 3500 organic vapor monitor
<b>OVM 3520</b>	3M 3520 organic vapor monitor with backup section
<b>OVM 3551</b>	3M 3551 organic vapor monitor for ethylene oxide
<b>OVS-2/GFF</b>	SKC XAD-2/glass fiber filter sorbent tube (226-30-16 OVS)
<b>OVS-2/QF</b>	SKC XAD-2/quartz filter sorbent tube (226-58 OVS)
<b>OVS-7</b>	SKC XAD-7/glass fiber filter sorbent tube (226-57 OVS)
<b>OVS-Tenax</b>	SKC tenax/glass fiber filter sorbent tube (226-56 OVS)
<b>Oxidizer</b>	SKC oxidizer with TEA-IMS for nitric oxide sampling
<b>PS</b>	SKC 520 inorganic mercury passive sampler (520-02A/03)
<b>Paint chips</b>	Paint chips sample in polyethylene bag
<b>Porapak-P</b>	SKC Porapak-P tube (226-114)
<b>Porapak-Q</b>	SKC Porapak-Q tube (226-115)
<b>PTFE</b>	25mm 0.5um Polytetrafluoroethylene (PTFE) filter (SKC 225-1708)
<b>PTFE1</b>	37mm 1.0um Polytetrafluoroethylene (PTFE) filter (Zefon FPTFE137)
<b>PTFE3</b>	37mm 5.0um Polytetrafluoroethylene (PTFE) filter (SKC 225-17A)
<b>PTFE4</b>	PALL Life Sciences 47mm, 0.45um polytetrafluoroethylene (Teflon) filter (TF-450, P/N 66149)
<b>PTFE5</b>	25mm 1.0um Polytetrafluoroethylene (PTFE) filter (SKC 225-2714)
<b>PTFE-SGT</b>	Teflon filter-silica gel tube in series
<b>PTFE/XAD2</b>	Teflon filter-XAD-2 tube in series
<b>PTFE-AgMF</b>	SKC polytetrafluoroethylene (Teflon) filter-silver membrane filter in series (225-1708) (225-1802)
<b>PVC</b>	37mm 5.0um pre-weighed polyvinyl chloride filter (SKC 225-5-37-P)
<b>PVC-KOH</b>	Polyvinyl chloride filter- Midget fritted glass bubbler containing 0.1 N potassium hydroxide in series
<b>QFF</b>	Quartz fiber filter (Millipore AQFA03700)



<b>QFF, titanium oxysulfate</b>	Quartz fiber filter (SKC 225-9030)
<b>QFF, Na<sub>2</sub>CO<sub>3</sub></b>	Quartz fiber filter (Millipore AQFA03700) treated with sodium carbonate (SKC 225-9032)
<b>Sep-Pak</b>	Waters dinitrophenylhydrazine (DNPH) treated cartridge for aldehydes (WAT047205)
<b>SGT</b>	Silica gel tube (SKC 226-10) with sorbent 75/150 mg
<b>SGT</b>	Silica gel tube (SKC 226-51) with sorbent 50/100 mg
<b>SGT**</b>	SKC specially cleaned silica gel tube (226-10-03)
<b>SGT,Acid</b>	SKC sulfuric acid treated silica gel tube (226-10-06)
<b>SGT, DNPH</b>	SKC dinitrophenylhydrazine impregnated silica gel tube (226-119)
<b>SGT/GFF-SGT/GFF</b>	2 specially washed and dried silica gel tubes in series (SKC 226-183)
<b>SGT, HgCl<sub>2</sub></b>	SKC mercuric chloride coated silica gel tube (226-10-02)
<b>Soda Lime</b>	SKC soda lime tube (226-210)
<b>TEA-IMS</b>	Triethanolamine-impregnated molecular sieve tube (SKC 226-40-02 / 226-40A with oxidizer)
<b>Tenax</b>	SKC tenax tube (226-35/226-35-01)
<b>Whatman 42</b>	Whatman 1442-070
<b>XAD-2</b>	SKC sorbent tube (226-30/226-30-04/226-30-06)
<b>XAD-2,NITC</b>	SKC sorbent tube coated with 10% 1-naphthylisothiocyanate (NITC) (226-30-18)
<b>XAD-2,p-An</b>	SKC sorbent tube (50/100) coated with 0.5mg of p-Anisidine (226-30-07)
<b>XAD-7</b>	SKC sorbent tube (226-95)
<b>XAD-7, Acid</b>	SKC XAD-7 sorbent (1,2-Dichloroethane) t tube treated with 10% phosphoric acid (226-98)

## Sampling rate (flow rate)

For passive monitors (OVM and AT monitor): cc/min. For all other compounds listed in the guide: liters per minute (lpm).

The sampling time in minutes for passive monitors and the sampling volume range indicates the minimum and maximum volume in liters. For bulk samples, the amount needed is listed in grams in this column. Minimum volumes are typically calculated to allow quantification at 10% of the occupational exposure limits unless otherwise stated.

## LOQ (Limit of Quantitation)/LOD (Limit of Detection)

Limit of Quantification and Limit of Detection are reported to two significant digits: %-Percent for bulk sample, fib/fld- Fibers per field, and µg- Micrograms/sample.

## Compatibility code

This code indicates analytes that can be collected and analyzed simultaneously on a single sample. The absence of a code means that the analyte is not compatible with other analytes and would need to be collected on separate samples. The codes are as follows:

<b>1%DMF/CS<sub>2</sub></b>	Desorption in 1% dimethylformamide in carbon disulfide
<b>1%IPA/CS<sub>2</sub></b>	Desorption in 1% isopropanol in carbon disulfide
<b>1%PRO/CS<sub>2</sub></b>	Desorption in 1% n-propanol in carbon disulfide
<b>5%IPA</b>	Desorption in 5% isopropanol in deionized water
<b>5% IPA/ CS<sub>2</sub></b>	Desorption in 5% isopropanol in carbon disulfide
<b>5%PRO/CS<sub>2</sub></b>	Desorption in 5% n-propanol in carbon disulfide
<b>95%EtOH</b>	Desorption in 95% ethanol in deionized water

<b>AC/CS<sub>2</sub></b>	Desorption in 2% acetone in carbon disulfide
<b>AC/MeOH</b>	Desorption in 1% methanol in acetone
<b>Acid1</b>	Inorganic acids group 1 (see List of Scans for Inorganic Acids)
<b>Acid2</b>	Organic acids group 2 (formic acid, acetic acid, butyric acid, and propionic acid)
<b>ACN</b>	Desorption in acetonitrile
<b>ACN/TOL</b>	Desorption in 50% acetonitrile/toluene
<b>ACN/DMSO</b>	Desorption in 90% acetonitrile/dimethylsulfoxide
<b>Acetone</b>	Desorption in acetone
<b>Aldehyde</b>	Aldehyde group
<b>Amine1</b>	Aliphatic amine group by GC-FID
<b>Amine2</b>	Aliphatic amine group by HPLC
<b>Amine3</b>	Aromatic amine group
<b>Benzene</b>	Desorption in benzene
<b>BUT/CS<sub>2</sub></b>	Desorption in 1% 2-butanol in carbon disulfide
<b>CCl<sub>4</sub></b>	Desorption in carbon tetrachloride
<b>Cl<sub>2</sub>&amp;Br<sub>2</sub></b>	Chlorine and bromine
<b>CS<sub>2</sub></b>	Desorption in carbon disulfide
<b>dil acid</b>	Diluted sulfuric acid
<b>DMF/CS<sub>2</sub></b>	Desorption in 50% dimethylformamide in carbon disulfide
<b>EA</b>	Ethanolamine/diethanolamine/triethanolamine
<b>Ethyl Acetate</b>	Desorption in ethyl acetate
<b>Ethyl Ether</b>	Desorption in ethyl ether
<b>FA</b>	Desorption in formic acid
<b>Isocyanate</b>	Isocyanate and diisocyanate group
<b>MC</b>	Desorption in methylene chloride
<b>MeOH</b>	Desorption in methanol
<b>MeOH/CS<sub>2</sub></b>	Desorption in 1% methanol in carbon disulfide
<b>%MeOH/MC</b>	Desorption in mixture of methanol and methylene chloride
<b>Metals</b>	Metal group by either NIOSH 7301 or OSHA ID-125G
<b>Metal1</b>	Compatible metal group by modified NIOSH 7300 method
<b>Metal2</b>	Compatible metal group by modified OSHA ID-121 method
<b>NO&amp;NO<sub>2</sub></b>	Nitric oxide and nitrogen dioxide
<b>Phenol&amp;cresol</b>	Phenol and cresol group
<b>PNAs</b>	Polynuclear aromatic hydrocarbons
<b>Silica</b>	Silica by XRD
<b>Tol</b>	Desorption in toluene

## Interferences

This column lists the possible interferences as stated in the analytical method.

## Comments

This column contains special instructions for sample collection and handling.

## Notes for gas analysis in comments

### Note 1, Breathing air: Grade D and Grade E, Collected in Cylinders

Use when sampling a compressed gas line at 50 psig (i.e., before it goes into the regulator of the SCBA). If the pressure is lower than 50 psig, be sure to record the line pressure. This method reports CO, CO<sub>2</sub>, %O<sub>2</sub>, total hydrocarbons as CH<sub>4</sub> (TH), total halogenated hydrocarbons as CCl<sub>4</sub> (HH) and dew point (DP) in ppm and degrees F at 0 psig. See Procedure for Removing the Cylinder and LOQ at end of notes. Gravimetric analysis for condensed oil (total particulate) requires a 47mm Teflon filter.

### Note 2, Medical gases – Contamination (NFPA 99)

Use when sampling a compressed medical gas line at 50 psig. When testing medical air, this method reports CO, CO<sub>2</sub>, %O<sub>2</sub>, total hydrocarbons as CH<sub>4</sub> (TH), total halogenated hydrocarbons as CCl<sub>4</sub> (HH) and dew point (DP) in ppm and degrees C at 50 psig. When testing nitrogen, this method reports CO, CO<sub>2</sub>, total hydrocarbons as CH<sub>4</sub> (TH), total halogenated hydrocarbons as CCl<sub>4</sub> (HH) and dew point (DP) in ppm and degrees C at 50 psig. For oxygen, this method reports the same as medical air, except dew point. For nitrous oxide, this method reports CO, % Air and total halogenated hydrocarbons as CCl<sub>4</sub>. See Procedure for Removing the Cylinder and LOQ at end of notes. Gravimetric analysis for condensed oil (total particulate) requires a 47-mm Teflon filter. The minimum sampling time for the gravimetric sample is 10 minutes in order to sample at least 1000 liters of gas, at 50 psig.

## Procedure for removing the cylinder

A common mistake is removing the cylinder and the “A” fixture together. The correct steps are: Remove fixture “B.” The pop-off valve will activate. Remove the cylinder from fixture “A”, then remove fixture “A” from the gas line. Please refer to sampling instructions. Failure to follow the sampling instructions will result in the cylinder not being pressurized. ***Samples that are not pressurized cannot be analyzed.***

## Limits of quantification

	Grav.	Air	CO	CO <sub>2</sub>	Dew Pt	O <sub>2</sub>	TH	HH
<b>GAS</b>	µg/s	%	ppm	ppm	see notes	%	ppm	ppm
<b>Grade D&amp;E</b>	50	NA	0.50	25	200ppm, -32.6°F	0.30	0.50	0.50
<b>Medical Air</b>	50	NA	0.50	25	200ppm, -21.5°C	0.30	0.50	0.50
<b>Nitrogen</b>	50	NA	0.50	25	200ppm, -21.5°C	NA	0.50	0.50
<b>Oxygen</b>	50	NA	0.50	25	NA	0.30	0.50	0.50
<b>Nitrous Oxide</b>	50	0.15	0.50	NA	NA	NA	NA	0.50

## Fee schedule for common analyses

This fee schedule gives a listing of prices for the most frequently requested analyses, effective April 1, 2019 through March 31, 2021. If the analysis you require is not listed, please contact us for a price estimate.

### Metals

	First	Additional
Common metals by ICP <sup>1,4</sup>	\$36	\$26
Beryllium oxide	\$55	
Chromium VI (paint) <sup>1</sup>	\$100	
Chromium VI (soluble) <sup>1</sup>	\$75	
13 metal scan (welding fume scan) <sup>1,2,4</sup>	\$150	
20 metal scan <sup>1,2,4</sup>	\$190	

### Microscopy

Asbestos fiber count (PCM)	\$35	
Asbestos fiber identification (PLM)	\$40	

### Gravimetric analysis

Dust (carbon black (NIOSH 5000), oil mist, respirable, welding fume)	\$27	
Inhalable dust <sup>3</sup>	\$35	
Condensed oil/particulate	\$50	

### High performance liquid chromatography (HPLC)

Common HPLC (aldehydes & isocyanates)	\$95	\$45
Anhydrides (maleic & trimellitic) <sup>3</sup>	\$100	\$50
Amines (if collected on NITC tubes) <sup>3</sup>	\$85	\$45
Aldehyde scan	\$230	
Isocyanate scan	\$190	
Pharmaceuticals <sup>3</sup>	\$110	\$50

### Ion chromatography (IC)

Common IC (inorganic acids, anion & cation) <sup>4</sup>	\$65	\$44
Chlorine dioxide <sup>3</sup>	\$70	
Chlorine & bromine <sup>3</sup>	\$75	
Fluorides (gaseous and particulate) <sup>3,4</sup>	\$115	
Organic acids (such as formic, acetic, butyric, lactic, propionic)	\$70	\$44
Inorganic and organic acid scan	\$165	

### Silica (by XRD)

	Air	Bulk
Quartz	\$70	\$100
Quartz and cristobalite	\$85	\$130
Quartz, Cristobalite and tridymite	\$105	\$140

## Solvents

### Common solvents

	First	Additional
Solid sorbent tubes (CT, SGT, etc.)	\$55	\$32
Organic vapor monitors (3M 3500/ 3520)	\$72	\$32
GC/MS Scan (15 analytes quantified)	\$250	

### Other solvents

Amines <sup>3</sup>	\$85	\$50
Diacetyl or acetoin	\$95	\$45
Ethylene oxide	\$105	
Glycols and phthalates	\$75	\$40
Organo sulfur compounds <sup>3</sup>	\$105	\$55
Mixtures on tubes or sorbent tubes in series	\$70	\$32
PCBs <sup>3</sup>	\$90	\$45
1,3,5 Triglycidyl isocyanurate <sup>3</sup>	\$95	

## Special methods analyses

Aminoethanols	\$88	\$45
Asphalt fume (NIOSH 5042)	\$75	
Carbon black (OSHA 196)	\$80	
Coal tar pitch volatile, BSF	\$90	
Cyanide (NIOSH 7904)	\$75	
Gas analysis (medical and Grade D Air)	\$165	
Hydrogen peroxide (OSHA VI-6)/filter <sup>4</sup>	\$85	
Metal working fluids (NIOSH 5524)	\$80	
Oil mist mineral by FTIR (NIOSH 5026) <sup>3</sup>	\$75	
PNA scan (NIOSH 5506)	\$280	
PNA scan (OSHA 58)	\$190	
Specialty metals by ICP	\$55	\$32
Specialty metals by AA (includes mercury)	\$65	\$32
Metals by ICP/MS <sup>3</sup>	\$50	\$32
Organic tin <sup>3</sup>	\$80	

**Prices are subject to change without notice.**

## Special charges

1. Preparation charge applies to bulk and wipes: \$10 per sample.
2. Customized reporting: \$40 per report.
3. Three sample minimum required.
4. Media charge applies to some specialty media such as PPI and treated filters.

### Method development/validation

We partner with our customers to develop new methods for analytical testing, as well as validate existing methods. Please contact the lab for quote on method development or validation costs.

### Sample blanks

Sample blanks are recommended for all sampling activities and are charged at the same rate as the sample.

## Sampling guide

Sampling Guide provides further information about air sampling and our lab services. The Sampling Guide may be accessed by contacting the lab (lmihlaboratory@libertymutual.com) or for our policy holders via Liberty Mutual SafetyNet™.

## Terms and conditions

The following statements describes the terms and conditions under which the Liberty Mutual Industrial Hygiene Laboratory (LMIHL) operates.

The client is responsible for the condition and custody of all samples prior to receipt, inspection and acceptance by LMIHL.

Methods used by LMIHL to analyze your samples are compliant with NIOSH, OSHA and other regulatory agencies. LMIHL reserves the right to interpret these methodologies when applying them to the analysis of client's samples based on reasonable, professional judgment of LMIHL personnel and recognized industry standards.

LMIHL reserves the right to use our standard template in reporting analytical results. Where reasonable, we will make every effort to honor requests for special hardcopy or electronic formats if requested in advance. LMIHL requires authorization in writing when requesting additional distribution of lab reports to other than the client

LMIHL may release reports upon the request of the client either verbally or by email. Such reports are considered tentative and may be subject to modification after completion of QA/QC review.

Report will only contain analytical data. LMIHL is not in the position to interpret data as they pertain to regulations, calculation of time-weighted average exposures from analytical results, etc.

## Equipment rental

Inhalable Samplers: \$20 each per week (two week maximum); Calibrator adaptor: \$10 per week additional.

## Sampling devices and media

Sampling pumps and compressed gas analysis media

Pre-calibrated air sampling pumps and gas analysis (Grade D, E or Medical Gas) air sampling equipment are loaned at no charge for two weeks. At least one week notice prior to your "need by date" is required by the lab. After two weeks, charges of \$15 per day, per pump, or per gas cylinder or filter will apply.

## Return policy

Media is provided at no cost when returned to laboratory for analysis within 30 days. After 30 days, unreturned media will be invoiced at cost plus shipping. Returned unused specialty media that cannot be reused (isocyanates filters, aldehyde badges, sorbent tubes, ozone filters, etc.) will be charged at cost plus shipping.

## Shipping

UPS ground shipping within the 48 states is provided for sampling media at no extra charge. Return shipping is the client's responsibility. Rush and international shipping charges will be added to analysis fees. Media orders for same day shipping must be submitted by 2 p.m. Eastern Time, Monday through Friday. Media requiring cold shipping cannot be delivered over the weekend. Order the media by contacting the lab.

## Common analytes

**Common aldehydes:** Acetaldehyde, acrolein, formaldehyde, benzaldehyde, butyraldehyde, valeraldehyde, and propionaldehyde.

**Common Isocyanates:** Hexamethylene diisocyanate (HDI), Isophorone diisocyanate (IPDI), Methylene bis-phenyl isocyanate (MDI), Toluene-2,4-diisocyanate (2,4-TDI), Toluene-2,6-diisocyanate (2,6-TDI).

**Common metals:** Aluminum, Antimony, Arsenic, Barium, Beryllium, Calcium, Cadmium, Chromium, Cobalt, Copper, Iron, Lanthanum, Lead, Lithium, Magnesium, Manganese, Molybdenum, Nickel, Phosphorus, Platinum, Potassium, Selenium, Sodium, Silver, Strontium, Tellurium, Tin, Titanium, Thallium, Vanadium, Zinc, Zirconium

**Specialty metals:** Arsenic Trioxide, Arsine, Beryllium Oxide, Bismuth, Boron, Gallium, Germanium, Germane, Gold, Indium, Palladium, Phosphine, Rhodium, Silane, Titanium dioxide, Tungsten, Yttrium

**Welding fume:** Includes antimony, beryllium, cadmium, chromium, cobalt, copper, iron, lead, manganese, molybdenum, nickel, vanadium, zinc.

**20 metal scan:** Aluminum, Antimony, Arsenic, Beryllium, Calcium, Cadmium, Chromium, Cobalt, Copper, Iron, Lead, Magnesium, Manganese, Nickel, Selenium, Tin, Titanium, Thallium, Vanadium, Zinc. Can be customized. Please contact lab.

**Inorganic acid scan:** hydrogen bromide, hydrogen chloride, hydrogen fluoride, nitric acid, phosphoric acid, sulfuric acid.

**Organic acid scan:** acetic acid, butyric acid, formic acid, propionic acid.

**Organo sulfur compounds:** carbon disulfide, dimethyl sulfide, etc.

**Solvent mixtures:** mineral spirits, naphthas, Stoddard solvents, TVOC, kerosene, etc.

## Alphabetical list of analytes

24	Acenaphthene (see PNA scan)	31	Amines, Aromatic (see scan for aromatic amines)
24	Acenaphthylene (see PNA scan)	31	Ammonia
24	Acetaldehyde	31	Ammonium Chloride Fume
24	Acetaldehyde	31	Amyl Acetate
25	Acetic Acid	32	Amyl Acetate
25	Acetic Acid	32	Anesthetic Gases scan
25	Acetic Anhydride	32	Aniline
25	Acetone	32	Anthracene (see PNA scan)
26	Acetone	33	Anthracene (see PNA scan)
26	Acetonitrile	33	Antimony and compounds as Sb
26	Acetonitrile	33	Antimony and compounds as Sb
26	Acrolein	33	Aromatic 100
27	Acrolein	34	Aromatic Amine scan
27	Acrolein	34	Arsenic and inorganic compounds, as As
27	Acrylamide	34	Arsenic and inorganic compounds, as As
27	Acrylic Acid	35	Arsenic and inorganic compounds, as As
28	Acrylonitrile (Vinyl Cyanide)	35	Arsenic Trioxide as As
28	Acrylonitrile (Vinyl Cyanide)	35	Arsine
28	Aldehyde scan	36	Arsine
28	Aliphatic Amine scan	36	Asbestos (bulk), all forms
29	Allyl Alcohol	36	Asbestos (Fibers)
29	Allyl Alcohol	36	Asbestos, all forms
29	Allyl Bromide	37	Asbestos, all forms
29	Allyl Bromide	37	Asphalt Fume
29	Allyl Chloride	37	Asphalt Fume as Benzene-Soluble Aerosol
30	Allyl Chloride	38	Barium and soluble compounds as Ba
30	Aluminum Metal and insoluble compounds	38	Benz[a]anthracene (see PNA scan)
30	Aluminum Oxide	38	Benzaldehyde
30	Amines, Aliphatic (see scan for aliphatic amines)	38	Benzaldehyde



39	Benzaldehyde	48	Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)
39	Benzene	48	Butyl Acrylate
39	Benzene	48	Butyl Acrylate
39	Benzo[a]pyrene (see PNA scan)	48	Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)
40	Benzo[a]pyrene (see PNA scan)	48	Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)
40	Benzo[b]fluoranthene (see PNA scan)	49	Butyl Cellosolve; (2-Butoxyethanol); (EGBE)
40	Benzo[e]pyrene	49	Butyl Cellosolve; (2-Butoxyethanol); (EGBE)
40	Benzo[ghi]perylene (see PNA scan)	49	Butyl(n-) Acetate
41	Benzo[k]fluoranthene (see PNA scan)	49	Butyl(n-) Acetate
41	Benzyl Alcohol	49	Butyl(n-) Alcohol
41	Benzyl Alcohol	50	Butyl(n-) Alcohol
41	Benzyl Chloride	50	Butyl(n-) Glycidyl Ether
41	Benzyl Chloride	50	Butyl(n-) Glycidyl Ether
42	Beryllium and compounds as Be	50	Butyl(sec-) Acetate
42	Beryllium and Compounds as Be	50	Butyl(sec-) Acetate
42	Beryllium and Compounds as Be	51	Butyl(sec-) Alcohol
43	Biphenyl (Diphenyl)	51	Butyl(sec-) Alcohol
43	Bismuth	51	Butyl(tert-) Acetate
43	Bisphenol A	51	Butyl(tert-) Acetate
43	Borate compounds, inorganic	51	Butyl(tert-) Alcohol
44	Breathing Air Grade D, Grade E	52	Butyl(tert-) Alcohol
44	Bromine	52	Butyraldehyde(n-)
44	Bromo(1-)-2-Chloroethane	52	Butyraldehyde(n-)
44	Bromoform	52	Butyraldehyde(n-)
45	Bromoform	53	Butyric Acid
45	Bromopropane(1-)	53	Cadmium and compounds as Cd
45	Bromopropane(1-)	53	Cadmium and compounds as Cd
45	Butadiene(1,3-)	54	Cadmium and compounds as Cd
46	Butadiene(1,3-)	54	Calcium
46	Butanedione(2,3-); (Butadione(2,3-), Diacetyl, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)	54	Calcium Carbonate
46	Butanone(2-); (Methyl Ethyl Ketone)	55	Calcium Carbonate
46	Butanone(2-); (Methyl Ethyl Ketone)	55	Calcium Hydroxide
47	Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)	55	Calcium Oxide
47	Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)	56	Calcium Oxide
47	Butoxyethoxy(2-(2-)) Ethanol	56	Calcium Silicate Synthetic Nonfibrous
47	Butoxyethoxy(2-(2-)) Ethyl Acetate	56	Calcium Sulfate (Gypsum)
47	Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)	56	Camphor

57	Camphor	66	Cobalt and Inorganic compounds as Co
57	Caprolactam	67	Cobalt and Inorganic compounds as Co
57	Carbaryl (SEVIN)	67	Cobalt and Inorganic compounds as Co
57	Carbon Black	67	Copper (Fume, Dusts and Mists) as Cu
58	Carbon Black	68	Copper (Fume, Dusts and Mists) as Cu
58	Carbon Disulfide	68	Copper (Fume, Dusts and Mists) as Cu
58	Carbon Disulfide	68	Cotton Dust, Raw
58	Carbon Tetrachloride (tetrachloromethane)	69	Cresol, all Isomers
59	Carbon Tetrachloride (tetrachloromethane)	69	Cumene
59	Cellosolve (2-Ethoxyethanol)	69	Cumene
59	Cellosolve (2-Ethoxyethanol)	69	Cyclohexane
59	Ceramic Fibers	70	Cyclohexane
60	Chlorine	70	Cyclohexanol
60	Chlorine Dioxide	70	Cyclohexanol
60	Chloro(2-)naphthalene	70	Cyclohexanone
60	Chlorobenzene	70	Cyclohexanone
61	Chlorobenzene	71	Cyclohexylamine
61	Chlorodiphenyl (Polychlorobiphenyl, 42% Chlorine)	71	Cyclopentane
61	Chlorodiphenyl (Polychlorobiphenyl, 54% Chlorine)	71	Cyclopentane
61	Chloroform (Trichloromethane)	71	Desflurane (Suprane)
62	Chloroform (Trichloromethane)	72	Desflurane (Suprane)
62	Chlorophenol(p-)	72	Diacetone Alcohol (4-Hydroxy-4-methyl-2-pentanone)
62	Chloroprene(beta-); (2-Chloro-1,3-butadiene)	72	Diacetone Alcohol (4-Hydroxy-4-methyl-2-pentanone)
62	Chloroprene(beta-); (2-Chloro-1,3-butadiene)	72	Diacetyl (Biacetyl, 2,3-Butadione, 2,3-Butanedione, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)
62	Chlorotoluene(o-)	73	Dibenzo[a,h]anthracene (see PNA scan)
63	Chlorotoluene(o-)	73	Dibromochloropropane (DBCP)
63	Chlorpyrifos (Dursban)	73	Dibutyl Ether
63	Chromium and Inorganic Compounds as Cr	73	Dibutyl Phthalate
64	Chromium and Inorganic compounds as Cr	74	Dichlorobenzene(o-)
64	Chromium and Inorganic compounds as Cr	74	Dichlorobenzene(o-)
64	Chromium, Hexavalent compounds as Cr	74	Dichlorobenzene(p-)
65	Chromium, Hexavalent Compounds as Cr	74	Dichlorobenzene(p-)
65	Chrysene (see PNA scan)	74	Dichloroethane(1,1)
65	Chrysene (see PNA scan)	74	Dichloroethane(1,1)
65	Coal Dust – Anthracite	75	Dichloroethane(1,1)
66	Coal Dust – Bituminous	75	Dichloroethylene(1,2-)(trans); (Acetylene dichloride)
66	Coal Tar Pitch Volatiles, as Benzene Soluble Aerosol		

75	Dichloroethylene(1,2-)(cis); (Acetylene dichloride)	82	Epichlorohydrin (1-Chloro- 2,3-epoxy propane)
75	Dichloroethylene(1,2-)(cis); (Acetylene dichloride)	82	Ethanolamine (2-Aminoethanol)
75	Dichloroethylene(1,2-)(trans); (Acetylene dichloride)	82	Ethoxyethanol(2-) (Cellosolve)
76	Dichloromethane (Methylene chloride)	83	Ethoxyethanol(2-) (Cellosolve)
76	Dichloromethane (Methylene chloride)	83	Ethoxyethyl(2-) Acetate (Cellosolve acetate)
76	Diesel Exhaust	83	Ethoxyethyl(2-) Acetate (Cellosolve acetate)
76	Diethanolamine	83	Ethyl 2-cyanoacrylate
77	Diethyl Ketone (3- Pentanone)	83	Ethyl 3-ethoxypropionate
77	Diethyl Ketone	84	Ethyl 3-ethoxypropionate
77	Diethyl Phthalate	84	Ethyl Acetate
77	Diethyl Sulfate	84	Ethyl Acetate
77	Diethylamine	84	Ethyl Acrylate
78	Diethylenetriamine	84	Ethyl Acrylate
78	Diglycidyl Ether of Bisphenol A	85	Ethyl Alcohol (Ethanol)
78	Dimethyl Acetamide	85	Ethyl Alcohol (Ethanol)
78	Dimethyl Acetamide	85	Ethyl Benzene
78	Dimethyl Disulfide	85	Ethyl Benzene
79	Dimethyl Sulfide	86	Ethyl Ether
79	Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)	86	Ethyl Ether
79	Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)	86	Ethyl Lactate
79	Dimethylformamide	86	Ethyl Methacrylate
79	Dimethylformamide	87	Ethylamine
80	Dioctyl Phthalate	87	Ethylene Chlorohydrin (2-Chloroethanol)
80	Dioxane(p-)	87	Ethylene Chlorohydrin (2-Chloroethanol)
80	Dioxane(p-)	87	Ethylene Dichloride (1,2-Dichloroethane)
80	Diphenyl (Biphenyl)	88	Ethylene Dichloride (1,2-Dichloroethane)
80	Dipropylene Glycol Methyl Ether (DPGME)	88	Ethylene Glycol
81	Dipropylene Glycol Methyl Ether (DPGME)	88	Ethylene Oxide
81	Dipropylene Glycol Methyl Ether Acetate (DPGMEA)	88	Ethylene Oxide
81	Divinyl Benzene	89	Ethylenediamine
81	Divinyl Benzene	89	Flour Dust
81	Enflurane (Ethrane)	89	Fluoranthene (see PNA scan)
82	Enflurane (Ethrane)	89	Fluorene (see PNA scan)
82	Epichlorohydrin (1-Chloro- 2,3-epoxy propane)	90	Fluorides, Particulate/Hydrogen Fluoride
		90	Forane (Isoflurane)
		90	Forane (Isoflurane)
		90	Formaldehyde
		91	Formaldehyde

91	Formaldehyde	101	Hydroquinone (Dihydroxybenzene)
91	Formamide	101	Hydroquinone (Dihydroxybenzene)
91	Formic Acid	101	Indium and Compounds as In
92	Furfural	101	Inorganic Acid Scan
92	Furfuryl Alcohol	102	Iodine and Iodides as I
92	Gasoline	102	Iron Oxide
92	Gasoline	102	Iron Oxide
93	Germanium	102	Iron
93	Glutaraldehyde	103	Isobutyl Acetate
93	Glutaraldehyde	103	Isobutyl Acetate
93	Glutaraldehyde	103	Isobutyl Alcohol
94	Gold	103	Isobutyl Alcohol
94	Grain Dust	104	Isocyanate Scan
94	Graphite	104	Isoflurane (Forane)
94	Halothane (Fluothane)	104	Isoflurane (Forane)
95	Halothane (Fluothane)	104	Isooctane
95	Heptane	105	Isooctane
95	Heptane	105	Isophorone
95	Heptanone(2-) (Methyl Amyl Ketone)	105	Isophorone
95	Heptanone(2-) (Methyl Amyl Ketone)	105	Isophorone Diisocyanate (IPDI)
96	Hexamethylene Diisocyanate (1,6- ) Homopolymer (HDI Homo)	106	Isophorone Diisocyanate (IPDI)
96	Hexamethylene Diisocyanate (1,6- ) Homopolymer (HDI Homo)	106	Isopropyl Acetate
96	Hexamethylene Diisocyanate (1,6-) (HDI)	106	Isopropyl Acetate
97	Hexamethylene Diisocyanate (1,6-) (HDI)	106	Isopropyl Alcohol (Isopropanol)
97	Hexane(n-)	107	Isopropyl Alcohol (Isopropanol)
97	Hexane(n-)	107	Kaolin
97	Hexyl Acrylate	107	Kerosene
98	Hexylene Glycol (2-Methyl-2,4-pentanediol)	107	Kerosene
98	Hydrazine	108	Lactic Acid
98	Hydrogen Bromide	108	Lanthanum
98	Hydrogen Bromide	108	Lead and Inorganic Compounds as Pb
99	Hydrogen Chloride	108	Lead and Inorganic Compounds as Pb
99	Hydrogen Chloride	109	Lead and Inorganic Compounds as Pb
99	Hydrogen Cyanide	109	Lead and Inorganic Compounds as Pb
100	Cyanide Salts as CN	109	Lead Chromate as Cr(VI)
100	Hydrogen Fluoride, as F	110	Limonene(d-)
100	Hydrogen Peroxide	110	Limonene(d-)
100	Hydrogen Sulfide	110	Lithium Salts
		110	Magnesium

111	Magnesium Oxide	118	Methyl Alcohol (Methanol)
111	Magnesium Oxide	118	Methyl Alcohol (Methanol)
111	Magnesium Oxide	119	Methyl Amyl Ketone (2-Heptanone)
111	Maleic Anhydride	119	Methyl Amyl Ketone (2-Heptanone)
112	Maleic Anhydride	119	Methyl Aniline
112	Manganese, Elemental and Inorganic compounds as Mn	119	Methyl Chloroform (1,1,1-Trichloroethane)
112	Manganese, Elemental and Inorganic compounds as Mn	119	Methyl Chloroform (1,1,1-Trichloroethane)
113	Manganese, elemental and Inorganic compounds as Mn	120	Methyl Cyclopentane
113	Medical Gases	120	Methyl Cyclopentane
113	Mercury as Hg (Elemental and inorganic forms)	120	Methyl Ethyl Ketone (2-Butanone, MEK)
114	Mercury as Hg (Elemental and inorganic forms)	120	Methyl Ethyl Ketone (2-Butanone, MEK)
114	Mercury as Hg Particulate	120	Methyl Isoamyl Ketone
114	Mesityl Oxide	121	Methyl Isoamyl Ketone
114	Metalworking Fluids	121	Methyl Isobutyl Ketone (MIBK)
115	Methanol (Methyl alcohol)	121	Methyl Isobutyl Ketone (MIBK)
115	Methanol (Methyl alcohol)	121	Methyl Isopropyl Ketone
115	Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)	121	Methyl Isopropyl Ketone
115	Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)	122	Methyl Methacrylate
116	Methoxyethanol(2-) (Methyl cellosolve, EGME)	122	Methyl Propyl Ketone (2-Pentanone)
116	Methoxyethanol(2-) (Methyl cellosolve, EGME)	122	Methyl Propyl Ketone (2-Pentanone)
116	Methoxyethoxy(2-(2-)) Ethanol (Diethylene glycol methyl ether)	122	Methyl Styrene(a-)
116	Methoxyethoxy(2-(2-)) Ethanol (Diethylene glycol methyl ether)	122	Methyl Styrene(a-)
116	Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)	123	Methyl Tert-butyl Ether (MTBE)
117	Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)	123	Methyl Tert-butyl ether (MTBE)
117	Methoxyethyl(2-) Ether (Diethylene glycol dimethyl ether)	123	Methyl Vinyl Ketone
117	Methoxyethyl(2-) Ether (Diethylene glycol dimethyl ether)	123	Methyl(1-)-2-pyrrolidinone
117	Methyl Acetate	124	Methyl(1-)-2-pyrrolidinone
117	Methyl Acetate	124	Methylacrylonitrile
118	Methyl Acrylate	124	Methylcyclohexane
118	Methyl Acrylate	124	Methylcyclohexane
		125	Methylene Bis(4-cyclohexylisocyanate) (HMDI)
		125	Methylene Bisphenyl Isocyanate (MDI)
		125	Methylene Bisphenyl Isocyanate (MDI)
		125	Methylene Chloride (Dichloromethane)
		126	Methylene Chloride (Dichloromethane)
		126	Methylene(4,4'-) Dianiline (MDA)
		126	Methylene(4,4'-)-bis(2 chloroaniline) (MOCA)
		126	Methylnaphthalene(2-)
		127	Mica

127	Mineral Oil (Oil mist)	136	Pentanedione(2,3-)
127	Mineral Oil, excluding Metal Working Fluids, Pure, highly and severely refined.	136	Pentanedione(2,4-)
127	Mineral Oil, used in metal working	137	Pentanone(2-) (Methyl propyl ketone)
128	Mineral Spirits (Stoddard Solvent)	137	Pentanone(2-) (Methyl propyl ketone)
128	Mineral Spirits (Stoddard Solvent)	137	Peracetic Acid
128	Mineral Wool Fiber	137	Perchloroethylene (Tetrachloroethylene)
128	Molybdenum as Mo	138	Perchloroethylene (Tetrachloroethylene)
129	Molybdenum as Mo	138	Perfluorooctanoic Acid
129	Molybdenum as Mo	138	Petroleum Ether
129	Morpholine	138	Petroleum Ether
129	Naphthalene	139	Phenanthrene (see PNA scan)
130	Naphthalene	139	Phenanthrene (see PNA scan)
130	Naphthalene (see PNA scan)	139	Phenol
130	Naproxen Sodium	139	Phenylcyclohexene(4-)
130	Naproxen Sodium	140	Phenylcyclohexene (4-)
131	Nickel and inorganic compounds as Ni	140	Phenylene(1,3-) diamine
131	Nickel and inorganic compounds as Ni	140	Phosphine
131	Nickel and inorganic compounds as Ni	140	Phosphoric Acid
131	Nicotine	141	Phosphoric Acid
132	Nitric Acid	141	Phosphorus (elements)
132	Nitric acid	141	Phthalic Anhydride
132	Nitric Oxide and Nitrogen Dioxide	141	Piperazine
132	Nitroethane	142	Platinum Metal and Soluble Salts as Pt
133	Nitrogen Dioxide	142	PNA Scan (NIOSH 5506)
133	Nitromethane	142	PNA Scan (OSHA 58)
133	Nitrous Oxide	142	Polychlorobiphenyl (Chlorodiphenyl, 54% Chlorine) (PCB)
133	Organic Solvent Scan	143	Polychlorobiphenyl (Chlorodiphenyl, 42% Chlorine) (PCB)
134	Oxalic Acid	143	Polyvinyl Chloride (PVC)
134	Ozone	143	Portland Cement
134	Palladium	143	Potassium Hydroxide
135	Paraffin Wax Fume	144	Propanol(n-)
135	Particles (insoluble or poorly soluble) Not otherwise specified; inhalable	144	Propanol(n-)
135	Particles (insoluble or poorly soluble) Not otherwise specified; respirable	144	Propionaldehyde
135	Particles (insoluble or poorly soluble) Not otherwise specified; total	144	Propionaldehyde
136	Pentane(n-)	145	Propionaldehyde
136	Pentane(n-)	145	Propionic Acid
		145	Propoxyethanol(2-) (Ethylene glycol monopropyl ether)

145	Propoxyethanol(2-) (Ethylene glycol monopropyl ether)	154	Silver Metal and Soluble Compounds as Ag
146	Propyl Bromide	154	Soapstone
146	Propyl Bromide	154	Soapstone
146	Propyl(n-) Acetate	155	Sodium
146	Propyl(n-) Acetate	155	Sodium Hydroxide
146	Propyl(n-) Alcohol	155	Starch
147	Propyl(n-) Alcohol	155	Stoddard Solvent
147	Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)	156	Stoddard Solvent
147	Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)	156	Strontium
147	Propylene Glycol Monomethyl Ether Acetate (PGMEA)	156	Strontium Chromate as Cr
147	Propylene Glycol Monomethyl Ether Acetate (PGMEA)	156	Styrene (Vinyl benzene)
148	Propylene Glycol (1,2-Propanediol)	157	Styrene (Vinyl benzene)
148	Propylene Oxide (1,2-Epoxypropane)	157	Sulfur Dioxide
148	Propylene Oxide (1,2-Epoxypropane)	157	Sulfuric Acid
148	Pyrene (see PNA scan)	157	Sulfuric Acid
149	Pyrene (see PNA scan)	157	Sulfuric Acid
149	Pyrethrum	158	Sulfuric Acid
149	Pyridine	158	Synthetic Vitreous Fibers
149	Resin Acids	158	Talc
150	Resorcinol	158	Tantalum and Tantalum Oxide Dust as Ta
150	Rhodium as Rh	159	Tellurium and Compounds as Te
150	Scan for Aldehydes	159	Tetrachloroethylene (Perchloroethylene)
150	Scan for Aliphatic Amines	159	Tetrachloroethylene (Perchloroethylene)
151	Scan for Anesthetic Gases	159	Tetrahydrofuran [THF]
151	Scan for Aromatic Amines	160	Tetrahydrofuran
151	Scan for Inorganic Acids	160	Thallium and Compounds, as Tl
151	Scan for Isocyanates	160	Thallium and Compounds, as Tl
152	Scan for Organic Solvents	160	Thiram
152	Scan for PNAs (NIOSH 5506)	161	Tin and Compounds as Sn
152	Scan for PNAs (OSHA 58)	161	Tin and Compounds as Sn
152	Selenium and Compounds as Se	161	Tin Organic Compounds as Sn
153	Selenium and Compounds as Se	162	Titanium
153	Sevoflurane (Sevofrane)	162	Titanium
153	Sevoflurane (Sevofrane)	162	Titanium Dioxide
153	Silica Cristobalite	163	Toluene
154	Silica Quartz	163	Toluene
		163	Toluene-2,4-diisocyanate (2,4-TDI)
		163	Toluene-2,4-diisocyanate (2,4-TDI)
		164	Toluene-2,6-diisocyanate (2,6-TDI)
		164	Toluene-2,6-diisocyanate (2,6-TDI)



164	Toluidine(o-)	170	Valeraldehyde
164	Tributyl Phosphate	170	Valeraldehyde
165	Trichloro(1,1,2-)-1,2,2-trifluoroethane	171	Vanadium Pentoxide as V
165	Trichloro(1,1,2-)-1,2,2-trifluoroethane	171	Vanadium Pentoxide as V
165	Trichlorobenzene(1,2,4-)	171	Vanadium
165	Trichlorobenzene(1,2,4-)	171	Vegetable Oil Mist
165	Trichloroethane(1,1,1-) (Methyl Chloroform)	172	Vinyl Acetate
166	Trichloroethane(1,1,1-) (Methyl Chloroform)	172	Vinyl Acetate
166	Trichloroethane(1,1,2-)	172	Vinyl Chloride (Chloroethylene)
166	Trichloroethane(1,1,2-)	172	Vinyl Chloride (Chloroethylene)
166	Trichloroethylene	173	Vinyl(1-)-2-pyrrolidinone
166	Trichloroethylene	173	Vinyl(1-)-2-pyrrolidinone
167	Triethanolamine	173	Vinylidene Chloride (1,1-Dichloroethylene)
167	Triethylamine	173	VM & P Naphtha
167	Triethylamine	174	VM & P Naphtha
167	Triethylenetetramine	174	Welding Fume Scan
168	Triglycidyl Isocyanurate(1,3,5)	174	Welding Fumes, Total
168	Trimellitic Anhydride	174	Wood Dust
168	Trimethylbenzene(1,2,4-)	175	Wood Dust
168	Trimethylbenzene(1,2,4-)	175	Xylene (Dimethyl benzene)
168	Trimethylbenzene(1,3,5-)	175	Xylene (Dimethyl benzene)
169	Trimethylbenzene(1,3,5-)	175	Yttrium and compounds, as Y
169	Tungsten and Compounds as W (in the absence of Cobalt)	176	Zinc
169	Tungsten, as W Soluble Compounds	176	Zinc
169	TVOC as n-Hexane	176	Zinc
170	TVOC as n-Hexane	177	Zinc Chloride Fume
170	Valeraldehyde	177	Zinc Oxide
		177	Zinc Oxide

**Acenaphthene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
83-32-9	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.33 µg	0.17 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

**Acenaphthylene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
208-96-8	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.45 µg	0.23 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

**Acetaldehyde**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-07-0	NIOSH 2016	HPLC	AT Monitor (N571AT)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
11.7	15-480	0.029 µg	0.015 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight.		

**Acetaldehyde**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-07-0	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	1-15	0.058 µg	0.029 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Sample at 1 lpm for STEL.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Acetic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media		
64-19-7	NIOSH 2011	IC	PTFE3-SGT** (SKC 225-17A, SKC 226-10-03)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.05-0.5	5-100	2.3 µg	1.2 µg	Acid2	
Interferences		Comments			
Particulate salts of the acid will give a positive interference.		Preferred method for STEL and also Lab-preferred method. Use 0.2 lpm for flow rate. Do not sample with inorganic acids.			

## Acetic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media		
64-19-7	OSHA ID-PV2119	IC	CT (SKC 226-01)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.2	5-100	2.3 µg	1.2 µg		
Interferences		Comments			

## Acetic Anhydride

CAS #	Analytical Method	Analytical Technique	Sampling Media		
108-24-7	OSHA 102	HPLC	GFF, vamine		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.05-0.5	20	14 µg	7.0 µg	ACN/DMSO	
Interferences		Comments			
		Sampling media has short shelf-life. Please contact the lab at least 5 days prior to survey date to order media. Media are prepared as requested.			

## Acetone

CAS #	Analytical Method	Analytical Technique	Sampling Media		
67-64-1	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2		0.5-3	0.91 µg	0.46 µg	CS <sub>2</sub>
Interferences		Comments			
		Store and ship cold overnight.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Acetone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-64-1	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
40.1	15-120	1.4 µg	0.70 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

## Acetonitrile

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-05-8	NIOSH 1606	GC-FID	CT (SKC 226-09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	3-25	1.0 µg	0.50 µg	15%MeOH/MC
Interferences		Comments		
Samples containing greater than 15% methanol or other alcohols		Large charcoal tubes are required for analyte collection since breakthrough volume is lower compared with smaller charcoal tubes. If also sampling for aldehydes, use AT monitors or Sep-Paks. Do not use DNPH tubes as they may off-gas acetonitrile.		

## Acetonitrile

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-05-8	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
48.2	15-120	1.5 µg	0.75 µg	DMF/CS <sub>2</sub> ; CS <sub>2</sub>
Interferences		Comments		
		Preferred method. Use 3M 3520. Separate front section of the monitor from the back section and cap immediately after sampling. If also sampling for aldehydes, use AT monitors or Sep-Paks. Do not use DNPH tubes as they may off-gas acetonitrile.		

## Acrolein

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-02-8	NIOSH 2016	HPLC	AT Monitor (N571AT)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
10.3	15-480	0.027 µg	0.014 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved		Refrigerate media before and after sampling. Ship samples cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Acrolein

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-02-8	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	1-15	0.054 µg	0.027 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use flow rate of 1 lpm.		

## Acrolein

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-02-8	NIOSH 2016	HPLC	Sep-Pak (WAT047205)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	10-100	0.14 µg	0.070 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use flow rate of 1 lpm.		

## Acrylamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-06-1	OSHA 21	GC-FID	GFF-SGT (SKC 225-16; SKC226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	120	0.62 µg	0.31 µg	MeOH
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents. 2019 NIC, TWA = 0.03ppm(IFV), Skin;DSEN;A2 adopted in 2020.		

## Acrylic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-10-7	OSHA PV2005	HPLC	Anasorb708/ Anasorb708 (SKC 226-30-08)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	24	0.33 µg	0.17 µg	
Interferences		Comments		
Propanoic acid will interfere.		Sample with 2 Anasorb 708 tubes in series. Separate and cap tubes after sampling.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Acrylonitrile (Vinyl Cyanide)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-13-1	OSHA 37	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	3.5-20	0.68 µg	0.34 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at 0.2 lpm.		

**Acrylonitrile (Vinyl Cyanide)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-13-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
43.8	15-480	1.0 µg	0.51 µg	DMF/CS <sub>2</sub> CS <sub>2</sub>
Interferences		Comments		

**Aldehyde scan**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 2016	HPLC		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
Interferences		Comments		
		See List of Scans for individual aldehydes. Refrigerate media before and after sampling. Ship samples cold overnight.		

**Aliphatic Amine scan**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 2010	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-1	5-30			Amine1
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.		Ethylamine, diethylamine and triethylamine are analyzed by this method. Please call Lab for other amines. High moisture will limit collection efficiency.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Allyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-18-6	NIOSH 1402	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	5-10	0.72 µg	0.36 µg	5% IPA/CS <sub>2</sub>
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents.		

### Allyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-18-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
40.4	15-480	1.1 µg	0.55 µg	MC or CS <sub>2</sub>
Interferences		Comments		

### Allyl Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-95-6	OSHA 1000	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-1.0	16-100	1.6 µg	0.80 µg	CS <sub>2</sub>
Interferences		Comments		

### Allyl Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-95-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.5	180-480	2.4 µg	1.2 µg	CS <sub>2</sub>
Interferences		Comments		

### Allyl Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-05-1	NIOSH 1000	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-1.0	16-100	0.96 µg	0.48 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 1 lpm.		

†(LPM) (CC/Min)

††(L) (Minutes)



## Allyl Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-05-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.1	15-480	1.4 µg	0.72 µg	CS <sub>2</sub>
Interferences		Comments		

## Aluminum Metal and insoluble compounds

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7429-90-5	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	200-1000	5.0 µg	2.5 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Aluminum Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1344-28-1	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	40-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 LPM, for area sampling up to 15 LPM.		

## Amines, Aliphatic (see scan for aliphatic amines)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 2010	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-1	5-30	Varies	Varies	Amine1
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.		Ethylamine, diethylamine and triethylamine are analyzed by this method. Please call Lab for other types of amines. High moisture will limit collection efficiency.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Amines, Aromatic (see scan for aromatic amines)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 2002	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.02-0.5	5-30	Varies	Varies	Amine3
Interferences		Comments		
Nitrogen compounds that co-elute will interfere		Aniline, methylaniline and o-toluidine are analyzed by this method. Please call Lab for other types of amines. High moisture will limit collection efficiency.		

### Ammonia

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7664-41-7	OSHA ID-188	IC	ORBO 77 (SUPELCO 20036) SKC 226-29	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-0.5	8-100	4.7 µg	2.4 µg	
Interferences		Comments		
Particulate ammonium salts will interfere.		Sample at a flow rate of 0.5 lpm for STEL.		

### Ammonium Chloride Fume

CAS #	Analytical Method	Analytical Technique	Sampling Media	
12125-02-9	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	40-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm; for area sampling up to 15 lpm.		

### Amyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
628-63-7	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	2.3 µg	1.2 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Amyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
628-63-7	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
26.0	15-480	3.5 µg	1.8 µg	CS <sub>2</sub>
Interferences		Comments		

## Anesthetic Gases scan

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05	12	Varies	Varies	
Interferences		Comments		
		Analyzes enflurane, halothane and isoflurane. See List of Scans for anesthetic gases. Please call lab for other anesthetic gases. Store and ship cold.		

## Aniline

CAS #	Analytical Method	Analytical Technique	Sampling Media	
62-53-3	NIOSH 2002	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.02-0.5	5-30	1.0 µg	0.50 µg	Amine3
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

## Anthracene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
120-12-7	OSHA 58	HPLC	GFF (SKC 225-7)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	960	0.25 µg	0.13 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Anthracene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media		
120-12-7	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
2		200-1000	0.25 µg	0.13 µg	PNAs
Interferences		Comments			
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.			

**Antimony and compounds as Sb**

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-36-0	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	50-1000	0.51 µg	0.26 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

**Antimony and compounds as Sb**

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-36-0	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		40-240	0.75 µg	0.38 µg	Metals
Interferences		Comments			
		As part of the Lab’s QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

**Aromatic 100**

CAS #	Analytical Method	Analytical Technique	Sampling Media		
64742-95-6	NIOSH 1550	GC-FID	CT (SKC 226-01, -09)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2		1.3-20	0.55 µg	0.28 µg	CS <sub>2</sub>
Interferences		Comments			
Other aromatic compounds.		Please send bulk sample. Ship bulk samples separately from air samples.			

†(LPM) (CC/Min)

††(L) (Minutes)

### Aromatic Amine scan

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 2002	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.02-0.5	5-30	Varies	Varies	Amine3
Interferences		Comments		
		Aniline, methylaniline and o-toluidine are analyzed by this method. Please call Lab for other types of amines.		

### Arsenic and inorganic compounds, as As

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-38-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	560-1000	0.56 µg	0.28 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Use ICP-MS analysis for lower detection limit at an additional cost. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

### Arsenic and inorganic compounds, as As

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-38-2	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	180-480	0.025 µg	0.012 µg	Metals
Interferences		Comments		
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. <b>Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)

## Arsenic and inorganic compounds, as As

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-38-2	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe (SKC 225-2414)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
NA		NA	6.5 µg	3.2 µg	Metals2
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

## Arsenic Trioxide as As

CAS #	Analytical Method	Analytical Technique	Sampling Media		
1327-53-3	NIOSH 7901	ICP-MS	MCE, carbonate		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-3	180-1000	0.033 µg	0.016 µg		
Interferences		Comments			
All forms of Arsenic are quantified.		<b>Media has one week shelf-life.</b> Media are prepared as requested. Please contact Lab in advance of survey date. As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. <b>Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.</b>			

## Arsine

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7784-42-1	NIOSH 6001	ICP	CT (226-01,-09)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2		10	0.052 µg	0.026 µg	
Interferences		Comments			
Other forms of Arsenic compounds (aerosol and gases) are quantified as Arsine.		High moisture may limit collection efficiency. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

†(LPM) (CC/Min)

††(L) (Minutes)

## Arsine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7784-42-1	NIOSH 6001	ICP-MS	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	10	0.026 µg	0.013 µg	
Interferences		Comments		
Other forms of Arsenic compounds (aerosol and gases) are quantified as Arsine.		High moisture may limit collection efficiency. As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. <b>Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.</b>		

## Asbestos (bulk), all forms

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1332-21-4	NIOSH 9002	PLM	Bulk	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
N.A.	1-10 grams	1.0%	1.0%	
Interferences		Comments		
Other fibers with optical properties similar to the asbestos minerals may give positive interferences.		Send bulk samples for asbestos analysis in double bagged ziplock bags with labels and chain of custody form outside the bag.		

## Asbestos (Fibers)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
Various	NIOSH 7400	PCM	MCE, 25 mm (ZEFON Z008BA)	
Sampling Rate†	Sampling Volume††	LOQ	LOQ	Compatibility Code
0.5-16	50-720	0.050 asb/field	0.01 asb/field	
Interferences		Comments		
Chain-like particles may appear fibrous leading to high count and high levels of non-fibrous dust particles may obscure fibers leading to low count.		Sample open faced. Adjust sampling flow rate and time to obtain optimum fiber loading on the filter. Do not overload filter. When shipping your samples, do not pack with untreated polystyrene as it can lead to fiber loss from electrostatic effect.		

## Asbestos, all forms

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1332-21-4	YAMATE LEVEL II	TEM	MCE 25-mm, 0.45-µm (ZEFON Z045BA)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2.48-17.62	74-529	0.050 asb/field	0.01 asb/field	
Interferences		Comments		
Other amphiboles similar to the asbestos minerals may give positive interferences. High dust background interferes with fiber identification.		Sample open faced. Adjust sampling flow rate and time to obtain optimum fiber loading on the filter. Do not overload filter. This analysis is sub-contracted to an AIHA-LAP, LLC accredited lab. <b>Standard turnaround time is 10 business days.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)



## Asbestos, all forms

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1332-21-4	NIOSH 7402	TEM	MCE, 25 mm (ZEFON Z008BA)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.5-16	50-720	0.050 asb/field	0.01 asb/field	
Interferences		Comments		
Other amphiboles similar to the asbestos minerals may give positive interferences. High dust background interferes with fiber identification.		Sample open faced. Adjust sampling flow rate and time to obtain optimum fiber loading on the filter. Do not overload filter. This method is designed for use with Method 7400 (phase contrast microscopy). This analysis is sub-contracted to an AIHA-LAP, LLC accredited lab. <b>Standard turnaround time is 10 business days.</b>		

## Asphalt Fume

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-42-4	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	40-7200	50 µg	10 µg	
Interferences		Comments		
All types of dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm; for area sampling up to 15 lpm.		

## Asphalt Fume as Benzene-Soluble Aerosol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-42-4	NIOSH 5042	GRAV	Pre-weighed PTFE1 (Zefon FPTFE137)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	176-960	130 µg	63 µg	
Interferences		Comments		
All substances soluble in benzene will interfere. Changes in humidity pre- and post weighing can affect accuracy.		Benzene extraction is done first. Scan is only done if BSF results are at or above the TLV. For inhalable sampling, please contact Lab for IOM sampler 1 week before intended use. The availability of IOM samplers is limited. Rental charges for use of IOM sampler apply.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Barium and soluble compounds as Ba

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-39-3	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	25-1000	0.50 µg	0.25 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

## Benz[a]anthracene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
56-55-3	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
2		200-1000	0.082 µg	0.041 µg	PNAs
Interferences		Comments			
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.			

## Benzaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media		
100-52-7	NIOSH 2016	HPLC	AT Monitor (N571AT)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
6.97	15-480	0.047 µg	0.024 µg	Aldehyde	
Interferences		Comments			
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight.			

## Benzaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media		
100-52-7	NIOSH 2016	HPLC	Sep-Pak (WAT047205)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.1-1.5	10-100	0.24 µg	0.12 µg	Aldehyde	
Interferences		Comments			
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use flow rate of 1 lpm for STEL			

†(LPM) (CC/Min)

††(L) (Minutes)

## Benzaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-52-7	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	1-15	0.094 µg	0.048 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use flow rate of 1 lpm for STEL		

## Benzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-43-2	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	5-30	0.40 µg	0.20 µg	CS <sub>2</sub>
Interferences		Comments		
		Sample at flow rate of 0.2 lpm for STEL.		

## Benzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-43-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.5	15-480	0.60 µg	0.30 µg	CS <sub>2</sub>
Interferences		Comments		

## Benzo[a]pyrene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-32-8	OSHA 58	HPLC	GFF (SKC 225-7)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	960	0.20 µg	0.10 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Benzo[a]pyrene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media		
50-32-8	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
2	200-1000	0.20 µg	0.10 µg	PNAs	
Interferences		Comments			
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.			

**Benzo[b]fluoranthene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media		
205-99-2	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
2	200-1000	0.19 µg	0.095 µg	PNAs	
Interferences		Comments			
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.			

**Benzo[e]pyrene**

CAS #	Analytical Method	Analytical Technique	Sampling Media		
192-97-2	NIOSH 5515	GC-MS	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
2	200-1000	0.57 µg	0.29 µg	PNAs	
Interferences		Comments			
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.			

**Benzo[ghi]perylene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media		
191-24-2	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
2	200-1000	0.27 µg	0.14 µg	PNAs	
Interferences		Comments			
Asphalt fumes will interfere		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.			

†(LPM) (CC/Min)

††(L) (Minutes)

**Benzo[k]fluoranthene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
207-08-9	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.25 µg	0.13 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

**Benzyl Alcohol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-51-6	NIOSH 1402	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	1-10	0.72 µg	0.36 µg	5%IPA/CS <sub>2</sub>
Interferences		Comments		

**Benzyl Alcohol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-51-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.1	15-480	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		

**Benzyl Chloride**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-44-7	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.73 µg	0.37 µg	CS <sub>2</sub>
Interferences		Comments		

**Benzyl Chloride**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-44-7	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.2	15-480	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Beryllium and compounds as Be

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-41-7	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	400-1000	0.010 µg	0.0050 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Minimum air volume required at 1/2 of TLV is 400 L. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b> Also indicate in your sample submission form if beryllium oxide is suspected to be present. Beryllium oxide requires analysis by OSHA ID-125G.		

## Beryllium and Compounds as Be

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-41-7	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	240	0.0012 µg	0.00059 µg	Metals
Interferences		Comments		
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. <b>Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.</b> Also indicate in your sample submission form if beryllium oxide is suspected to be present. Beryllium oxide requires analysis by OSHA ID-125G.		

## Beryllium and Compounds as Be

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-41-7	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe (SKC 225-2414)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
NA	NA	0.082 µg	0.041 µg	Metals2
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b> Also indicate in your sample submission form if beryllium oxide is suspected to be present. Beryllium oxide requires analysis by OSHA ID-125G.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Biphenyl (Diphenyl)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
92-52-4	OSHA PV2022	GC-FID	XAD-7 (SKC-226-95)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.20	20	0.46 µg	0.23 µg	CS <sub>2</sub>
Interferences		Comments		

**Bismuth**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-69-9	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	0.50 µg	0.25 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

**Bisphenol A**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
80-05-7	OSHA 1018	HPLC	GFF	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.0	240	0.32 µg	0.16 µg	
Interferences		Comments		
		Samples must be stored and shipped cold.		

**Borate compounds, inorganic**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
varies	NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25	varies	varies	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		LOQ for boron is 1.1 µg/sample. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b> OEL is as inhalable. Order IOM samplers a week before survey date. Rental charges for IOM samplers apply.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Breathing Air Grade D, Grade E

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	NIOSH 0500 LMI-A5	GRAV GC-FID GC-ECD GC-TCD GC-XSD	PTFE4 (PALL TF-450) Cylinder		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
NA	40-1000 (gravimetric) 300cc/25 psig (Cylinder)	Note 1	–	BA	
Interferences		Comments			
At high levels argon interferes with oxygen and nitrous oxide interferes with carbon dioxide.		A particulate blank is required.			

## Bromine

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7726-95-6	NIOSH 6011	IC	PTFE-AgMF (SKC 225-1708, SKC 225-1802)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.3-1.0	70-360	4.5 µg	3.0 µg	Cl <sub>2</sub> &Br <sub>2</sub>	
Interferences		Comments			
Hydrobromic acid will interfere.		Preferred method.			

## Bromo(1-)-2-Chloroethane

CAS #	Analytical Method	Analytical Technique	Sampling Media		
107-04-0	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2		1-10	1.0 µg	0.50 µg	CS <sub>2</sub>
Interferences		Comments			

## Bromoform

CAS #	Analytical Method	Analytical Technique	Sampling Media		
75-25-2	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2		1-50	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments			

†(LPM) (CC/Min)

††(L) (Minutes)



**Bromoform**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-25-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
29.3	15-480	1.7 µg	0.85 µg	CS <sub>2</sub>
Interferences		Comments		

**Bromopropane(1-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-94-5	OSHA PV2061	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	12	1.3 µg	0.65 µg	CS <sub>2</sub>
Interferences		Comments		

**Bromopropane(1-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-94-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.7	15-480	2.0 µg	1.0 µg	CS <sub>2</sub>
Interferences		Comments		

**Butadiene(1,3-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-99-0	NIOSH 1024	GC-MS	CT-CT (SKC 226-01)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	10-25	1.1 µg	0.55 µg	MC
Interferences		Comments		
High humidity (>80%) or other hydrocarbons present at permissible levels decrease sampler's capacity.		Use two large charcoal tubes in series. Separate and cap tubes after sampling. Ship cold overnight. If unable to ship overnight, store cold and then ship cold the following day.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Butadiene(1,3-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-99-0	3M Method	GC-MS	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
42.8	15-480	1.7 µg	0.85 µg	MC CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate front section of the monitor from the back section and cap immediately after sampling. Ship cold overnight. If unable to ship cold overnight, refrigerate, and then ship cold the following day.		

**Butanedione(2,3-); (Butadione(2,3-), Diacetyl, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
431-03-8	OSHA 1013	GC-MS	SGT/GFF-SGT/GFF (SKC 226-183)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	9 (TWA) 3 (15-min short term)	0.16 µg	0.080 µg	95% EtOH
Interferences		Comments		
		Samples are collected on two specially washed silica gel tubes in series. Samples should be protected from the light during and after sampling. Separate and cap tubes after sampling. Sample separately from CS <sub>2</sub> compatible solvents. Order aluminum foil for wrapping the samples. For STEL sampling, a minimum of 3.0L is required. For TWA sampling, a minimum of 9.0L is required. Store and ship cold overnight.		

**Butanone(2-); (Methyl Ethyl Ketone)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-93-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
36.3	15-480	1.2 µg	0.59 µg	CS <sub>2</sub>
Interferences		Comments		
Isopropyl acetate may co-elute with MEK.		Ship and store cold.		

**Butanone(2-); (Methyl Ethyl Ketone)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-93-3	NIOSH 2500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.25-12	0.79 µg	0.40 µg	CS <sub>2</sub>
Interferences		Comments		
Isopropyl acetate may co-elute with MEK.		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-76-2	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	2-10	1.2 µg	0.60 µg	5%MeOH/MC
Interferences		Comments		
		Sample separately from other solvents.		

**Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-76-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.2	15-480	1.8 µg	0.90 µg	MC CS <sub>2</sub>
Interferences		Comments		
		Sample separately from other solvents.		

**Butoxyethoxy(2-(2-)) Ethanol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-34-5	OSHA PV2095	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	2-10	0.58 µg	0.29 µg	5%MeOH/MC
Interferences		Comments		
		Sample separately from other solvents.		

**Butoxyethoxy(2-(2-)) Ethyl Acetate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
124-17-4	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.87 µg	0.44 µg	CS <sub>2</sub>
Interferences		Comments		

**Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-07-2	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.94 µg	0.47 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-07-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.3	15-480	1.4 µg	0.70 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl Acrylate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-32-3	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.69 µg	0.35 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl Acrylate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-32-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.3	15-480	1.0 µg	0.50 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-07-2	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.94 µg	0.47 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-07-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.3	15-480	1.4 µg	0.70 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Butyl Cellosolve; (2-Butoxyethanol); (EGBE)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-76-2	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	2-10	1.2 µg	0.60 µg	5%MeOH/MC
Interferences		Comments		

**Butyl Cellosolve; (2-Butoxyethanol); (EGBE)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-76-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.2	15-480	1.8 µg	0.90 µg	MC
Interferences		Comments		
		Sample separately from other solvents.		

**Butyl(n-) Acetate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-86-4	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.84 µg	0.42 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl(n-) Acetate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-86-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.6	15-480	1.3 µg	0.65 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl(n-) Alcohol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-36-3	NIOSH 1401	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	2-10	0.67 µg	0.34 µg	1%IPA/CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Butyl(n-) Alcohol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-36-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
34.3	15-480	1.0 µg	0.50 µg	MC
Interferences		Comments		
		Sample separately from other solvents.		

**Butyl(n-) Glycidyl Ether**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2426-08-6	NIOSH 1616	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	15-30	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl(n-) Glycidyl Ether**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2426-08-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.0	15-480	1.7 µg	0.85 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl(sec-) Acetate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
105-46-4	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.71 µg	0.36 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl(sec-) Acetate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
105-46-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.6	15-480	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Butyl(sec-) Alcohol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-92-2	NIOSH 1401	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	2-10	1.0 µg	0.50 µg	1%IPA/CS <sub>2</sub>
Interferences		Comments		

**Butyl(sec-) Alcohol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-92-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
34.8	15-480	1.5 µg	0.75 µg	MC CS <sub>2</sub>
Interferences		Comments		

**Butyl(tert-) Acetate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-88-5	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.84 µg	0.42 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl(tert-) Acetate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-88-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
29.4	15-480	1.3 µg	0.65 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyl(tert-) Alcohol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-65-0	NIOSH 1400	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	2-10	0.55 µg	0.28 µg	BUT/CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Butyl(tert-) Alcohol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-65-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.2	15-480	0.83 µg	0.42 µg	CS <sub>2</sub>
Interferences		Comments		

**Butyraldehyde(n-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-72-8	NIOSH 2016	HPLC	AT Monitor (N571AT)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
8.20	15-480	0.037 µg	0.019 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight.		

**Butyraldehyde(n-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-72-8	NIOSH 2016	HPLC	Sep-Pak (WAT047205)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	10-100	0.19 µg	0.095 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. <b>Preferred for STEL sampling.</b> Sample at 1.0 lpm for STEL.		

**Butyraldehyde(n-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-72-8	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	1-15	0.073 µg	0.037 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. <b>Preferred for STEL sampling.</b> Sample at 1.0 lpm for STEL.		

†(LPM) (CC/Min)

††(L) (Minutes)



## Butyric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media		
107-92-6	NIOSH 2011	IC	PTFE3-SGT** (SKC 225-17A, SKC 226-10-03)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.5		15-100	2.3 µg	1.2 µg	Acid2
Interferences		Comments			
		Preferred method for STEL and also Lab-preferred method. Use 0.2 lpm flow rate. DO NOT sample with inorganic acids.			

## Cadmium and compounds as Cd

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-43-9	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	500-1000	0.10 µg	0.050 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

## Cadmium and compounds as Cd

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-43-9	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe (SKC 225-2414)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
NA		NA	0.50 µg	0.25 µg	Metals2
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab’s QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

†(LPM) (CC/Min)

††(L) (Minutes)

## Cadmium and compounds as Cd

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-43-9	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		160-240	0.032 µg	0.016 µg	Metals
Interferences		Comments			
		As part of the Lab’s QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

## Calcium

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-70-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	100-1000	7.9 µg	4.0 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of calcium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.			

## Calcium Carbonate

CAS #	Analytical Method	Analytical Technique	Sampling Media		
1317-65-3	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	100-1000	20 µg	10 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of calcium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. If client wants different calcium salts speciated, please indicate in sample submission form or call/email Lab.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Calcium Carbonate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1317-65-3	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	20-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

## Calcium Hydroxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1305-62-0	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	100-1000	15 µg	7.5 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of calcium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. If client wants different calcium salts speciated, please indicate in sample submission form or call/email Lab.</b>		

## Calcium Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1305-78-8	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	100-1000	11 µg	5.5 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of calcium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. If client wants different calcium salts speciated, please indicate in sample submission form or call/email Lab.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)

## Calcium Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media		
1305-78-8	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	40-240	4.6 µg	2.3 µg	Metals	
Interferences		Comments			
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

## Calcium Silicate Synthetic Nonfibrous

CAS #	Analytical Method	Analytical Technique	Sampling Media		
1344-95-2	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-15		40-7200	50 µg	10 µg	
Interferences		Comments			

## Calcium Sulfate (Gypsum)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7778-18-9	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-15	20-7200	50 µg	10 µg		
Interferences		Comments			
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.			

## Camphor

CAS #	Analytical Method	Analytical Technique	Sampling Media		
76-22-2	NIOSH 1301	GC-FID	CT (SKC 226-01, -09)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2		1-25	0.53 µg	0.27 µg	MeOH/CS <sub>2</sub>
Interferences		Comments			
		Preferred method. Sample at 0.2 lpm for STEL. Sample separately from CS <sub>2</sub> compatible solvents.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Camphor

CAS #	Analytical Method	Analytical Technique	Sampling Media	
76-22-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
21.4	15-480	0.80 µg	0.40 µg	CS <sub>2</sub>
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents.		

## Caprolactam

CAS #	Analytical Method	Analytical Technique	Sampling Media	
105-60-2	OSHA PV2012	HPLC	OVS 7 (SKC 226-57)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	100	2.1 µg	1.1 µg	
Interferences		Comments		

## Carbaryl (SEVIN)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
63-25-2	OSHA 63	HPLC	OVS-2 (SKC 226-30-16)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	10-60	0.45 µg	0.23 µg	
Interferences		Comments		

## Carbon Black

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1333-86-4	NIOSH 5000	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-2	30-570	50 µg	30 µg	
Interferences		Comments		
All other dusts will interfere.		Preferred method. For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Carbon Black

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1333-86-4	OSHA ID-196	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	460-960	850 µg	420 µg	
Interferences		Comments		
Particulates that are insoluble in THF and that either vaporize or lose weight between 150°C and 600°C will interfere.		Please notify lab prior to sample collection.		

## Carbon Disulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-15-0	NIOSH 1600	GC-MS	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	2-10	0.52 µg	0.26 µg	Tol
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents. Store and ship cold.		

## Carbon Disulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-15-0	3M Method	GC-MS	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
42.8	15-480	0.78 µg	0.39 µg	Tol
Interferences		Comments		
		Use 3M 3520. Separate front section of the monitor from the back section and cap immediately after sampling. Sample separately from CS <sub>2</sub> compatible solvents. Store and ship cold.		

## Carbon Tetrachloride (tetrachloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
56-23-5	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1.5-40	11 µg	5.5 µg	CS <sub>2</sub>
Interferences		Comments		
		Use a flow rate of 0.2 lpm for STEL. Preferred method for STEL sampling.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Carbon Tetrachloride (tetrachloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
56-23-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.2	15-480	17 µg	8.5 µg	CS <sub>2</sub>
Interferences		Comments		

### Cellosolve (2-Ethoxyethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-80-5	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1-6	0.54 µg	0.27 µg	5%MeOH/MC
Interferences		Comments		
		Sample separately from other solvents.		

### Cellosolve (2-Ethoxyethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-80-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.4	15-480	0.81 µg	0.41 µg	MC CS <sub>2</sub>
Interferences		Comments		
		Sample separately from other solvents.		

### Ceramic Fibers

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 7400	PCM	MCE, 25 mm (ZEFON Z008BA)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.5-16	50-720	0.050 fiber/ field	0.01 fiber/field	
Interferences		Comments		
Chain-like particles which may appear fibrous cause positive interference. High levels of non-fibrous dust particles may obscure fibers.		Adjust sampling flow rate and time to obtain optimum fiber loading on the filter. Do not overload filter. Sample open faced. When shipping your samples, do not pack with untreated polystyrene as can lead to fiber loss from electrostatic effect.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Chlorine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7782-50-5	NIOSH 6011	IC	PTFE-AgMF (SKC 225-1708, SKC 225-1802)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.3-1.0	80-360	2.3 µg	1.2 µg	Cl <sub>2</sub> &Br <sub>2</sub>
Interferences		Comments		
Hydrochloric acid will interfere.		Preferred method. Use a flow rate of 1.0 lpm for STEL. Order media one week ahead, media is prepared when ordered.		

## Chlorine Dioxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10049-04-4	OSHA ID-202	IC	Impinger 4	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.5	35-120	0.75 µg	0.38 µg	
Interferences		Comments		
		Shelf-life of impinger solution is 3 months. Transfer sample solution to labeled vials before shipping. To avoid "freezing" of glass to glass, rinse impinger with distilled water before returning to the lab.		

## Chloro(2-)naphthalene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
91-58-7	NIOSH 5515	GC-MS	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.60 µg	0.30 µg	PNAs
Interferences		Comments		
		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

## Chlorobenzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-90-7	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1.5-40	0.57 µg	0.29 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)



## Chlorobenzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-90-7	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
29.3	15-480	0.86 µg	0.43 µg	CS <sub>2</sub>
Interferences		Comments		

## Chlorodiphenyl (Polychlorobiphenyl, 42% Chlorine)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
53469-21-9	NIOSH 5503	GC-MS	GFF-Florisil (Millipore SX0001300/01 /AP2001300 SKC 226-39)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	1-50	0.97 µg	0.49 µg	
Interferences		Comments		
Other chlorinated pesticides may interfere in the quantification of PCB.				

## Chlorodiphenyl (Polychlorobiphenyl, 54% Chlorine)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
11097-69-1	NIOSH 5503	GC-MS	GFF-Florisil(Millipore SX0001300/01/ AP2001300 SKC 226-39)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	1-50	1.1 µg	0.55 µg	
Interferences		Comments		
Other chlorinated pesticides may interfere in the quantification of PCB.				

## Chloroform (Trichloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-66-3	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-50	5.9 µg	3.0 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Chloroform (Trichloromethane)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-66-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
33.5	15-480	8.9 µg	4.5 µg	CS <sub>2</sub>
Interferences		Comments		

**Chlorophenol(p-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-48-9	NIOSH 2014	HPLC	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	10-40	0.10 µg	0.050 µg	
Interferences		Comments		

**Chloroprene(beta-); (2-Chloro-1,3-butadiene)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
126-99-8	NIOSH 1002	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.1	1.5-8	0.43 µg	0.22 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold.		

**Chloroprene(beta-); (2-Chloro-1,3-butadiene)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
126-99-8	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.2	15-480	0.65 µg	0.33 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold.		

**Chlorotoluene(o-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-49-8	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1.5-40	0.47 µg	0.24 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

### Chlorotoluene(o-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-49-8	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.3	15-480	0.70 µg	0.35 µg	CS <sub>2</sub>
Interferences		Comments		

### Chlorpyrifos (Dursban)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2921-88-2	NIOSH 5600	GC-NPD	OVS-2/QF (SKC 226-58)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2-1	12-240	0.039 µg	0.020 µg	
Interferences		Comments		
Organophosphate compounds may interfere.				

### Chromium and Inorganic Compounds as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-47-3	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	670-1000	1.0 µg	0.50 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Minimum air volume required at ½ of TLV is 670 L. This applies to Cr III compounds. See chromium VI and chromate entries for special instructions on collection for this form of chromium. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)

## Chromium and Inorganic compounds as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-47-3	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		240-500	0.75 µg	0.38 µg	Metals
Interferences		Comments			
		Minimum air volume required at ½ of TLV is 500 L. This applies to Cr III compounds. See chromium VI and chromate entries for special instructions on collection for this form of chromium. As part of the Lab’s QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. <b>Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.</b>			

## Chromium and Inorganic compounds as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-47-3	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe (SKC 225-2414)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
NA		NA	1.0 µg	0.50 µg	Metals2
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

## Chromium, Hexavalent compounds as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-47-3	OSHA ID-215	IC	PVC (SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
2		500-960	0.050 µg	0.020 µg	
Interferences		Comments			
Fe (II) appears to cause a negative interference during sampling and storage.		Please indicate in sample submission form the operation sampled (e.g., spray paint, chrome plating, welding, etc.) as the extraction method is different for spray paint samples. Refrigerate samples and ship <i>overnight</i> as soon as possible. 500 liters is the minimum air volume at 50% TLV. Samples from plating operations must be analyzed within 6 days from date of sampling. Samples from welding operations must be analyzed within 8 days from date of sampling. If sampling CrVI as inhalable, order IOM samplers a week ahead of survey date. Remember IOM samplers are not stock items to be kept indefinitely.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Chromium, Hexavalent Compounds as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-47-3	OSHA W4001	IC	QFF (Millipore AQFA03700) Wipes	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
NA	NA	0.050 µg	0.020 µg	
Interferences		Comments		
Fe (II) appears to cause a negative interference during sampling and storage.		Please follow "Hexavalent Chromium Wipe Sampling Instructions." Do not use ghost wipes, Whatman, mixed cellulose ester (MCE) or glass fiber filters as they convert Cr(VI) to Cr(III). Please indicate in sample submission form the operation sampled as the extraction method is different for spray paint samples. Refrigerate samples and ship <i>overnight</i> as soon as possible. If sampling a chromium plating operation, request vials containing 10% Na <sub>2</sub> CO <sub>3</sub> with 2%NaHCO <sub>3</sub> to stabilize the samples. Samples from welding operations must be analyzed within 8 days from date of sampling.		

## Chrysene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
218-01-9	OSHA 58	HPLC	GFF (SKC 225-7)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	960	0.11 µg	0.055 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.		

## Chrysene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
218-01-9	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.11 µg	0.055 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

## Coal Dust – Anthracite

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0600	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.7	50-816	50 µg	10 µg	
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Coal Dust – Bituminous

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0600	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.7	50-816	50 µg	10 µg	
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

### Coal Tar Pitch Volatiles, as Benzene Soluble Aerosol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
65996-93-2	OSHA 58	GRAV	GFF (SKC 225-7)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	960	130 µg	28 µg	
Interferences		Comments		
The method is non-specific and measures all substances soluble in benzene.		This method is used to air monitor coke oven emissions, petroleum combustion products & asphalt fumes. If the BSF exceeds the PEL, then the sample is analyzed by HPLC to determine the presence of selected polynuclear aromatic hydrocarbons (PAHs). Wrap sample cassettes with aluminum foil.		

### Cobalt and Inorganic compounds as Co

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-48-4	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	140-1000	0.10 µg	0.050 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)

### Cobalt and Inorganic compounds as Co

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-48-4	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		75-240	0.025 µg	0.013 µg	Metals
Interferences		Comments			
		As part of the Lab’s QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

### Cobalt and Inorganic compounds as Co

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-48-4	OSHA ID-121	ICP	Ghost wipe ( 225-2414)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
NA		NA	0.5 µg	0.25 µg	Metals 2
Interferences		Comments			
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab’s QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.			

### Copper (Fume, Dusts and Mists) as Cu

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-50-8	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	70-1000	1.4 µg	0.70 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

†(LPM) (CC/Min)

††(L) (Minutes)

### Copper (Fume, Dusts and Mists) as Cu

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-50-8	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		40-240	0.10 µg	0.050 µg	Metals
Interferences		Comments			
		As part of the Lab’s QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

### Copper (Fume, Dusts and Mists) as Cu

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-50-8	OSHA ID-121	ICP	Ghost wipe ( 225-2414)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
NA		NA	5.0 µg	2.5 µg	Metals 2
Interferences		Comments			
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab’s QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.			

### Cotton Dust, Raw

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	NIOSH 0500	GRAV	PVC See Comment		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
See Comment	1000-3600	50 µg	10 µg		
Interferences		Comments			
All other dusts will interfere.		OSHA: Open-faced sampling with 3-piece pre-weighed PVC cassette on a vertical elutriator, cotton-dust sampler at 7.4 lpm. ACGIH: Use BGI GK2.69 cyclone at 1.6 lpm; sample 768 liters for 65% of new thoracic TLV.			

†(LPM) (CC/Min)

††(L) (Minutes)



### Cresol, all Isomers

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1319-77-3 95-48-7 108-39-4 106-44-5	OSHA 32	HPLC	XAD-7 (SKC 226-95)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	5-24	0.39 µg	0.20 µg	Phenol & cresol
Interferences		Comments		
		This method is applicable for all isomers of cresol (ortho-, meta-, and para-).		

### Cumene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-82-8	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01 -0.2	2-30	0.55 µg	0.28 µg	CS <sub>2</sub>
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%. 2018 NIC of TWA to 1ppm, A3 was changed to TWA = 5ppm, A3 in 2020.		

### Cumene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-82-8	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.5	15-480	0.83 µg	0.42 µg	CS <sub>2</sub>
Interferences		Comments		
		2018 NIC of TWA to 1ppm, A3 was changed to TWA = 5ppm, A3 in 2020.		

### Cyclohexane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-82-7	NIOSH 1500	GC-MS	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	2.5-5	0.17 µg	0.083 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Cyclohexane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-82-7	3M Method	GC-MS	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.4	15-360	0.25 µg	0.13 µg	CS <sub>2</sub>
Interferences		Comments		

## Cyclohexanol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-93-0	NIOSH 1402	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	1.1 µg	0.55 µg	5%PRO/CS <sub>2</sub>
Interferences		Comments		

## Cyclohexanol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-93-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
29.5	15-480	1.7 µg	0.85 µg	MC
Interferences		Comments		
		Sample separately from other solvents.		

## Cyclohexanone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-94-1	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.72 µg	0.36 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

## Cyclohexanone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-94-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.9	15-480	1.1 µg	0.54 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Cyclohexylamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-91-8	NIOSH 2010	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-1	5-30	5.9 µg	3.0 µg	Amine1
Interferences		Comments		
Methanol could interfere in low level analysis.				

## Cyclopentane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
287-92-3	NIOSH 1500	GC-FID	CT (SKC 226-01)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	2.5-5.0	0.40 µg	0.20 µg	CS <sub>2</sub>
Interferences		Comments		

## Cyclopentane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
287-92-3	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
36.2	15-60	0.60 µg	0.30 µg	CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate front section of the monitor from the back section and cap immediately after sampling.		

## Desflurane (Suprane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
57041-67-5	OSHA 106	GC-FID	Anasorb 747 (SKC 226-81A)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05	3	6.1 µg	3.1 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Desflurane (Suprane)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
57041-67-5	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.1	15-480	9.1 µg	4.6 µg	CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate front section of the monitor from the back section and cap immediately after sampling. Store and ship cold overnight.		

**Diacetone Alcohol (4-Hydroxy-4-methyl-2-pentanone)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-42-2	NIOSH 1402	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.98 µg	0.49 µg	5%PRO/CS <sub>2</sub>
Interferences		Comments		

**Diacetone Alcohol (4-Hydroxy-4-methyl-2-pentanone)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-42-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.2	15-480	1.5 µg	0.75 µg	MC
Interferences		Comments		
		Sample separately from other solvents.		

**Diacetyl (Biacetyl, 2,3-Butadione, 2,3-Butanedione, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
431-03-8	OSHA 1013	GC-MS	SGT/GFF-SGT/GFF (SKC 226-183)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	9 (TWA) 3 (15-min short term)	0.16 µg	0.080 µg	95%EtOH
Interferences		Comments		
		Samples are collected on two specially washed silica gel tubes in series. Samples should be protected from the light during and after sampling. Order aluminum foil for wrapping the samples during and after sampling. Separate and cap tubes after sampling. Sample separately from CS <sub>2</sub> compatible solvents. For STEL sampling, a minimum of 3.0L is required. For TWA sampling, a minimum of 9.0L is required. Store and ship cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Dibenzo[a,h]anthracene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
53-70-3	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.93 µg	0.47 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

**Dibromochloropropane (DBCP)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-12-8	OSHA In-house	GC-FID	Anasorb 747/ Anasorb 747 (SKC 226-81A)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.20	20	0.41 µg	0.205 µg	CS <sub>2</sub>
Interferences		Comments		
		Sample with 2 Anasorb 747 tubes in series. Separate and cap tubes after sampling. Ship cold overnight.		

**Dibutyl Ether**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
142-96-1	NIOSH 1610	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	1-10	0.85 µg	0.43 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold.		

**Dibutyl Phthalate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
84-74-2	NIOSH 5020	GC-FID	MCE (SKC 225-5)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-3	10-200	0.46 µg	0.23 µg	CS <sub>2</sub>
Interferences		Comments		
		Sample at 1.0 lpm for STEL.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Dichlorobenzene(o-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-50-1	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.50 µg	0.25 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

**Dichlorobenzene(o-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-50-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.8	15-480	0.75 µg	0.38 µg	CS <sub>2</sub>
Interferences		Comments		

**Dichlorobenzene(p-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-46-7	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-8	0.45 µg	0.23 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

**Dichlorobenzene(p-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-46-7	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.8	15-480	0.68 µg	0.34 µg	CS <sub>2</sub>
Interferences		Comments		

**Dichloroethane(1,1)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-34-3	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.5-15	1.7 µg	0.83 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Dichloroethane(1,1)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-34-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
33.2	15-480	2.5 µg	1.3 µg	CS <sub>2</sub>
Interferences		Comments		

**Dichloroethylene(1,2-)(trans); (Acetylene dichloride)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-59-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.2	15-360	2.1 µg	1.1 µg	CS <sub>2</sub>
Interferences		Comments		

**Dichloroethylene(1,2-)(cis); (Acetylene dichloride)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-59-0	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.2-5	1.6 µg	0.80 µg	CS <sub>2</sub>
Interferences		Comments		

**Dichloroethylene(1,2-)(cis); (Acetylene dichloride)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-59-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.2	15-480	2.4 µg	1.2 µg	CS <sub>2</sub>
Interferences		Comments		

**Dichloroethylene(1,2-)(trans); (Acetylene dichloride)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-59-0	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	15-360	1.4 µg	0.70 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Dichloromethane (Methylene chloride)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-09-2	NIOSH 1005	GC-FID	CT-CT (SKC 226-01)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.5-2.5	3.7 µg	1.9 µg	CS <sub>2</sub>
Interferences		Comments		
		Separate tubes and cap immediately after sampling. Ship cold immediately.		

**Dichloromethane (Methylene chloride)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-09-2	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
37.9	15-240	5.6 µg	2.8 µg	CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate front section of the monitor from the back section and cap immediately after sampling. Ship cold immediately.		

**Diesel Exhaust**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 5040	EGA-TDA	QFF	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2-4	142-19000	1.3 µg	NA	
Interferences		Comments		
		This analysis is sub-contracted to an AIHA-LAP, LLC accredited lab. We require a week's notice to procure the media. Media has short shelf-life. For underground mines, use diesel particulate matter (DPM) cassettes (SKC 225-317). Turnaround time is 10 working days.		

**Diethanolamine**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-42-2	NIOSH 2007	IC	ORBO 53 or SGT** (SUPELCO 20265) (SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.5	30	3.0 µg	1.5 µg	EA
Interferences		Comments		
		Store in freezer after sampling. Ship cold.		

†(LPM) (CC/Min)

††(L) (Minutes)



### Diethyl Ketone (3- Pentanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-22-0	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.67 µg	0.34 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

### Diethyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-22-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.7	15-480	1.0 µg	0.50 µg	CS <sub>2</sub>
Interferences		Comments		

### Diethyl Phthalate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
84-66-2	OSHA 104	GC-FID	OVS-Tenax (SKC 226-56)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	240	0.58 µg	0.29 µg	Tol
Interferences		Comments		

### Diethyl Sulfate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-67-5	NIOSH 2524	GC-FID	Porapak-P (SKC 226-114)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	0.25-12	7.2 µg	3.6 µg	Ethyl Ether
Interferences		Comments		

### Diethylamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-89-7	NIOSH 2010	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-1.0	5-30	0.31 µg	0.16 µg	Amine1
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

†(LPM) (CC/Min)

††(L) (Minutes)

## Diethylenetriamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-40-0	OSHA 60	HPLC	XAD-2, NITC (SKC 226-30-18)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	10	0.10 µg	0.050 µg	Amine2
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

## Diglycidyl Ether of Bisphenol A

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1675-54-3	OSHA 1018	HPLC	GFF	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.5	270	0.50 µg	0.25 µg	
Interferences		Comments		

## Dimethyl Acetamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-19-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.0	15-480	0.99 µg	0.50 µg	MC
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents.		

## Dimethyl Acetamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-19-5	NIOSH 2004	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-1	15-80	0.66 µg	0.33 µg	MeOH
Interferences		Comments		
		Silica gel has a high affinity for water; high relative humidity may adversely affect the efficiency of analyte adsorption.		

## Dimethyl Disulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
624-92-0	LM-GC-59	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	5	1.6 µg	0.80 µg	CS <sub>2</sub>
Interferences		Comments		
		Turnaround is 10 business days.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Dimethyl Sulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-18-3	LM-GC-59	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	5	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		
		Turnaround is 10 business days.		

### Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-83-8	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.59 µg	0.30 µg	CS <sub>2</sub>
Interferences		Comments		

### Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-83-8	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.6	15-480	0.89 µg	0.45 µg	CS <sub>2</sub>
Interferences		Comments		

### Dimethylformamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
68-12-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.5	15-480	1.7 µg	0.83 µg	MC
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents.		

### Dimethylformamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
68-12-2	NIOSH 2004	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-1	15-80	1.5 µg	0.75 µg	MeOH
Interferences		Comments		
		Silica gel has a high affinity for water; high relative humidity may adversely affect the efficiency of analyte adsorption.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Diethyl Phthalate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
117-84-1	OSHA 104	GC-FID	OVS-Tenax (SKC 226-56)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.0	1-240	0.27 µg	0.14 µg	Tol
Interferences		Comments		

**Dioxane(p-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-91-1	NIOSH 1602	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	1-10	1.3 µg	0.65 µg	CS <sub>2</sub>
Interferences		Comments		

**Dioxane(p-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-91-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
34.5	15-480	2.0 µg	1.0 µg	CS <sub>2</sub>
Interferences		Comments		

**Diphenyl (Biphenyl)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
92-52-4	OSHA PV2022	GC-FID	XAD-7 (SKC 226-95)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.20	20	0.46 µg	0.23 µg	CS <sub>2</sub>
Interferences		Comments		

**Dipropylene Glycol Methyl Ether (DPGME)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
34590-94-8	NIOSH 2554	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	2-10	3.0 µg	1.5 µg	MC
Interferences		Comments		
		2020 NIC, TWA = 50ppm, STEL = 100ppm		

†(LPM) (CC/Min)

††(L) (Minutes)

### Dipropylene Glycol Methyl Ether (DPGME)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
34590-94-8	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
25.3	15-480	4.5 µg	2.3 µg	CS <sub>2</sub>
Interferences		Comments		
		2020 NIC, TWA = 50ppm, STEL = 100ppm		

### Dipropylene Glycol Methyl Ether Acetate (DPGMEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
88917-22-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
22.8	15-480	1.5 µg	0.75 µg	CS <sub>2</sub>
Interferences		Comments		

### Divinyl Benzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1321-74-0	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.57 µg	0.29 µg	CS <sub>2</sub>
Interferences		Comments		

### Divinyl Benzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1321-74-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
23.3	15-480	0.86 µg	0.43 µg	CS <sub>2</sub>
Interferences		Comments		

### Enflurane (Ethrane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
13838-16-9	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	1-10	2.6 µg	1.3 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Enflurane (Ethrane)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
13838-16-9	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.3	15-480	4.0 µg	2.0 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

**Epichlorohydrin (1-Chloro-2,3-epoxy propane)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-89-8	NIOSH 1010	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	2-30	1.5 µg	0.75 µg	CS <sub>2</sub>
Interferences		Comments		

**Epichlorohydrin (1-Chloro-2,3-epoxy propane)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-89-8	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
29.6	15-480	2.3 µg	1.2 µg	CS <sub>2</sub>
Interferences		Comments		

**Ethanolamine (2-Aminoethanol)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-43-5	NIOSH 2007	IC	ORBO 53 or SGT** (SUPELCO 20265) (SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.5	5	3.0 µg	1.5 µg	EA
Interferences		Comments		
		Store in freezer after sampling. Ship cold.		

**Ethoxyethanol(2-) (Cellosolve)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-80-5	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1-6	0.54 µg	0.27 µg	5%MeOH/MC
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Ethoxyethanol(2-) (Cellosolve)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-80-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.4	15-480	0.81 µg	0.41 µg	MC
Interferences		Comments		
		Sample separately from other solvents.		

**Ethoxyethyl(2-) Acetate (Cellosolve acetate)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-15-9	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	2.8 µg	1.4 µg	CS <sub>2</sub>
Interferences		Comments		

**Ethoxyethyl(2-) Acetate (Cellosolve acetate)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-15-9	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
26.6	15-480	4.2 µg	2.1 µg	CS <sub>2</sub>
Interferences		Comments		

**Ethyl 2-cyanoacrylate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7085-85-0	OSHA 55	HPLC	XAD-7, Acid (SKC 226-98)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	12	0.70 µg	0.35 µg	
Interferences		Comments		
		Ship and store cold.		

**Ethyl 3-ethoxypropionate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
763-69-9	OSHA PV2025	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	1-10	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

### Ethyl 3-ethoxypropionate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
763-69-9	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
25.8	15-480	1.7 µg	0.85 µg	CS <sub>2</sub>
Interferences		Comments		

### Ethyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-78-6	NIOSH 1457	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	1.2 µg	0.60 µg	CS <sub>2</sub>
Interferences		Comments		

### Ethyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-78-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
34.5	15-360	1.8 µg	0.90 µg	CS <sub>2</sub>
Interferences		Comments		

### Ethyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
140-88-5	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.98 µg	0.49 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm. Store and ship cold overnight.		

### Ethyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
140-88-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.2	15-480	1.5 µg	0.75 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)



**Ethyl Alcohol (Ethanol)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-17-5	NIOSH 1400	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.1-1	1.2 µg	0.60 µg	BUT/CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

**Ethyl Alcohol (Ethanol)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-17-5	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
43.7	15-120	1.8 µg	0.90 µg	ACN
Interferences		Comments		
		Use 3M 3520. Separate the front section of the monitor from the back section and cap immediately after sampling. Store and ship cold overnight. Sample separately from other solvents.		

**Ethyl Benzene**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-41-4	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-30	0.5 µg	0.25 µg	CS <sub>2</sub>
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%.		

**Ethyl Benzene**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-41-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.3	15-480	0.75 µg	0.36 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Ethyl Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
60-29-7	NIOSH 1610	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.75 µg	0.38 µg	Ethyl Acetate
Interferences		Comments		
Hexane may co-elute with the analyte of interest.		High humidity may greatly decrease breakthrough volume. Store and ship cold.		

## Ethyl Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
60-29-7	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
36.8	15-240	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate the front section of the monitor from the back section and cap immediately after sampling. Store and ship cold.		

## Ethyl Lactate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
687-47-8	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.85 µg	0.43 µg	CS <sub>2</sub>
Interferences		Comments		

## Ethyl Methacrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
97-63-2	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.65 µg	0.33 µg	CS <sub>2</sub>
Interferences		Comments		
		Ship and store cold.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Ethylamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-04-7	NIOSH 2010	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-1	5-30	NA	NA	Amine1
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

## Ethylene Chlorohydrin (2-Chloroethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-07-3	NIOSH 2513	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.51 µg	0.26 µg	5%IPA/CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at flow rate of 0.2 lpm. High humidity may greatly decrease the breakthrough volume.		

## Ethylene Chlorohydrin (2-Chloroethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-07-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
33.9	15-480	0.77 µg	0.39 µg	MC CS <sub>2</sub>
Interferences		Comments		

## Ethylene Dichloride (1,2-Dichloroethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-06-2	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-50	1.5 µg	0.74 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Ethylene Dichloride (1,2-Dichloroethane)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-06-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
33.2	15 -480	2.3 µg	1.1 µg	CS <sub>2</sub>
Interferences		Comments		

**Ethylene Glycol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-21-1	NIOSH 5523	GC-FID	OVS 7 (SKC 226-57)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.5-2	5-60	2.2 µg	1.1 µg	MeOH
Interferences		Comments		

**Ethylene Oxide**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-21-8	ASTM D5578-04	GC-FID	ORBO 78 (SUPELCO 20355)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.15	1-24	0.64 µg	0.32 µg	ACN/TOL
Interferences		Comments		
		Store and ship cold. Sample separately from CS <sub>2</sub> compatible solvents.		

**Ethylene Oxide**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-21-8	3M Method	GC-FID	OVM (3M 3551)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
49.3	15 -480	0.96 µg	0.48 µg	ACN/TOL
Interferences		Comments		
		Store and ship cold. Sample separately from CS <sub>2</sub> compatible solvents.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Ethylenediamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-15-3	OSHA 60	HPLC	XAD-2, NITC (SKC 226-30-18)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	10	0.080 µg	0.040 µg	Amine2
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

## Flour Dust

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	HSE MDHS-14	GRAV	PVC, IOM	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	960	100 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		Use IOM sampler with pre-weighed PVC. Contact Lab one week before intended use. The availability of IOM samplers is limited. Rental charge for the IOM samplers applies.		

## Fluoranthene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
206-44-0	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.30 µg	0.15 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

## Fluorene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
86-73-7	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.32 µg	0.16µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Fluorides, Particulate/Hydrogen Fluoride

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	NIOSH 7906	IC	Cellulose Nitrate,Na <sub>2</sub> CO <sub>3</sub> (SKC 225-9031)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-2	70-1000	13 µg	6.5 µg	Acid1	
Interferences		Comments			
Recovery of gaseous HF is reduced at high humidity.		Treated filter is stable for 14 days. Order media one week ahead of survey date. Media are prepared when ordered. Store and ship cold. <b>Specialty filter. Media charge applies.</b>			

## Forane (Isoflurane)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
26675-46-7	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.05	12	4.9 µg	2.5 µg	CS <sub>2</sub>	
Interferences		Comments			
		Store and ship cold overnight.			

## Forane (Isoflurane)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
26675-46-7	3M Method	GC-FID	OVM (3M 3500)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
28.3	15-480	7.1 µg	3.6 µg	CS <sub>2</sub>	
Interferences		Comments			
		Store and ship cold overnight.			

## Formaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media		
50-00-0	NIOSH 2016	HPLC	AT Monitor (N571AT)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
16.2	15-480	0.018 µg	0.0090 µg	Aldehyde	
Interferences		Comments			
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Formaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-00-0	NIOSH 2016	HPLC	Sep-Pak (WAT047205)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	10-100	0.090 µg	0.045 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use a flow rate of 1.0 lpm for STEL.		

## Formaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-00-0	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	1-15	0.036 µg	0.018 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use a flow rate of 1.0 lpm for STEL. For IAQ and LEED sampling where TVOC and formaldehyde are collected, do not sample with Charcoal Tubes and DNPH tubes in tandem. The sorbent in the DNPH tubes may off-gas acetonitrile which can cause a positive interference in the TVOC results.		

## Formamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-12-7	NIOSH 2004	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-1	15-80	3.1 µg	1.6 µg	MeOH
Interferences		Comments		
		Silica gel has a high affinity for water; high relative humidity may adversely affect the efficiency of analyte adsorption. 2019 NIC, TWA = 1ppm, Skin; A3 adopted in 2020.		

## Formic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-18-6	NIOSH 2011	IC	PTFE3-SGT** (SKC 225-17A, SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.5	5-100	1.5 µg	0.75 µg	Acid2
Interferences		Comments		
		Do not sample with inorganic acids.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Furfural

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-01-1	OSHA 72	GC-FID	SKC 226-81A	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1	180	0.96 µg	0.48 µg	1%DMF/CS <sub>2</sub>
Interferences		Comments		
Furfuryl Alcohol is a sampling interference.				

## Furfuryl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-00-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.6	15-480	0.96 µg	0.48 µg	MC
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents.		

## Gasoline

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8006-61-9	NIOSH 1550	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1.3-20	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Please ship bulk sample separately from the air samples. <b>Preferred for STEL sampling.</b> Use a flow rate of 0.2 lpm for STEL.		

## Gasoline

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8006-61-9	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.5	15-480	1.7 µg	0.85 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Please ship bulk sample separately from the air samples.		

†(LPM) (CC/Min)

††(L) (Minutes)



## Germanium

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-56-4	NIOSH 7300	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		100-1000	0.38 µg	0.19 µg	Metals
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

## Glutaraldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media		
111-30-8	OSHA 64	HPLC	AT Monitor (N571AT)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
7.24		15-480	0.016 µg	0.0080 µg	
Interferences		Comments			
		Refrigerate media before and after sampling. Ship cold overnight.			

## Glutaraldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media		
111-30-8	NIOSH 2016	HPLC	Sep-Pak (WAT047205)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5		10-100	0.080 µg	0.040 µg	
Interferences		Comments			
		Refrigerate media before and after sampling. Ship cold overnight. Preferred for STEL sampling.			

## Glutaraldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media		
111-30-8	NIOSH 2532	HPLC	SGT, DNPH (SKC 226-119)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.5		1-30	0.032 µg	0.016 µg	
Interferences		Comments			
		Refrigerate media before and after sampling. Ship cold overnight. Preferred for STEL sampling.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Gold

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-57-5	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	50-1000	0.28 µg	0.14 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

## Grain Dust

GRAV Dust					
CAS #	Analytical Method	Analytical Technique	Sampling Media		
	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-15	50-7200	50 µg	10 µg		
Interferences		Comments			
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.			

## Graphite

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7782-42-5	NIOSH 0600	GRAV	PVC (SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1.7	100-816	50 µg	10 µg		
Interferences		Comments			
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.			

## Halothane (Fluothane)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
151-67-7	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.05	12	3.4 µg	1.7 µg	CS <sub>2</sub>	
Interferences		Comments			
		Store and ship cold overnight.			

†(LPM) (CC/Min)

††(L) (Minutes)

### Halothane (Fluothane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
151-67-7	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.2	15-480	5.1 µg	2.5 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

### Heptane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
142-82-5	NIOSH 1500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.45 µg	0.23 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at 0.2 lpm.		

### Heptane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
142-82-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.9	15-480	0.68 µg	0.34 µg	CS <sub>2</sub>
Interferences		Comments		

### Heptanone(2-) (Methyl Amyl Ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-43-0	NIOSH 1301	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-25	0.64 µg	0.32 µg	MeOH/CS <sub>2</sub>
Interferences		Comments		

### Heptanone(2-) (Methyl Amyl Ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-43-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.9	15-480	0.96 µg	0.48 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

### Hexamethylene Diisocyanate (1,6-) Homopolymer (HDI Homo)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
28182-81-2	OSHA PV2125	HPLC	GFF, 1-2PP	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	15	0.30 µg	0.15 µg	Isocyanate
Interferences		Comments		
Any compound having the same retention time as the analyte is a possible interference. However, chromatographic conditions can be altered to separate an interference.		Sample open-faced. Media must be stored cold before and after sampling. Filter is stable for 6 months if kept cold. Ship cold overnight.		

### Hexamethylene Diisocyanate (1,6-) Homopolymer (HDI Homo)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
28182-81-2	OSHA PV2125	HPLC	GFF Wipes	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
N/A	N/A	0.30 µg	0.15 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 1 month if kept cold.		

### Hexamethylene Diisocyanate (1,6-) (HDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
822-06-0	OSHA 42	HPLC	GFF, 1-2PP	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	15-240	0.020 µg	0.010 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Sample open-faced. Media must be stored cold before and after sampling. Filter is stable for 6 months if kept cold. Ship cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Hexamethylene Diisocyanate (1,6-) (HDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
822-06-0	OSHA 42	HPLC	GFF Wipes	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
N/A	N/A	0.020 µg	0.010 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 1 month if kept cold.		

### Hexane(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-54-3	NIOSH 1500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.45 µg	0.23 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at 0.2 lpm.		

### Hexane(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-54-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.0	15-480	0.72 µg	0.36 µg	CS <sub>2</sub>
Interferences		Comments		

### Hexyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2499-95-8	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.90 µg	0.45 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Hexylene Glycol (2-Methyl-2,4-pentanediol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-41-5	NIOSH 5523	GC-FID	OVS 7 (SKC 226-57)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.5-2	5-60	2.5 µg	1.3 µg	MeOH
Interferences		Comments		

### Hydrazine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
302-01-2	OSHA 108	HPLC	GFF, Acid	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	240	0.025 µg	0.013 µg	
Interferences		Comments		
		Media have short shelf-life. Media are prepared on request. Please contact the Lab 5 days before survey.		

### Hydrogen Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10035-10-6	NIOSH 7903	IC	MCE2, SGT** (SKC 225-19, SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2-0.5	10-100	4.6 µg	2.3 µg	Acid1
Interferences		Comments		
Particulate salts of the acid will give a positive interference.		Use a flow rate of 0.5 lpm for STEL sampling.		

### Hydrogen Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10035-10-6	NIOSH 7907	IC	SKC 225-9032	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	35-600	3.8 µg	1.9 µg	
Interferences		Comments		
Inorganic acids can react with co-sampled particulate matter on the pre-filter, leading to low results.		Order media one week ahead of survey. Ship and store cold. <b>Specialty Filter. Media charge applies.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)

## Hydrogen Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7647-01-0	NIOSH 7903	IC	MCE2, SGT** (SKC 225-19, SKC 226-10-03)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.2-0.5	10-100	2.3 µg	1.2 µg	Acid1	
Interferences		Comments			
Particulate salts of the acid will give a positive interference.		Use a flow rate of 0.5 lpm for STEL sampling.			

## Hydrogen Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7647-01-0	NIOSH 7907	IC	SKC 225-9032		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
2	40-600	13 µg	6.4 µg		
Interferences		Comments			
Inorganic acids can react with co-sampled particulate matter on the pre-filter, leading to low results. Potentially interfering particulate chlorides and nitrates removed by the pre-filter can react with the sampled acids and liberate HCl and HNO <sub>3</sub> , which gets collected on the sampling filter, leading to high results.		Order media one week ahead of survey. Ship and store cold. <b>Specialty Filter. Media charge applies.</b>			

## Hydrogen Cyanide

CAS #	Analytical Method	Analytical Technique	Sampling Media		
74-90-8	NIOSH 7904	ISE	Soda Lime Tube (SKC 226-210)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.2	10-90	5.2 µg	2.6 µg		
Interferences		Comments			
Sulfide, chloride, iodide, bromide, cadmium, zinc, silver, nickel, cuprous iron and mercury interfere.		Use a flow rate of 0.2 lpm for STEL sampling. Method is not covered under our AIHA-LAP, LLC scope of accreditation.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Cyanide Salts as CN

CAS #	Analytical Method	Analytical Technique	Sampling Media	
74-90-8	NIOSH 7904	ISE	PVC	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	15-100	5.2 µg	2.6 µg	
Interferences		Comments		
Sulfide, chloride, iodide, bromide, cadmium, zinc, silver, nickel, cuprous iron and mercury interfere.		Preferred for STEL sampling. Method is not covered under our AIHA-LAP, LLC scope of accreditation.		

## Hydrogen Fluoride, as F

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7664-39-3	NIOSH 7903	IC	MCE2, SGT** (SKC 225-19, SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2-0.5	30-100	3.0 µg	1.5 µg	Acid1
Interferences		Comments		
Particulate salts of the acid will give a positive interference.		Use the maximum flow rate for STEL sampling.		

## Hydrogen Peroxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7722-84-1	OSHA ID-1019	UV/VIS	QFF, titanium oxysulfate (SKC 225-9030)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.0	240	10.0 µg	5.0 µg	
Interferences		Comments		
Any compound with a response, or that reacts with the titanium oxysulfate to produce a response, at 410nm is a potential interferent.		After sampling seal the cassette with the end plugs and wrap each cassette with tin foil. Order tin foil for wrapping the samples. Filters have limited shelf-life. This method is not covered under our AIHA-LAP, LLC scope of accreditation. <b>Specialty Filter. Media charge applies.</b>		

## Hydrogen Sulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7783-06-4	NIOSH 6013	IC	ORBO 34 (SUPELCO 20211)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	30-70	8.9 µg	4.4 µg	
Interferences		Comments		
Sulfur dioxide gas may give a positive interference for hydrogen sulfide.		Use a flow rate of 0.5 lpm for STEL.		

†(LPM) (CC/Min)

††(L) (Minutes)



### Hydroquinone (Dihydroxybenzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-31-9	NIOSH 5004	HPLC	MCE (SKC 225-5)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	30-180	10 µg	5 µg	
Interferences		Comments		
		Hydroquinone is unstable on the collecting media. Stabilize immediately after collecting by transferring filter into a vial containing 1% acetic acid.		

### Hydroquinone (Dihydroxybenzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-31-9	OSHA PV2094	HPLC	XAD-7, Acid (SKC 226-98)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2	20	0.30 µg	0.15 µg	
Interferences		Comments		
		Preferred method.		

### Indium and Compounds as In

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-74-6	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	150-1000	1.5 µg	0.75 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

### Inorganic Acid Scan

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 7903	IC	MCE2, SGT** (SKC 225-19, SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2-0.5	50-100	N/A	N/A	
Interferences		Comments		
		See List of Scans for list of individual inorganic acid/anion.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Iodine and Iodides as I

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7553-56-2	NIOSH 6005	IC	CT, KOH (SKC 226-67)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.5-1.0	50-200	5.4 µg	2.7 µg	
Interferences		Comments		
Particulate iodide salts, hydrogen iodide or organic iodides may give a positive interference.		Use a flow rate of 1.0 lpm for STEL.		

## Iron Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1309-37-1	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	4.3 µg	2.2 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of iron are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Iron Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1309-37-1	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	3.6 µg	1.8 µg	Metals
Interferences		Comments		
		All forms of iron are quantified. As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. <b>Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.</b>		

## Iron

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1309-37-1	OSHA ID-121	ICP	Ghost wipe (225-2414)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
NA	NA	30 µg	15 µg	Metals 2
Interferences		Comments		
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Isobutyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-19-0	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.81 µg	0.41 µg	CS <sub>2</sub>
Interferences		Comments		

### Isobutyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-19-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.0	15-480	1.2 µg	0.60 µg	CS <sub>2</sub>
Interferences		Comments		

### Isobutyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-83-1	NIOSH 1401	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.59 µg	0.30 µg	1%PRO/CS <sub>2</sub>
Interferences		Comments		

### Isobutyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-83-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.9	15-480	0.89 µg	0.45 µg	MC CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Isocyanate Scan

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	OSHA 42	HPLC	GFF, 1-2PP	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		See List of Scans for individual isocyanates. Sample open-faced. Media must be stored cold before and after sampling. Filter is stable for 6 months if kept cold. Ship cold overnight.		

## Isoflurane (Forane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
26675-46-7	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05	12	4.9 µg	2.5 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight. 2019 NIC, TWA = 5 ppm, A4 was retained in 2020.		

## Isoflurane (Forane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
26675-46-7	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.3	15-480	7.4 µg	3.2 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight. 2019 NIC, TWA = 5 ppm, A4 was retained in 2020.		

## Isooctane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-84-1	NIOSH 1500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.52 µg	0.26 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Isooctane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-84-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.1	15-480	0.78 µg	0.39 µg	CS <sub>2</sub>
Interferences		Comments		

## Isophorone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-59-1	NIOSH 2508	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.90 µg	0.45 µg	CS <sub>2</sub>
Interferences		Comments		
		Use a flow rate of 0.2 lpm for STEL. High humidity may greatly decrease the breakthrough volume.		

## Isophorone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-59-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
21.7	15-480	1.4 µg	0.70 µg	CS <sub>2</sub>
Interferences		Comments		

## Isophorone Diisocyanate (IPDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
4098-71-9	OSHA 42	HPLC	GFF, 1-2PP	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	15-240	0.034 µg	0.017 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Sample open-faced. Media must be stored cold before and after sampling. Filter is stable for 6 months if kept cold. Ship cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Isophorone Diisocyanate (IPDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
4098-71-9	OSHA 42	HPLC	GFF Wipes	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
N/A	N/A	0.034 µg	0.017 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 1 month if kept cold.		

### Isopropyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-21-4	NIOSH 1454	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.02-0.2	0.1-9	1.0 µg	0.50 µg	CS <sub>2</sub>
Interferences		Comments		
		Use a flow rate of 0.2 lpm for STEL.		

### Isopropyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-21-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.7	15-480	1.5 µg	0.75 µg	CS <sub>2</sub>
Interferences		Comments		

### Isopropyl Alcohol (Isopropanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-63-0	NIOSH 1400	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.3-3	1.0 µg	0.51 µg	BUT/CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Isopropyl Alcohol (Isopropanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-63-0	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
39.4	15-240	1.5 µg	0.77 µg	ACN CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate the front section of the monitor from the back section and cap immediately after sampling. Ship cold.		

## Kaolin

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1332-58-7	NIOSH 0600	GRAV	PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.7	100-816	50 µg	10 µg	
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

## Kerosene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8008-20-6	NIOSH 1550	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1.3-20	2.9 µg	1.5 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample; ship it separately from the air samples.		

## Kerosene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8008-20-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.3	15-480	4.4 µg	2.2 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample; ship it separately from the air samples.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Lactic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-21-5	NIOSH 2011	IC	SGT**(SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2-0.5	5-100	3.0 µg	1.2 µg	Acid2
Interferences		Comments		
		Do not sample for inorganic acids using the same tube.		

## Lanthanum

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-91-0	NIOSH 7301	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	120-1000	0.10 µg	0.050 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Lead and Inorganic Compounds as Pb

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-92-1	OSHA ID-121	ICP	ghost wipe (SKC 225-2414)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
NA	NA	1.0 µg	0.50 µg	Metals2
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Lead and Inorganic Compounds as Pb

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-92-1	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	120-1000	0.25 µg	0.13 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)



## Lead and Inorganic Compounds as Pb

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7439-92-1	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		40-240	0.032 µg	0.016 µg	Metals
Interferences		Comments			
		As part of the Lab’s QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

## Lead and Inorganic Compounds as Pb

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7439-92-1	OSHA ID-121 NIOSH 7301	ICP	Paint chips		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
NA		NA	0.25 µg	0.13 µg	Metals2
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Needs at least 1 gram of bulk sample. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.			

## Lead Chromate as Cr(VI)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7758-97-6	OSHA ID-215	IC	PVC (SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
2	500-960	0.05 µg	0.02 µg		
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Please indicate in sample submission form the operation sampled (e.g., spray paint, chrome plating, welding, etc.) as the extraction method is different for spray paint samples. Refrigerate samples and ship overnight as soon as possible. 500 liters is the minimum air volume at 50% TLV. Samples from plating operations must be analyzed within 6 days from date of sampling. Samples from welding operations must be analyzed within 8 days from date of sampling. Hexavalent chromium compounds, including Chromite ore processing as Chrome (VI).			

†(LPM) (CC/Min)

††(L) (Minutes)

**Limonene(d-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
138-86-3	NIOSH 1552	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.52 µg	0.26 µg	CS <sub>2</sub>
Interferences		Comments		

**Limonene(d-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
138-86-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.7	15-480	0.78 µg	0.39 µg	CS <sub>2</sub>
Interferences		Comments		

**Lithium Salts**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-93-1	NIOSH 7301	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	250-1000	0.10 µg	0.050µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of lithium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

**Magnesium**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-95-4	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	1.0 µg	0.50 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of magnesium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Magnesium Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1309-48-4	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	1.7 µg	0.85 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of magnesium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

## Magnesium Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1309-48-4	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	40-240	2.5 µg	1.3 µg	Metals
Interferences		Comments		
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.		

## Magnesium Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1309-48-4	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	40-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

## Maleic Anhydride

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-31-6	OSHA 86	HPLC	GFF, Vamine	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.5	60	0.74 µg	0.37 µg	ACN/DMSO
Interferences		Comments		
Both phthalic and trimellitic anhydride should be considered as potential sampling interferences.		Sampling media has short shelf-life, so media is prepared when ordered. Please order filters at least 5 days prior to survey date.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Maleic Anhydride

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-31-6	OSHA 25	HPLC	XAD-2,p-An/XAD-2 (SKC 226-30-07, SKC 226-30)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	20	0.030 µg	0.015 µg	
Interferences		Comments		
Both phthalic and trimellitic anhydride should be considered as potential sampling interferences.		Sample with XAD-2 treated tube and XAD-2 untreated tube in series. Separate and cap tubes after sampling. Media has short shelf-life. Please contact the Lab before survey date for more info. Preferred method.		

## Manganese, Elemental and Inorganic compounds as Mn

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-96-5	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	125-1000	0.15 µg	0.075 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. Manganese may be sampled using respirable PPI. Use a flow rate of 2 lpm. Please indicate in your request if you want the filter inside the PPI pre-weighed. Please contact lab a week before your survey if using respirable PPI. Media charge applies for PPIs.		

## Manganese, Elemental and Inorganic compounds as Mn

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-96-5	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	75-240	0.15 µg	0.075µg	Metals
Interferences		Comments		
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples. Manganese may be sampled using respirable PPI. Use a flow rate of 2 lpm. Please indicate in your request if you want the filter inside the PPI pre-weighed. Please contact lab a week before your survey if using respirable PPI. Media charge applies for PPIs.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Manganese, elemental and Inorganic compounds as Mn

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-96-5	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe (SKC 225-2414)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
NA	NA	0.60 µg	0.30 µg	Metals2
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Medical Gases

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0500 LMI-A5	GRAV, GC-FID GC-ECD, GC-TCD GC-XSD, GC-DID	PTFE4 (PALL TF-450) Cylinder	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
Note 2* see Comments	1000 (gravimetric) 300cc/25 psig (Cylinder)	Note 2* see Comments	Note 2* see Comments	
Interferences		Comments		
At high levels argon interferes with oxygen and carbon dioxide interferes with nitrous oxide.		*See Note 2 in the Comments section of "IH Lab Sampling Guide Analyte Descriptions and Abbreviations"		

## Mercury as Hg (Elemental and inorganic forms)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-97-6	NIOSH 6009	AA-CV	Carulite/ Anasorb C300 (SKC 226-17-1A/3A)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.15-0.25	50-120	0.12 µg	0.060 µg	
Interferences		Comments		
Inorganic and organic mercury compounds may cause a positive interference. Oxidizing gases, including chlorine, do not interfere.		A 37-mm, cellulose ester membrane filter in a cassette preceding the sorbent may be used if particulate mercury is to be determined separately.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Mercury as Hg (Elemental and inorganic forms)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-97-6	OSHA ID-140	AA-CV	PS (SKC 520-02A/03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.020	9.6	0.12 µg	0.060 µg	
Interferences		Comments		
Particulate mercury compounds are a positive interference.		Call Lab one week before sampling. Lab has limited number of passive mercury monitor holders. Refer to OSHA ID-145 if sampling in workplaces containing both mercury vapor and particulate. Specialty media. Media charge applies.		

## Mercury as Hg Particulate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-97-6	OSHA ID-145	AA-CV	Air MCE (SKC225-5) Wipe (Whatman#42) Bulk	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	50	0.12 µg	0.060 µg	
Interferences		Comments		
		If mercury vapor is suspected to be present, please refer to OSHA ID-140 or NIOSH 6009 for additional sampling information.		

## Mesityl Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-79-7	NIOSH 1301	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2	25	1.3 µg	0.65 µg	CS <sub>2</sub>
Interferences		Comments		

## Metalworking Fluids

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 5524	GRAV	Pre-weighed PTFE1 (Zefon FPTFE137)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	768-960	75 µg	15 µg	
Interferences		Comments		
The method is non-specific and measures all substances extractable by organic solvents.		NIOSH 5524 recommends submitting one bulk sample of each type of fluid for solubility testing. Refrigerate samples if unable to ship immediately. Ship cold. Please note that samples must be analyzed within 2 weeks after collection.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Methanol (Methyl alcohol)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-56-1	OSHA 5001	GC-FID	Anasorb747/Anasorb747 (SKC 226-82)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05	1-5	2.3 µg	1.1 µg	DMF/CS <sub>2</sub>
Interferences		Comments		
		Preferred method. Sample with 2 Anasorb 747 tubes in series ("Part A" as the front section and "Part B" as the back section. Please order as a set.) Recommended air volume is 5 liters when relative humidity is >50% and 3 liters when relative humidity is <50% at 25°C. Separate and cap tubes after sampling. Store and ship cold overnight.		

**Methanol (Methyl alcohol)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-56-1	NIOSH 2000	GC-FID	SGT-SGT (SKC 226-51)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.02-0.2	1-5	2.5 µg	1.3 µg	5%IPA
Interferences		Comments		
		Sample using 2 SGT tubes in series. Separate and cap tubes after sampling. Use a flow rate of 0.2 lpm for STEL. Do not use in area where humidity is high. Store and ship cold overnight.		

**Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-98-2	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1-10	1.3 µg	0.65 µg	MeOH/MC
Interferences		Comments		

**Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-98-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.4	15-480	2.0 µg	1.0 µg	MC CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Methoxyethanol(2-) (Methyl cellosolve, EGME)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-86-4	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1-50	1.7 µg	0.85 µg	%MeOH/MC
Interferences		Comments		
		Preferred method for STEL sampling. Sample at 0.05 lpm.		

**Methoxyethanol(2-) (Methyl cellosolve, EGME)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-86-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
36.3	120-480	2.6 µg	1.3 µg	MC CS <sub>2</sub>
Interferences		Comments		

**Methoxyethoxy(2-(2-)) Ethanol (Diethylene glycol methyl ether)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-77-3	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1-10	1.3 µg	0.65 µg	5%MeOH/MC
Interferences		Comments		

**Methoxyethoxy(2-(2-)) Ethanol (Diethylene glycol methyl ether)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-77-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.8	15-480	2.0 µg	1.0 µg	MC CS <sub>2</sub>
Interferences		Comments		

**Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-49-6	NIOSH 1451	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.2-20	0.46 µg	0.23 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred method for STEL sampling.		

†(LPM) (CC/Min)

††(L) (Minutes)



**Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-49-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
29.0	120-480	0.67 µg	0.34 µg	CS <sub>2</sub>
Interferences		Comments		

**Methoxyethyl(2-) Ether (Diethylene glycol dimethyl ether)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-96-6	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1-10	0.58 µg	0.29 µg	5%MeOH/MC
Interferences		Comments		

**Methoxyethyl(2-) Ether (Diethylene glycol dimethyl ether)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-96-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
26.7	15-480	0.87 µg	0.44 µg	MC CS <sub>2</sub>
Interferences		Comments		

**Methyl Acetate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-20-9	NIOSH 1458	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.2-10	1.5 µg	0.75 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred method for STEL sampling. Sample at 0.2 lpm.		

**Methyl Acetate**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-20-9	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
37.0	15-480	2.3 µg	1.2 µg	CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate the front section of the monitor from the back section and cap immediately after sampling. Ship cold.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Methyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-33-3	NIOSH 1459	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-5	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		

## Methyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-33-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.8	15-480	1.7 µg	0.85 µg	CS <sub>2</sub>
Interferences		Comments		

## Methyl Alcohol (Methanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-56-1	OSHA 5001	GC-FID	Anasorb747/Anasorb747 (SKC 226-82)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05	1-5	2.3 µg	1.1 µg	DMF/CS <sub>2</sub>
Interferences		Comments		
		Preferred method. Sample Anasorb 747 tubes in series (with "Part A" in front and "PART B" in back section. Recommended air volume is 5 liters when relative humidity is >50% and 3 liters when relative humidity is <50% at 25° C. Separate and cap tubes after sampling. Store and ship cold overnight.		

## Methyl Alcohol (Methanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-56-1	NIOSH 2000	GC-FID	SGT-SGT (SKC 226-51)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.02-0.2	1-5	2.3 µg	1.7 µg	5%IPA
Interferences		Comments		
		Sample using 2 SGT tubes in series. Use a flow rate of 0.2 lpm for STEL. Separate and cap tubes immediately after sampling. Do not use in area where humidity is high. Store and ship cold overnight		

†(LPM) (CC/Min)

††(L) (Minutes)

**Methyl Amyl Ketone (2-Heptanone)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-43-0	NIOSH 1301	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-25	0.61 µg	0.31 µg	MeOH/CS <sub>2</sub>
Interferences		Comments		
		Use a flow rate of 0.2 lpm for STEL.		

**Methyl Amyl Ketone (2-Heptanone)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-43-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.9	15-480	0.92 µg	0.46 µg	CS <sub>2</sub>
Interferences		Comments		

**Methyl Aniline**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-61-8	NIOSH 2002	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.02-0.5	5-30	NA	NA	Amine3
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

**Methyl Chloroform (1,1,1-Trichloroethane)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-55-6	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	2.2 µg	1.1 µg	CS <sub>2</sub>
Interferences		Comments		

**Methyl Chloroform (1,1,1-Trichloroethane)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-55-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.9	15-480	3.3 µg	1.7 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

### Methyl Cyclopentane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-37-7	NIOSH 1500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.38 µg	0.19 µg	CS <sub>2</sub>
Interferences		Comments		

### Methyl Cyclopentane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-37-7	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.5	15-480	0.57 µg	0.29 µg	CS <sub>2</sub>
Interferences		Comments		

### Methyl Ethyl Ketone (2-Butanone, MEK)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-93-3	NIOSH 2500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.25-12	0.86 µg	0.43 µg	CS <sub>2</sub>
Interferences		Comments		
Isopropyl acetate may co-elute with MEK.		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm. Ship and store cold.		

### Methyl Ethyl Ketone (2-Butanone, MEK)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-93-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
36.3	15-480	1.3 µg	0.65 µg	CS <sub>2</sub>
Interferences		Comments		
Isopropyl acetate may co-elute with MEK.		Ship and store cold.		

### Methyl Isoamyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-12-3	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.66 µg	0.33 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

### Methyl Isoamyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-12-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.0	15-480	0.99 µg	0.50 µg	CS <sub>2</sub>
Interferences		Comments		

### Methyl Isobutyl Ketone (MIBK)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-10-1	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.54 µg	0.27 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at 0.2 lpm.		

### Methyl Isobutyl Ketone (MIBK)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-10-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.0	15-480	0.80 µg	0.40 µg	CS <sub>2</sub>
Interferences		Comments		

### Methyl Isopropyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
563-80-4	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	1.0 µg	0.50 µg	CS <sub>2</sub>
Interferences		Comments		
		Ship and store cold.		

### Methyl Isopropyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
563-80-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.8	15-480	1.5 µg	0.75 µg	CS <sub>2</sub>
Interferences		Comments		
		Ship and store cold.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Methyl Methacrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
80-62-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.8	15-480	1.5 µg	0.75 µg	CS <sub>2</sub>
Interferences		Comments		
		Ship cold.		

### Methyl Propyl Ketone (2-Pentanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-87-9	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.73 µg	0.37 µg	CS <sub>2</sub>
Interferences		Comments		

### Methyl Propyl Ketone (2-Pentanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-87-9	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
33.0	15-480	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		

### Methyl Styrene(a-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-83-9	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-30	0.41 µg	0.21 µg	CS <sub>2</sub>
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%. Store and ship cold overnight.		

### Methyl Styrene(a-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-83-9	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
25.0	15-480	0.62 µg	0.31 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Methyl Tert-butyl Ether (MTBE)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1634-04-4	NIOSH 1615	GC-FID	CT-CT (SKC 226-01)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-0.2	2-96	0.67 µg	0.34 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at 0.2 lpm. Store and ship cold immediately.		

### Methyl Tert-butyl ether (MTBE)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1634-04-4	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.8	15-480	1.0 µg	0.50 µg	CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate the front section of the monitor from the back section and cap immediately after sampling. Store and ship cold immediately.		

### Methyl Vinyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-94-4	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.27 µg	0.14 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

### Methyl(1-)-2-pyrrolidinone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
872-50-4	NIOSH 1302	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	0.5-125	1.1 µg	0.55 µg	5%MeOH/MC
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Methyl(1-)-2-pyrrolidinone**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
872-50-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.8	15-480	1.7 µg	0.85 µg	MC CS <sub>2</sub>
Interferences		Comments		

**Methylacrylonitrile**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
126-98-7	NIOSH 1604	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	3.5-20	1.0 µg	0.50 µg	AC/CS <sub>2</sub>
Interferences		Comments		

**Methylcyclohexane**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-87-2	NIOSH 1500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	4	0.45 µg	0.23 µg	CS <sub>2</sub>
Interferences		Comments		

**Methylcyclohexane**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-87-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.9	15-420	0.68 µg	0.34 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)



### Methylene Bis(4-cyclohexylisocyanate) (HMDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
5124-30-1	OSHA 47	HPLC	GFF, 1-2PP	
Sampling Rate†	Sampling Volume††	LOQ	LOQ	Compatibility Code
1	15-240	0.070 µg	0.035 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months if kept cold.		

### Methylene Bisphenyl Isocyanate (MDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
101-68-8	OSHA 47	HPLC	GFF, 1-2PP	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	15-240	0.014 µg	0.0070 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months if kept cold.		

### Methylene Bisphenyl Isocyanate (MDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
101-68-8	OSHA 47	HPLC	GFF Wipes	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
N/A	N/A	0.014 µg	0.0070 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 1 month if kept cold.		

### Methylene Chloride (Dichloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-09-2	NIOSH 1005	GC-FID	CT-CT (SKC 226-01)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	0.5-2.5	2.8 µg	1.4 µg	CS <sub>2</sub>
Interferences		Comments		
		Sample using two CT in series. Separate and cap tubes after sampling. Ship and store cold immediately.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Methylene Chloride (Dichloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-09-2	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
37.9	15-240	4.2 µg	2.1 µg	CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate the front section of the monitor from the back section and cap immediately after sampling. Ship and store cold immediately.		

### Methylene(4,4'-) Dianiline (MDA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
101-77-9	NIOSH 5029	HPLC	GFF, Acid	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-2	10-1000	0.68 µg	0.34 µg	
Interferences		Comments		
4,4'-Diphenyl methane diisocyanate (MDI) will interfere.		Media has short shelf life so it is not kept in stock. Please contact the lab one week before sampling to order media. Within 4 hours of sampling, transfer the filter to a glass vial containing 4 ml 0.1 N methanolic KOH. This method can also be used for wipe sampling.		

### Methylene(4,4')-bis(2 chloroaniline) (MOCA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
101-14-4	NIOSH P&CAM 236	HPLC	GFF-SGT (SKC 225-16, SKC226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2-1	50	0.50 µg	0.25 µg	
Interferences		Comments		

### Methylnaphthalene(2-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
91-57-6	NIOSH 5515	GC-MS	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.66 µg	0.33 µg	PNAs
Interferences		Comments		
		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Mica

CAS #	Analytical Method	Analytical Technique	Sampling Media	
12001-26-2	NIOSH 0600	GRAV	PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.7	500-1200	50 µg	10 µg	
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and BGI-4L at 2.2 lpm. 2020 NIC, TWA = 0.1mg/m <sup>3(R)</sup>		

## Mineral Oil (Oil mist)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8012-95-1	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	40-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		This method is not for oil mist containing PNAs; it does not collect vapor.		

## Mineral Oil, excluding Metal Working Fluids, Pure, highly and severely refined.

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8012-95-1	NIOSH 5026	FTIR	PVC, IOM	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	100-500	32 µg	16 µg	
Interferences		Comments		
Any aerosol (e.g. tobacco smoke) which absorbs infrared radiation near 2950 cm <sup>-1</sup> interferes.		Concentrated bulk oil sample required for analysis. Request IOMs one week before survey date. Rental charges for use of IOM samplers apply. This method is not covered under our AIHA-LAP, LLC scope of accreditation.		

## Mineral Oil, used in metal working

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8012-95-1	NIOSH 5524	GRAV	Pre-weighed PTFE1 (Zefon FPTFE137)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	768-960	75 µg	38 µg	
Interferences		Comments		
The method is non-specific and measures all substances extractable by organic solvents.		Bulk is recommended to test solubility. Refrigerate samples if unable to ship immediately. Ship cold. Samples must be analyzed within two weeks of collection.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Mineral Spirits (Stoddard Solvent)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-41-3	NIOSH 1550	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1.3-20	1.4 µg	0.71 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

### Mineral Spirits (Stoddard Solvent)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-41-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.3	15-480	2.1 µg	1.1 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

### Mineral Wool Fiber

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	40-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

### Molybdenum as Mo

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-98-7	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	0.10 µg	0.050 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)

## Molybdenum as Mo

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7439-98-7	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		40-240	0.025 µg	0.013 µg	Metals
Interferences		Comments			
		As part of the Lab’s QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

## Molybdenum as Mo

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7439-98-7	OSHA ID-121	ICP	Ghost wipe (225-2414)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
NA		NA	1.0 µg	0.5 µg	Metals 2
Interferences		Comments			
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab’s QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.			

## Morpholine

CAS #	Analytical Method	Analytical Technique	Sampling Media		
110-91-8	NIOSH S-150	GC-FID	SGT (SKC 226-10)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.1		10	2.8 µg	1.4 µg	dil acid
Interferences		Comments			

## Naphthalene

CAS #	Analytical Method	Analytical Technique	Sampling Media		
91-20-3	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2		1-10	0.47 µg	0.24 µg	CS <sub>2</sub>
Interferences		Comments			
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Naphthalene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
91-20-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.6	15-480	0.70 µg	0.35 µg	CS <sub>2</sub>
Interferences		Comments		

## Naphthalene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
91-20-3	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.32 µg	0.16 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

## Naproxen Sodium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
22204-53-1	LM-Pharma-4	LC-MS	PTFE5	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200	0.0010 µg	0.00050 µg	
Interferences		Comments		
		It is critical that the PTFE filter used for sampling must have a pore size of 1.0µm. Contact the lab to determine the minimum air volume required.		

## Naproxen Sodium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
22204-53-1	LM-Pharma-4	LC-MS	GFF Wipes	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
N/A	N/A	0.0020 µg	0.0010 µg	
Interferences		Comments		
		Follow "Naproxen Sodium Wipe Sampling Procedure. Order media one week ahead of survey. Media are prepared when ordered.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Nickel and inorganic compounds as Ni

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-02-0	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	50-1000	0.10 µg	0.05 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Nickel and inorganic compounds as Ni

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-02-0	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	40-240	0.16 µg	0.080 µg	Metals
Interferences		Comments		
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. <b>Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.</b>		

## Nickel and inorganic compounds as Ni

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-02-0	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe (SKC 225-2414)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
NA	NA	0.50 µg	0.25 µg	Metals2
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Nicotine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
54-11-5	NIOSH 2544	HPLC	XAD-2 (SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	60-400	2.0 µg	1.0 µg	
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Nitric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7697-37-2	NIOSH 7903	IC	MCE2, SGT** (SKC 225-19, SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2-0.5	5-100	2.3 µg	1.2 µg	
Interferences		Comments		
Particulate salts of the acid will give a positive interference.		Use the maximum flow rate for STEL sampling.		

## Nitric acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7697-37-2	NIOSH 7907	IC	SKC 225-9032	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	35-600	3.8 µg	1.9 µg	Acid1
Interferences		Comments		
Inorganic acids can react with co-sampled particulate matter on the pre-filter, leading to low results. Potentially interfering particulate chlorides and nitrates removed by the pre-filter can react with the sampled acids and liberate HCl and HNO <sub>3</sub> , which gets collected on the sampling filter, leading to high results.		Order media one week ahead of survey. Media are prepared when ordered. Ship and store cold. <b>Specialized filter. Media charge applies.</b>		

## Nitric Oxide and Nitrogen Dioxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10102-43-9	OSHA ID-190 (NO), OSHA ID-182 (NO <sub>2</sub> )	IC	SKC 226-40A	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	3-24	0.78 µg	0.39 µg	NO&NO <sub>2</sub>
Interferences		Comments		
		Use SKC 226-40A (2TEA coated tubes + oxidizer) to sample NO only or NO <sub>2</sub> and NO. Store and ship cold.		

## Nitroethane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-24-3	NIOSH 2526	GC-FID	XAD-2 (SKC 226-3002A) (Part A + Part B)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1.5-3	1.8 µg	0.90 µg	Ethyl Acetate
Interferences		Comments		
		Sample using 2 XAD-2, (front 600 mg and backup 300 mg), tubes in series. After sampling, separate and cap sorbent tubes.		

†(LPM) (CC/Min)

††(L) (Minutes)



## Nitrogen Dioxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10102-44-0	OSHA ID-182	IC	TEA-IMS (SKC 226-40-02)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2	35	1.2 µg	0.60 µg	
Interferences		Comments		
		Sample with TEA IMS (226-40-02) for NO <sub>2</sub> only. Store and ship cold.		

## Nitromethane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-52-5	NIOSH 2527	GC-FID	CS 106 (SKC 226-111A)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	5-10	2.3 µg	1.2 µg	Ethyl Actate
Interferences		Comments		

## Nitrous Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10024-97-2	OSHA ID-166	GC-ECD	AT N <sub>2</sub> O Monitor (X575AT)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.75	15-480	0.20 µg	NA	
Interferences		Comments		
Halogenated anesthetic gases, CFC's and HCFC's do not interfere when tested at their PELs.		Monitors must be received by Lab within one week after sampling and stored at controlled room temperature. Do not use after expiration date. This analysis is sub-contracted to an AIHA-LAP, LLC accredited lab.		

## Organic Solvent Scan

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	LM-GCMS-13			
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
Interferences		Comments		
		See List of Scans for list of individual organic solvents.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Oxalic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
144-62-7	OSHA ID-PV2115	IC	GFF (SKC 225-7)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.0	30-100	3.0 µg	1.5 µg	
Interferences		Comments		

## Ozone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10028-15-6	OSHA ID-214	IC	GFF, NaNO <sub>2</sub>	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.25-0.5	See Comments section.	4.8 µg	2.4 µg	
Interferences		Comments		
Particulate salts of nitrate and nitric acid give positive interferences for ozone. SO <sub>2</sub> will cause a negative interference.		Treated filter is stable for 1 month. Use flow rate of 1.5 lpm for short term sampling. For longer sampling time, use a flow rate of 0.5 lpm. For 25% of TLV, you will need 125L. Order media one week ahead of survey. Media are prepared when ordered. Sulfur dioxide (SO <sub>2</sub> ) has a negative interference on ozone. If SO <sub>2</sub> is suspected, check for its presence using a SO <sub>2</sub> detector tube. If present, request oxidizer tubes in addition to the ozone filters. Oxidizer tube precedes the ozone filter in the sampling train. After sampling seal the cassette with the end plugs and wrap each cassette with aluminum foil.		

## Palladium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-05-3	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	50-1000	0.32 µg	0.16 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Paraffin Wax Fume

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8002-74-2	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	100-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

### Particles (insoluble or poorly soluble) Not otherwise specified; inhalable

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	HSE MDHS 14	GRAV	PVC, IOM	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	40-960	100 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		Use IOM sampler with pre-weighed PVC. Contact Lab 1 week before intended use. The availability of IOM samplers is limited. Rental charge for the IOM samplers applies.		

### Particles (insoluble or poorly soluble) Not otherwise specified; respirable

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0600	GRAV	PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.7	100-816	50 µg	10 µg	
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

### Particles (insoluble or poorly soluble) Not otherwise specified; total

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	40-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Pentane(n-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-66-0	NIOSH 1500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	4	0.48 µg	0.24 µg	CS <sub>2</sub>
Interferences		Comments		

**Pentane(n-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-66-0	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.3	15-180	0.73 µg	0.37 µg	CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate front section of the monitor from the back section and cap immediately after sampling. Ship and store cold immediately.		

**Pentanedione(2,3-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
600-14-6	OSHA 1016	GC-MS	SGT/GFF-SGT/GFF (SKC 226-183)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	10L (TWA); 3L (short term)	0.59 µg	0.30 µg	95% EtOH
Interferences		Comments		
		Samples are collected on two specially washed silica gel tubes in series. Samples should be protected from the light during and after sampling. Separate and cap tubes after sampling. Sample separately from CS <sub>2</sub> compatible solvents. Order aluminum foil for wrapping the samples.		

**Pentanedione(2,4-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-54-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.7	15-300	2.7 µg	1.3 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Pentanone(2-) (Methyl propyl ketone)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-87-9	NIOSH 1300	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.70 µg	0.35 µg	CS <sub>2</sub>
Interferences		Comments		

**Pentanone(2-) (Methyl propyl ketone)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-87-9	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
33.0	15-480	1.1 µg	0.53 µg	CS <sub>2</sub>
Interferences		Comments		

**Peracetic Acid**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-21-0	NON 57	HPLC	SKC 225-9030 (Hydrogen Peroxide) + SKC-226-193-UC (Peracetic Acid)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1 LPM	15 Liters Max	0.50 µg/sample	N/A	Sample with hydrogen peroxide
Interferences		Comments		
Must sample for hydrogen peroxide		This analysis is sub-contracted to another laboratory. This method is not covered under the laboratory's AIHA-LAP, LLC scope of accreditation. Turnaround time is 10 business days.		

**Perchloroethylene (Tetrachloroethylene)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-18-4	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-40	1.8 µg	0.90 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

### Perchloroethylene (Tetrachloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-18-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.3	15-480	2.7 µg	1.4 µg	CS <sub>2</sub>
Interferences		Comments		

### Perfluorooctanoic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
335-67-1	DuPont-PFOA	LC-MS	OVS-2 (SKC 226-30-16)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.0	480	1.4 µg	0.70 µg	
Interferences		Comments		
		Sample 480 L if possible; 280 L will give 50% of recommended exposure limit.		

### Petroleum Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8032-32-4	NIOSH 1550	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1.3-20	0.50 µg	0.25 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

### Petroleum Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8032-32-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
33.1	15-480	0.75 µg	0.38 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Phenanthrene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
85-01-8	OSHA 58	HPLC	GFF (SKC 225-7)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	960	0.12 µg	0.060 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.		

**Phenanthrene (see PNA scan)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
85-01-8	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.12 µg	0.060 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

**Phenol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-95-2	OSHA 32	HPLC	XAD-7 (SKC 226-95)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	5-24	0.20 µg	0.10 µg	Phenol and cresol
Interferences		Comments		

**Phenylcyclohexene(4-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
4994-16-5	NIOSH 1500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-60	0.38 µg	0.18 µg	CS <sub>2</sub>
Interferences		Comments		
		Please indicate if samples are for "LEED" compliance. Minimum sample volume for "LEED" samples is 60 L. Preferred air sampling method.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Phenylcyclohexene (4-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
4994-16-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
20.3	15-480	0.57 µg	0.28 µg	CS <sub>2</sub>
Interferences		Comments		
		Please indicate if samples are for "LEED" compliance. Minimum sample volume for "LEED" samples is 60 L.		

### Phenylene(1,3-) diamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-45-2	OSHA 87	HPLC	GFF, Acid	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	10-1000	0.83 µg	0.42 µg	
Interferences		Comments		
		Sampling media has a short shelf-life so it is not kept in stock. Contact lab one week before sampling to order media.		

### Phosphine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7803-51-2	OSHA 1003	ICP	GFF- PE,HgCl <sub>2</sub> (SKC 225-9018)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	250	3.6 µg	1.8 µg	
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Minimum air volume required at ¼ of TLV is 250. Sampling media has a shelf-life of only two weeks so it is not kept in stock. Contact lab one week before sampling to order media. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

### Phosphoric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7664-38-2	NIOSH 7903	IC	MCE2, SGT** (SKC 225-19, SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2-0.5	25-100	2.3 µg	1.2 µg	Acid1
Interferences		Comments		
Particulate salts of the acid will give a positive interference.		Sample at a flow rate of 0.5 lpm for STEL.		

†(LPM) (CC/Min)

††(L) (Minutes)



## Phosphoric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7664-38-2	NIOSH 7908	IC	SKC 225-9033	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-5	50-1000	3.8 µg	1.9 µg	
Interferences		Comments		
Particulate salts of sulfate or phosphate will give positive interference.		Ship and store cold. Specialty filter. Media charge applies.		

## Phosphorus (elements)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7723-14-0	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	330-1000	3.3 µg	1.7 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Phthalic Anhydride

CAS #	Analytical Method	Analytical Technique	Sampling Media	
85-44-9	OSHA 90	HPLC	GFF, Vamine	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.0	75	0.13 µg	0.065 µg	
Interferences		Comments		
Isocyanates, acid chlorides and other anhydrides will give a positive interference.		Sample open-faced. Order filters one week ahead, filters are prepared when ordered and have one month shelf-life.		

## Piperazine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-85-0	OSHA In-house IMIS P250	HPLC	XAD-2, NITC (SKC 226-30-18)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	10	0.12 µg	0.060 µg	
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Platinum Metal and Soluble Salts as Pt

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-06-4	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	500-1000	0.25 µg	0.13 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.			

## PNA Scan (NIOSH 5506)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
2	200-1000				
Interferences		Comments			
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold. Analysis price is for the filter and tube together.			

## PNA Scan (OSHA 58)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	OSHA 58	HPLC	GFF (SKC 225-7)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
2	960				
Interferences		Comments			
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.			

## Polychlorobiphenyl (Chlorodiphenyl, 54% Chlorine) (PCB)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
11097-69-1	NIOSH 5503	GC-MS	GFF-Florisil (Millipore SX0001300/01/ AP2001300 SKC 226-39)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2		1-50	1.1 µg	0.55 µg	
Interferences		Comments			
Other chlorinated pesticides may interfere in the quantification of PCB.					

†(LPM) (CC/Min)

††(L) (Minutes)

### Polychlorobiphenyl (Chlorodiphenyl, 42% Chlorine) (PCB)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
53469-21-9	NIOSH 5503	GC-MS	GFF-Florisil (Millipore SX0001300/01/ AP2001300 SKC 226-39)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2		1-50	0.97 µg	0.49 µg	
Interferences		Comments			
Other chlorinated pesticides may interfere in the quantification of PCB.					

### Polyvinyl Chloride (PVC)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
9002-86-2	NIOSH 0600	GRAV	PVC (SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1.7	100-816	50 µg	10 µg		
Interferences		Comments			
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.			

### Portland Cement

CAS #	Analytical Method	Analytical Technique	Sampling Media		
65997-15-1	NIOSH 0600	GRAV	PVC (SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1.7 (See comment)	500-816	50 µg	10 µg		
Interferences		Comments			
All other dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.			

### Potassium Hydroxide

CAS #	Analytical Method	Analytical Technique	Sampling Media		
1310-58-3	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	85-100	17 µg	8.5 µg	Metals	
Interferences		Comments			
All forms of potassium are quantified. Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

†(LPM) (CC/Min)

††(L) (Minutes)

**Propanol(n-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-23-8	NIOSH 1401	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.83 µg	0.42 µg	1%IPA/CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

**Propanol(n-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-23-8	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
39.7	15-420	1.2 µg	0.60 µg	MC CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

**Propionaldehyde**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-38-6	NIOSH 2016	HPLC	AT Monitor (N571AT)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
9.58	15-480	0.029 µg	0.015 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight.		

**Propionaldehyde**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-38-6	NIOSH 2016	HPLC	Sep-Pak (WAT047205)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	10-100	0.15 µg	0.075 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight. <b>Preferred for STEL sampling.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)

## Propionaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-38-6	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	1-15	0.059 µg	0.030 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight. Preferred for STEL sampling.		

## Propionic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-09-4	NIOSH 2011	IC	PTFE3-SGT** (SKC 225-17A, SKC 226-10-03)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.5	15-100	2.3 µg	1.2 µg	Acid2
Interferences		Comments		
		Do not sample with inorganic acids.		

## Propoxyethanol(2-) (Ethylene glycol monopropyl ether)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2807-30-9	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1-10	1.1 µg	0.55 µg	5%MeOH/MC
Interferences		Comments		

## Propoxyethanol(2-) (Ethylene glycol monopropyl ether)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2807-30-9	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
29.4	15-480	1.7 µg	0.85 µg	MC CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Propyl Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-94-5	OSHA PV2061	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	12	1.3 µg	0.65 µg	1%DMF/CS <sub>2</sub>
Interferences		Comments		

## Propyl Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-94-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.7	15-480	2.0 µg	1.0 µg	CS <sub>2</sub>
Interferences		Comments		

## Propyl(n-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-60-4	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	1-10	0.78 µg	0.39 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Use a flow rate of 0.2 lpm for STEL.		

## Propyl(n-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-60-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.1	15-480	1.2 µg	0.60 µg	CS <sub>2</sub>
Interferences		Comments		

## Propyl(n-) Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-23-8	NIOSH 1401	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.83 µg	0.42 µg	1%IPA/CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Propyl(n-) Alcohol**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-23-8	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
39.7	15-420	1.2 µg	0.60 µg	MC CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold overnight.		

**Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-98-2	NIOSH 1403	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1-10	1.3 µg	0.65 µg	5%MeOH/MC
Interferences		Comments		

**Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-98-2	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.4	15-480	2.0 µg	1.0 µg	MC CS <sub>2</sub>
Interferences		Comments		

**Propylene Glycol Monomethyl Ether Acetate (PGMEA)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-65-6	NIOSH 1450	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	1.1 µg	0.55 µg	CS <sub>2</sub>
Interferences		Comments		

**Propylene Glycol Monomethyl Ether Acetate (PGMEA)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-65-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
25.2	15-480	1.7 µg	0.85 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

### Propylene Glycol (1,2-Propanediol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
57-55-6	NIOSH 5523	GC-FID	OVS 7 (SKC 226-57)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.5-2	5-60	6.2 µg	3.1 µg	MeOH
Interferences		Comments		

### Propylene Oxide (1,2-Epoxypropane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-56-9	ASTM D5578-04	GC-FID	ORBO 78 (SUPELCO 20355)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.15	1-24	0.51 µg	0.25 µg	ACN/TOL
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents. Store and ship cold overnight.		

### Propylene Oxide (1,2-Epoxypropane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-56-9	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
37.7	15-480	0.76 µg	0.38 µg	
Interferences		Comments		
		Use 3M 3520. Separate front section of the monitor from the back section and cap immediately after sampling. Sample separately from CS <sub>2</sub> compatible solvents. Store and ship cold overnight.		

### Pyrene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
129-00-0	OSHA 58	HPLC	GFF (SKC 225-7)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	960	0.60 µg	0.30 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.		

†(LPM) (CC/Min)

††(L) (Minutes)



## Pyrene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
129-00-0	NIOSH 5506	HPLC	PTFE2/XAD-2 (PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000	0.60 µg	0.30 µg	PNAs
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

## Pyrethrum

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8003-34-7	NIOSH 5008	HPLC	GFF (SKC 225-7)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	20-400	0.13 µg	0.065 µg	
Interferences		Comments		

## Pyridine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-86-1	NIOSH 1613	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-1.0	18-150	0.23 µg	0.12 µg	MC
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents.		

## Resin Acids

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8050-09-7	LM-LC-30	HPLC	GFF, IOM	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200 (minimum)	0.091 µg	0.045 µg	
Interferences		Comments		
		2019 NIC Resin acids as total resin acids, TWA = 0.001mg/cu m (l), DSEN;RSEN adopted in 2020. Resin Acids include Abietic Acid and Dehydroabietic Acid. 200-L will give you 46% (0.00046 mg/cu m) of TLV. 400-L will give you 23% (0.00023 mg/cu m) of TLV. 900-L will give you 10% (0.00010 mg/cu m) of TLV.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Resorcinol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-46-3	OSHA PV2053	GC-FID	OVS 7 (SKC 226-57)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	20-400	5.0 µg	2.5 µg	MeOH
Interferences		Comments		

## Rhodium as Rh

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-16-6	OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	500-1000	0.50 µg	0.25 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

## Scan for Aldehydes

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 2016	HPLC		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
Interferences		Comments		
		See List of Scans for individual aldehydes. Ship and store cold.		

## Scan for Aliphatic Amines

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 2010	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-1	5-30			Amine1
Interferences		Comments		
		Ethylamine, diethylamine and triethylamine can also be analyzed using this method. Please call Lab for other types of amines.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Scan for Anesthetic Gases

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.05		12			
Interferences		Comments			
		Analyzes enflurane, halothane and isoflurane. See individual anesthetic gases. Please call lab for additional anesthetic gases.			

### Scan for Aromatic Amines

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	NIOSH 2002	GC-FID	SGT (SKC 226-10)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.02-0.5	5-30			Amine3	
Interferences		Comments			
		Analyze for aniline, methyl aniline and o-toluidine. Please call Lab for other types of amines.			

### Scan for Inorganic Acids

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	NIOSH 7903	IC	MCE2, SGT** (SKC 225-19, SKC 226-10-03)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.2-0.5	50-100				
Interferences		Comments			
		See List of Scans for individual inorganic acids/anions.			

### Scan for Isocyanates

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	OSHA 42	HPLC	GFF, 1-2PP		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
Interferences		Comments			
		Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months if kept cold. See List of Scans at the back of the guide for list of individual isocyanate.			

†(LPM) (CC/Min)

††(L) (Minutes)

### Scan for Organic Solvents

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	LM-GCMS-13	GC-FID GC-MS		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
Interferences		Comments		
		See List of Scans for individual organic solvents.		

### Scan for PNAs (NIOSH 5506)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 5506	HPLC	PTFE2/XAD-2 (SKC PALL P5PJ037, SKC 226-30-04)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	200-1000			
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold. Price is for the analysis of filter and tube together.		

### Scan for PNAs (OSHA 58)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	OSHA 58	HPLC	GFF (SKC 225-7)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	960			
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.		

### Selenium and Compounds as Se

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7782-49-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	120-1000	1.0 µg	0.50 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)

## Selenium and Compounds as Se

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7782-49-2	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		75-240	0.75 µg	0.38 µg	Metals
Interferences		Comments			
		As part of the Lab’s QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

## Sevoflurane (Sevofrane)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
28523-86-6	OSHA 106	GC-FID	CT (SKC 226-01, -09)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.05		3	2.9 µg	1.5 µg	CS <sub>2</sub>
Interferences		Comments			
		Store and ship cold overnight.			

## Sevoflurane (Sevofrane)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
28523-86-6	3M Method	GC-FID	OVM (3M 3500)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
27.3		15-480	4.4 µg	2.2 µg	CS <sub>2</sub>
Interferences		Comments			
		Store and ship cold overnight.			

## Silica Cristobalite

CAS #	Analytical Method	Analytical Technique	Sampling Media		
14464-46-1	NIOSH 7500	XRD	PVC (SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
see comments		600-1200	7.5 µg	5.0 µg	Silica
Interferences		Comments			
Mica, potash, feldspars, zircon, graphite and aluminosilicates will interfere. Bulk sample is required for interference check.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm; 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm. Dusty atmospheres require much smaller sample volumes (<600 L) to obtain optimum dust loading of 2 mg on filter. 600-L will give you 52% (0.013 mg/cu m) of TLV and 900-L for 33% (0.0083 mg/cu m) of TLV. Sample 600-L to quantify at 25% of the new OSHA PEL.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Silica Quartz

CAS #	Analytical Method	Analytical Technique	Sampling Media	
14808-60-7	NIOSH 7500	XRD	PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
see comments	600-1200	7.5 µg	5.0 µg	Silica
Interferences		Comments		
Mica, potash, feldspars, zircon, graphite and aluminosilicates will interfere. Bulk sample is required for interference check.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm; 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm. Dusty atmospheres require much smaller sample volumes (<600 L) to obtain optimum dust loading of 2 mg on filter. Sample 600-L for 52% (0.013 mg/cu m) of the TLV and 900-L for 33% (0.0083 mg/cu m) of TLV. Sample 600-L to quantify at 25% of the new OSHA PEL.		

## Silver Metal and Soluble Compounds as Ag

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-22-4	NIOSH 7300	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	250-1000	0.50 µg	0.25 µg	Metals1
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Soapstone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	34-7200	10 µg	5.0 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

## Soapstone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0600	GRAV	PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.7	67-816	50 µg	25 µg	
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (British Medical Research Council-SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Sodium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-23-5	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	60-1000	7.5 µg	3.8 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Sodium Hydroxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1310-73-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	60-1000	13 µg	6.5 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of sodium are quantified. Use a flow rate of 2 lpm for STEL. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Starch

CAS #	Analytical Method	Analytical Technique	Sampling Media	
9005-25-8	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	20-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 LPM, for area sampling up to 15 LPM.		

## Stoddard Solvent

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-41-3	NIOSH 1550	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1.3-20	2.3 µg	1.1 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Stoddard Solvent

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-41-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.3	15-480	3.4 µg	1.7 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

## Strontium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-24-6	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	620-1000	0.052 µg	0.026 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Strontium Chromate as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7789-06-2	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	620-1000	0.031 µg	0.016 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Samples are analyzed as strontium and calculated as strontium chromate as Cr. Other forms of strontium will interfere. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Styrene (Vinyl benzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-42-5	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.02-1	1-10	0.43 µg	0.22 µg	CS <sub>2</sub>
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%. Store and ship cold overnight. 2019 NIC, TWA = 10ppm, STEL = 20 ppm, OTO, A3, BEI adopted in 2020.		

†(LPM) (CC/Min)

††(L) (Minutes)



## Styrene (Vinyl benzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-42-5	3M Method	GC-FID	OVM (3M 3500) SKC 575-006	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.9 13.55	15-480 15-240	0.65 µg	0.33 µg	CS <sub>2</sub> Tol
Interferences		Comments		
		Store and ship cold overnight. 2019 NIC, TWA = 10ppm, STEL = 20 ppm, OTO, A3, BEI adopted in 2020.		

## Sulfur Dioxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7446-09-5	OSHA ID-200	IC	IABC (SKC 226-80)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1	10-50	3.1 µg	1.6 µg	
Interferences		Comments		
Particulate salts of sulfate, sulfur trioxide and sulfuric acid may give positive interferences for sulfur dioxide.				

## Sulfuric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7664-93-9	OSHA ID-113	IC	PPI Thoracic-MCE (225-3861)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2.0	90-230	4.6 µg	2.3 µg	Acid1
Interferences		Comments		
Particulate salts of sulfate will give a positive interference.		The TLV for sulfuric acid is as thoracic particulate mass. Use 37-mm MCE 3 piece cassette for BGI GK2.69 cyclone at 1.6 lpm or request thoracic PPI pre-loaded with 0.8µm MCE filter. Media charge applies.		

## Sulfuric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7664-93-9	NIOSH 7903	IC	MCE2, SKC 226-10-03	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.2-0.5	50-150	4.6 µg	2.3 µg	
Interferences		Comments		
Particulate salts of sulfate will give a positive interference.		If using 226-10-03 in a very humid environment, use two sorbent tubes in series.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Sulfuric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7664-93-9	NIOSH 7908	IC	SKC 225-9033	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-5	50-500	4.6 µg	2.3 µg	
Interferences		Comments		
Particulate salts of sulfate or phosphate will give positive interference.		Ship and store cold. Specialty filter. Media charge applies.		

## Synthetic Vitreous Fibers

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 7400	PCM	MCE, 25 mm (ZEFON Z008BA)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.5-16	50-720	0.050 fiber/fld	0.01 fib/fld	
Interferences		Comments		
Chain-like particles may appear fibrous and high levels of non-fibrous dust particles may obscure fibers.		Adjust sampling flow rate and time to obtain optimum fiber loading on the filter. Do not overload filter. Sample open faced. When shipping your samples, do not pack with untreated polystyrene as can lead to fiber loss from electrostatic effect.		

## Talc

CAS #	Analytical Method	Analytical Technique	Sampling Media	
14807-96-6	NIOSH 0600	GRAV	PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1.7	100-816	50 µg	10 µg	
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (British Medical Research Council-SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm. TWA of 2 mg/cu m as respirable fraction is for particulate matter containing no asbestos and <1% crystalline silica.		

## Tantalum and Tantalum Oxide Dust as Ta

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-25-7	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	100-1000	50 µg	25 µg	
Interferences		Comments		
All other dusts will interfere.				

†(LPM) (CC/Min)

††(L) (Minutes)

## Tellurium and Compounds as Te

CAS #	Analytical Method	Analytical Technique	Sampling Media	
13494-80-9	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	160-1000	0.54 µg	0.27 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Excludes hydrogen telluride. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

## Tetrachloroethylene (Perchloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-18-4	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-40	1.8 µg	0.90 µg	CS <sub>2</sub>
Interferences		Comments		

## Tetrachloroethylene (Perchloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-18-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
28.3	15-480	2.7 µg	1.4 µg	CS <sub>2</sub>
Interferences		Comments		

## Tetrahydrofuran [THF]

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-99-9	NIOSH 1609	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-9	0.83 µg	0.42 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm. High humidity may greatly decrease breakthrough volume.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Tetrahydrofuran

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-99-9	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
37.2	15-480	1.2 µg	0.62 µg	CS <sub>2</sub>
Interferences		Comments		

## Thallium and Compounds, as Tl

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-28-0	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	PVC, IOM MCE, IOM	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	500-1000	1.0 µg	0.50 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. Contact Lab 1 week before intended use. The availability of IOM samplers is limited. rental charge for IOM samplers applies.		

## Thallium and Compounds, as Tl

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-28-0	NIOSH 7301 NIOSH 7303	ICP-MS	PVC, IOM MCE, IOM	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	750-1000	1.3 µg	0.65 µg	Metals
Interferences		Comments		
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples. Contact Lab one week before intended use. The availability of IOM samplers is limited. Rental charge for IOM samplers applies.		

## Thiram

CAS #	Analytical Method	Analytical Technique	Sampling Media	
137-26-8	NIOSH 5005	HPLC	PTFE1 (SKC 225-1705)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	20-400	0.79 µg	0.40 µg	
Interferences		Comments		
		Return to Lab immediately after sampling.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Tin and Compounds as Sn

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-31-5	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	0.52 µg	0.26 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Includes metal, metal oxide, inorganic tin compounds (except tin hydride) and organic compounds. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Tin and Compounds as Sn

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-31-5	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	40-240	0.025 µg	0.013 µg	Metals
Interferences		Comments		
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. <b>Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples..</b>		

## Tin Organic Compounds as Sn

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-31-5	NIOSH 5504	GFAA	GFF/XAD-2 (SKC 225-7) (SKC 226-30)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-1.5	50-500	0.050 µg	0.025 µg	Metals
Interferences		Comments		
All forms of tin organic compounds are quantified.		Ship assembled sampler in dry ice overnight. This method is not covered under our AIHA-LAP, LLC scope of accreditation.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Titanium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-32-6	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	0.25 µg	0.13 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Titanium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-32-6	NIOSH 7301	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-500	0.15 µg	0.075 µg	Metals
Interferences		Comments		
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.		

## Titanium Dioxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
13463-67-7	OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	0.42 µg	0.21 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of titanium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Toluene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-88-3	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.44 µg	0.22 µg	CS <sub>2</sub>
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%. 2020 NIC, TWA = 20ppm, OTO;A4;BEi.		

## Toluene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-88-3	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.4	15-480	0.67 µg	0.33 µg	CS <sub>2</sub>
Interferences		Comments		
		2020 NIC, TWA = 20ppm, OTO;A4;BEi.		

## Toluene-2,4-diisocyanate (2,4-TDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
584-84-9	OSHA 42	HPLC	GFF, 1-2PP	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	15-240	0.010 µg	0.0050 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months if kept cold.		

## Toluene-2,4-diisocyanate (2,4-TDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
584-84-9	OSHA 42	HPLC	GFF Wipes	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
N/A	N/A	0.010 µg	0.0050 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 1 month if kept cold.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Toluene-2,6-diisocyanate (2,6-TDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
91-08-7	OSHA 42	HPLC	GFF, 1-2PP	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	15-240	0.010 µg	0.0050 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months if kept cold.		

### Toluene-2,6- diisocyanate (2,6-TDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
91-08-7	OSHA 42	HPLC	GFF Wipes	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
N/A	N/A	0.010 µg	0.0050 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead, media is prepared when ordered. Derivatizing solution has a shelf life of 1 month if kept cold.		

### Toluidine(o-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-53-4	NIOSH 2002	GC-FID	SGT (SKC 226-10)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.02-0.5	5-30	1.4 µg	0.70 µg	Amine3
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

### Tributyl Phosphate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
126-73-8	NIOSH 5034	GC-FID	MCE (SKC 225-5)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-3	2-100	0.41 µg	0.20 µg	Ethyl Ether
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents.		

†(LPM) (CC/Min)

††(L) (Minutes)



**Trichloro(1,1,2)-1,2,2-trifluoroethane**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
76-13-1	NIOSH 1020	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.05	1-10	2.9 µg	1.5 µg	CS <sub>2</sub>
Interferences		Comments		

**Trichloro(1,1,2)-1,2,2-trifluoroethane**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
76-13-1	3M Method	GC-FID	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
29.1	15-480	4.4 µg	2.2 µg	CS <sub>2</sub>
Interferences		Comments		

**Trichlorobenzene(1,2,4-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
120-82-1	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.29 µg	0.15 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

**Trichlorobenzene(1,2,4-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
120-82-1	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.8	15-480	0.44 µg	0.22 µg	CS <sub>2</sub>
Interferences		Comments		
		Use 3M 3520. Separate the front section of the monitor from the back section and cap immediately after sampling.		

**Trichloroethane(1,1,1-) (Methyl Chloroform)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-55-6	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	2.2 µg	1.1 µg	CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

†(LPM) (CC/Min)

††(L) (Minutes)

**Trichloroethane(1,1,1-) (Methyl Chloroform)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-55-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
30.9	15-480	3.3 µg	1.7 µg	CS <sub>2</sub>
Interferences		Comments		

**Trichloroethane(1,1,2-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-00-5	NIOSH 1003	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	2-60	1.6 µg	0.80 µg	CS <sub>2</sub>
Interferences		Comments		

**Trichloroethane(1,1,2-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-00-5	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
29.7	15-480	2.9 µg	1.4 µg	CS <sub>2</sub>
Interferences		Comments		

**Trichloroethylene**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-01-6	NIOSH 1022	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	2.0 µg	1.0 µg	CS <sub>2</sub>
Interferences		Comments		

**Trichloroethylene**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-01-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
31.1	15-480	3.0 µg	1.5 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

## Triethanolamine

CAS #	Analytical Method	Analytical Technique	Sampling Media		
102-71-6	NIOSH 2007	IC	ORBO 53 or SGT** (SUPELCO 20265) (SKC 226-10-03)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.5		15-25	7.5 µg	3.8 µg	EA
Interferences		Comments			
		Store in freezer after sampling. Ship and store cold.			

## Triethylamine

CAS #	Analytical Method	Analytical Technique	Sampling Media		
121-44-8	NIOSH 2010	GC-FID	SGT (SKC 226-10)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-1		5-30	0.29 µg	0.15 µg	Amine1
Interferences		Comments			
Nitrogen compounds that co-elute will interfere.					

## Triethylamine

CAS #	Analytical Method	Analytical Technique	Sampling Media		
121-44-8	OSHA PV2060	GC-FID	XAD-7, Acid (SKC 226-98)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
0.1		10	4.0 µg	2.0 µg	
Interferences		Comments			
Nitrogen compounds that co-elute will interfere.					

## Triethylenetetramine

CAS #	Analytical Method	Analytical Technique	Sampling Media		
112-24-3	OSHA 60	HPLC	XAD-2, NITC (SKC 226-30-18)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
0.1	10	0.24 µg	0.12 µg	Amine2	
Interferences		Comments			
Nitrogen compounds that co-elute will interfere.					

†(LPM) (CC/Min)

††(L) (Minutes)

**Triglycidyl Isocyanurate(1,3,5)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2451-62-9	Ciba-Geigy C321A	GC-MS	PTFE1 or PTFE5	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1	100-240	4.0 µg	2.0 µg	Acetone
Interferences		Comments		
		Samples can be collected using either 25mm or 37mm filters.		

**Trimellitic Anhydride**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
552-30-7	OSHA 98	HPLC	GFF, Vamine	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	480	0.050 µg	0.025 µg	
Interferences		Comments		
		Order media one week ahead of survey. Media are prepared when ordered. Sample open-faced.		

**Trimethylbenzene(1,2,4-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-63-6	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.44 µg	0.22 µg	CS <sub>2</sub>
Interferences		Comments		

**Trimethylbenzene(1,2,4-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-63-6	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
24.4	15-480	0.66 µg	0.33 µg	CS <sub>2</sub>
Interferences		Comments		

**Trimethylbenzene(1,3,5-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-67-8	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.45 µg	0.23 µg	CS <sub>2</sub>
Interferences		Comments		

†(LPM) (CC/Min)

††(L) (Minutes)

**Trimethylbenzene(1,3,5-)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-67-8	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
26.3	15-480	0.68 µg	0.34 µg	CS <sub>2</sub>
Interferences		Comments		

**Tungsten and Compounds as W (in the absence of Cobalt)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-33-7	OSHA ID-213	ICP	MCE (SKC 225-5)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	30-480	0.25 µg	0.12 µg	
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Recommended air volume for STEL is 30L. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

**Tungsten, as W Soluble Compounds**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-33-7	OSHA ID-213	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
2	30-480	0.25 µg	0.12 µg	
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Recommended air volume for STEL is 30L. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

**TVOC as n-Hexane**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 1500	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-10	0.40 µg	0.20 µg	CS <sub>2</sub>
Interferences		Comments		
		For IAQ and LEED sampling where TVOC and formaldehyde are collected, do not sample with Charcoal Tubes and DNPH tubes in tandem. The sorbent in the DNPH tubes may off-gas acetonitrile which can cause a positive interference in the TVOC results.		

†(LPM) (CC/Min)

††(L) (Minutes)

**TVOC as n-Hexane**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
32.0	15-480	0.60 µg	0.30 µg	CS <sub>2</sub>
Interferences		Comments		

**Valeraldehyde**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-62-3	NIOSH 2016	HPLC	AT Monitor (N571AT)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
7.21	15-480	0.043 µg	0.022 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight.		

**Valeraldehyde**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-62-3	NIOSH 2016	HPLC	Sep-Pak (WAT047205)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	10-100	0.22 µg	0.11 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight. Preferred for STEL sampling. Sample at 1.5 lpm for STEL.		

**Valeraldehyde**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-62-3	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-1.5	1-15	0.087 µg	0.044 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight. Preferred for STEL sampling. Sample at 1.5 lpm for STEL.		

†(LPM) (CC/Min)

††(L) (Minutes)

### Vanadium Pentoxide as V

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-62-2	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		40-480	0.20 µg	0.10 µg	Metals
Interferences		Comments			
		As part of the Lab’s QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

### Vanadium Pentoxide as V

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-62-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	120-1000	0.10 µg	0.050 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of vanadium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.			

### Vanadium

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-62-2	OSHA ID-121	ICP	Ghost wipe ( 225-2414)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
NA		NA	0.22 µg	0.11 µg	Metals 2
Interferences		Comments			
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.			

### Vegetable Oil Mist

CAS #	Analytical Method	Analytical Technique	Sampling Media		
	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-15	20-7200	50 µg	10 µg		
Interferences		Comments			
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm. TLV withdrawn.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Vinyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-05-4	NIOSH 1453	GC-FID	CT(SKC 226-01,-09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.1-0.2	0.75-24	1.2 µg	0.60 µg	5%MeOH/MC CS <sub>2</sub>
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

## Vinyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-05-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.8	15-480	1.8 µg	0.90 µg	CS <sub>2</sub>
Interferences		Comments		

## Vinyl Chloride (Chloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-01-4	NIOSH 1007	GC-FID GC-MS	CT-CT (SKC 226-01)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05	0.7-5	0.29 µg	0.14 µg	CS <sub>2</sub>
Interferences		Comments		
		Sample using 2 charcoal tubes in series. Separate and cap tubes before shipping. Possible loss of sample after 2 or more weeks of storage at room temperature. Ship and store cold immediately.		

## Vinyl Chloride (Chloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-01-4	3M Method	GC-FID GC-MS	OVM (3M 3520)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
40.8	15-480	0.43 µg	0.22 µg	CS <sub>2</sub>
Interferences		Comments		
		Use OVM 3520. Separate front and back section of the monitor immediately after sampling. Ship cold immediately.		

†(LPM) (CC/Min)

††(L) (Minutes)



**Vinyl(1-)-2-pyrrolidinone**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
88-12-0	NIOSH 1302	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.05-0.2	3-125	0.39 µg	0.20 µg	5%MeOH/MC
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents.		

**Vinyl(1-)-2-pyrrolidinone**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
88-12-0	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
26.7	120-480	0.59 µg	0.30 µg	MC
Interferences		Comments		
		Sample separately from CS <sub>2</sub> compatible solvents. Sampling for 480 minutes allows quantification at 40% of the TLV.		

**Vinylidene Chloride (1,1-Dichloroethylene)**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-35-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
35.1	15-480	0.60 µg	0.30 µg	CS <sub>2</sub>
Interferences		Comments		
		Store and ship cold.		

**VM & P Naphtha**

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8032-32-4	NIOSH 1550	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1.3-28	1.8 µg	0.90 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

†(LPM) (CC/Min)

††(L) (Minutes)

## VM & P Naphtha

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8032-32-4	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
33.2	15-480	2.7 µg	1.4 µg	CS <sub>2</sub>
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

## Welding Fume Scan

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	400-1000			
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		See List of Scans for individual metals. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

## Welding Fumes, Total

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	40-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

## Wood Dust

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	NIOSH 0500	GRAV	Pre-weighed PVC (SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-15	40-7200	50 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Wood Dust

CAS #	Analytical Method	Analytical Technique	Sampling Media	
	HSE MDHS-14	GRAV	PVC, IOM	
Sampling Rate†	Sampling Volume††	LOQ	LOQ	Compatibility Code
2	960	100 µg	10 µg	
Interferences		Comments		
All other dusts will interfere.		Use IOM sampler with pre-weighed PVC. Contact Lab one week before intended use. The availability of IOM samplers is limited. Rental charge for the IOM samplers applies.		

## Xylene (Dimethyl benzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1330-20-7	NIOSH 1501	GC-FID	CT (SKC 226-01, -09)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
0.01-0.2	1-30	1.0 µg	0.50 µg	CS <sub>2</sub>
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%.		

## Xylene (Dimethyl benzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1330-20-7	3M Method	GC-FID	OVM (3M 3500)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
27.3	15-480	1.5 µg	0.75 µg	CS <sub>2</sub>
Interferences		Comments		

## Yttrium and compounds, as Y

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-65-5	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	25-1000	0.025 µg	0.013 µg	Metals
Interferences		Comments		
		Spectral interferences are the primary interferences encountered in ICP-AES analysis.		

†(LPM) (CC/Min)

††(L) (Minutes)

## Zinc

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-66-6	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code	
1-4	120-1000	3.0 µg	1.5 µg	Metals	
Interferences		Comments			
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab’s QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>			

## Zinc

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-66-6	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
1-4		40-240	0.75 µg	0.38 µg	Metals
Interferences		Comments			
		As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.			

## Zinc

CAS #	Analytical Method	Analytical Technique	Sampling Media		
7440-66-6	OSHA ID-121	ICP	Ghost wipe ( 225-2414)		
Sampling Rate†		Sampling Volume††	LOQ	LOD	Compatibility Code
NA		NA	160 ug	80 ug	Metals 2
Interferences		Comments			
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.			

†(LPM) (CC/Min)

††(L) (Minutes)

## Zinc Chloride Fume

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7646-85-7	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	120-1000	6.3 µg	3.1 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Analysis is for water soluble zinc compounds. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Zinc Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1314-13-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	120-1000	3.7 µg	1.9 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Samples are analyzed for zinc and calculated as zinc oxide. Current TLV is for respirable sampling. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. <b>Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.</b>		

## Zinc Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1314-13-2	NIOSH 7301 NIOSH 7303	ICP-MS	MCE or PVC (SKC 225-5 or SKC 225-5-37-P)	
Sampling Rate†	Sampling Volume††	LOQ	LOD	Compatibility Code
1-4	40-240	0.93 µg	0.46 µg	Metals
Interferences		Comments		
		Samples are analyzed for zinc and calculated as zinc oxide. Current TLV is for respirable sampling. As part of the Lab's QC protocol, yttrium, rhodium, and lutetium are used as internal standards in ICP-MS analysis. <b>Please indicate in your sample submission form if yttrium, rhodium, and/or lutetium are present in the area where you collected your samples.</b>		

†(LPM) (CC/Min)

††(L) (Minutes)

## List of Analytes by CAS #

50-00-0	Formaldehyde	75-35-4	Vinylidene chloride (1,1-Dichloroethylene)
50-21-5	Lactic acid	75-52-5	Nitromethane
50-32-8	Benzo[a]pyrene	75-56-9	Propylene oxide
53-70-3	Dibenzo[ah]anthracene (see PNA scan)	75-65-0	Butyl(tert-) alcohol
54-11-5	Nicotine	76-13-1	Trichloro(1,1,2-)-1,2,2- trifluoroethane
56-23-5	Carbon tetrachloride	76-22-2	Camphor
56-55-3	Benz[a]anthracene	77-73-6	Dicyclopentadiene
57-55-6	Propylene glycol (1,2-Propanediol)	78-59-1	Isophorone
60-29-7	Ethyl ether	78-83-1	Isobutyl alcohol
62-53-3	Aniline	78-92-2	Butyl(sec-) alcohol
63-25-2	Carbaryl (SEVIN)	78-93-3	Butanone(2-) (Methyl ethyl ketone)
64-17-5	Ethyl alcohol (Ethanol)	78-94-4	Methyl vinyl ketone
64-18-6	Formic acid	79-00-5	Trichloroethane(1,1,2-)
64-19-7	Acetic acid	79-01-6	Trichloroethylene
67-56-1	Methyl alcohol (Methanol)	79-06-1	Acrylamide
67-63-0	Isopropyl alcohol	79-08-3	Bromoacetic acid
67-64-1	Acetone	79-09-4	Propionic acid
67-66-3	Chloroform	79-10-7	Acrylic Acid
68-12-2	Dimethylformamide	79-11-8	Chloroacetic acid
71-23-8	Propyl(n-) alcohol	79-20-9	Methyl acetate
71-36-3	Butyl(n-) alcohol	79-21-0	Peracetic Acid
71-43-2	Benzene	79-24-3	Nitroethane
71-55-6	Trichloroethane(1,1,1-) (Methyl Chloroform)	79-43-6	Dichloroacetic acid
74-90-8	Hydrogen cyanide	80-05-7	Bisphenol A
75-01-4	Vinyl chloride	80-62-6	Methyl methacrylate
75-04-7	Ethyl amine	83-32-9	Acenaphthene
75-05-8	Acetonitrile	84-66-2	Diethyl phthalate
75-07-0	Acetaldehyde	84-74-2	Dibutyl phthalate
75-09-2	Methylene chloride (Dichloromethane)	85-01-8	Phenanthrene
75-12-7	Formamide	85-44-9	Phthalic anhydride
75-15-0	Carbon disulfide	86-73-7	Fluorene
75-18-3	Dimethyl sulfide	88-12-0	Vinyl(1-)-2-pyrrolidinone
75-21-8	Ethylene oxide	91-08-7	Toluene-2,6-diisocyanate (2,6-TDI)
75-25-2	Bromoform	91-20-3	Naphthalene
75-34-3	Dichloroethane(1,1-)	91-57-6	Methylnaphthalene(2-)
		91-58-7	Chloro(2-)naphthalene

92-52-4	Biphenyl	107-06-2	Ethylene dichloride (1,2-Dichloroethane)
95-49-8	Chlorotoluene(o-)	107-07-3	Ethylene chlorohydrin
95-50-1	Dichlorobenzene(o-)	107-13-1	Acrylonitrile
95-53-4	Toluidine(o-)	107-15-3	Ethylenediamine
95-63-6	Trimethylbenzene(1,2,4-)	107-18-6	Allyl alcohol
96-12-8	Dibromochloropropane	107-21-1	Ethylene glycol
96-22-0	Diethyl ketone	107-41-5	Hexylene glycol (2-Methyl- 2,4-pentanediol)
96-33-3	Methyl acrylate	107-87-9	Pentanone(2-) (Methyl propyl ketone)
96-37-7	Methyl cyclopentane	107-98-2	Propylene glycol monomethyl ether
97-63-2	Ethyl methacrylate	108-05-4	Vinyl acetate
98-00-0	Furfuryl alcohol	108-10-1	Methyl isobutyl ketone
98-01-1	Furfural	108-21-4	Isopropyl acetate
98-82-8	Cumene	108-24-7	Acetic anhydride
98-83-9	Methyl styrene(a-)	108-45-2	Phenylene(1,3-) diamine
98-86-2	Acetophenone	108-46-3	Resorcinol
100-41-4	Ethyl benzene	108-65-6	Propylene glycol monomethyl ether acetate
100-42-5	Styrene	108-67-8	Trimethylbenzene(1,3,5-)
100-44-7	Benzyl chloride	108-83-8	Dimethyl(2,6-)-4-heptanone
100-51-6	Benzyl alcohol	108-87-2	Methylcyclohexane
100-52-7	Benzaldehyde	108-88-3	Toluene
100-61-8	Methyl aniline	108-90-7	Chlorobenzene
101-14-4	Methylene(4,4')-bis(2 chloroaniline) (MOCA)	108-91-8	Cyclohexylamine
101-68-8	Methylene bisphenyl isocyanate (MDI)	108-93-0	Cyclohexanol
101-77-9	Methylene(4,4') dianiline (MDA)	108-94-1	Cyclohexanone
102-71-6	Triethanolamine	108-95-2	Phenol
105-46-4	Butyl(sec-) acetate	109-60-4	Propyl(n-) acetate
105-60-2	Caprolactam	109-66-0	Pentane(n-)
106-46-7	Dichlorobenzene(p-)	109-86-4	Methoxyethanol(2-) (Methyl cellosolve)
106-48-9	Chlorophenol(p-)	109-89-7	Diethylamine
106-89-8	Epichlorohydrin	109-99-9	Tetrahydrofuran
106-94-5	Bromopropane(1-)	110-12-3	Methyl isoamyl ketone
106-95-6	Allyl bromide	110-19-0	Isobutyl acetate
106-99-0	Butadiene(1,3-)	110-43-0	Heptanone(2-)
107-02-8	Acrolein	110-49-6	Methoxyethyl(2-) acetate (Methyl cellosolve acetate)
107-04-0	Bromo(1-)-2-chloroethane		
107-05-1	Allyl chloride		

110-54-3	Hexane(n-)	141-32-3	Butyl acrylate
110-62-3	Valeraldehyde	141-43-5	Ethanolamine (2-Aminoethanol)
110-80-5	Ethoxyethanol(2-) (Cellosolve)	141-78-6	Ethyl acetate
110-82-7	Cyclohexane	141-79-7	Mesityl oxide
110-85-0	Piperazine	142-82-5	Heptane
110-86-1	Pyridine	142-96-1	Dibutyl ether
110-91-8	Morpholine	151-67-7	Halothane (Fluothane)
111-15-9	Ethoxyethyl(2-) acetate	191-24-2	Benzo[ghi]perylene
111-30-8	Glutaraldehyde	192-97-2	Benzo[e]pyrene
111-40-0	Diethylene triamine	193-39-5	Indeno[1,2,3-cd]pyrene
111-42-2	Diethanolamine	205-99-2	Benzo[b]fluoranthene
111-76-2	Butoxyethanol(2-) (Butyl cellosolve)	206-44-0	Fluoranthene
111-77-3	Methoxyethoxy(2-(2-)) ethanol	207-08-9	Benzo[k]fluoranthene
111-96-6	Methoxyethyl(2-) ether	208-96-8	Acenaphthylene
112-07-2	Butoxyethyl(2-) acetate	218-01-9	Chrysene
112-24-3	Triethylenetetramine	287-92-3	Cyclopentane
112-34-5	Butoxyethoxy(2-(2-)) ethanol	302-01-2	Hydrazine
117-84-1	Dioctyl phthalate	335-67-1	Perfluorooctanoic acid
120-12-7	Anthracene	431-03-8	Diacetyl
120-82-1	Trichlorobenzene(1,2,4-)	540-59-0	Dichloroethylene(1,2-)
121-44-8	Trimethylamine	540-84-1	Isooctane
123-31-9	Hydroquinone	540-88-5	Butyl(tert-) acetate
123-38-6	Propionaldehyde	552-30-7	Trimellitic anhydride
123-42-2	Diacetone alcohol	563-80-4	Methyl isopropyl ketone
123-54-6	Pentanedione(2,4-)	584-84-9	Toluene-2,4-diisocyanate (2,4-TDI)
123-72-8	n-Butyraldehyde	624-92-0	Dimethyl disulfide
123-86-4	Butyl(n-) acetate	628-63-7	Amyl acetate
123-91-1	Dioxane(p-)	631-64-1	Dibromoacetic acid
124-17-4	Butoxyethoxy(2-(2-)) ethyl acetate	687-47-8	Ethyl lactate
126-73-8	Tributyl phosphate	763-69-9	Ethyl 3-ethoxypropionate
126-98-7	Methylacrylonitrile	822-06-0	Hexamethylene diisocyanate (HDI)
126-99-8	Chloroprene(b-)	872-50-4	Methyl(1-)-2-pyrrolidinone
127-18-4	Perchloroethylene	1304-82-1	Bismuth telluride
127-19-5	Dimethyl acetamide	1305-62-0	Calcium hydroxide
129-00-0	Pyrene	1305-78-8	Calcium oxide
137-26-8	Thiram	1309-37-1	Iron oxide
138-86-3	Limonene(d-)	1309-48-4	Magnesium oxide (fume)
140-88-5	Ethyl acrylate	1310-58-3	Potassium hydroxide



1310-73-2	Sodium hydroxide	7440-06-4	Platinum
1314-13-2	Zinc oxide	7440-16-6	Rhodium as Rh
1314-62-1	Vanadium pentoxide	7440-22-4	Silver
1317-65-3	Calcium carbonate	7440-23-5	Sodium
1319-77-3	Cresol, all isomers	7440-24-6	Strontium
1321-74-0	Divinyl benzene	7440-28-0	Thallium
1327-53-3	Arsenic Trioxide as As	7440-31-5	Tin
1330-20-7	Xylene	7440-32-6	Titanium
1332-21-4	Asbestos	7440-33-7	Tungsten and Compounds as W (in the absence of Cobalt)
1332-58-7	Kaolin	7440-36-0	Antimony
1333-86-4	Carbon black	7440-38-2	Arsenic and inorganic compounds, as As
1344-28-1	Aluminum oxide	7440-39-3	Barium and soluble compounds as Ba
1344-95-2	Calcium silicates synthetic nonfibrous	7440-41-7	Beryllium and compounds as Be
1634-04-4	Methyl tert-butyl ether (MTBE)	7440-42-8	Boron, sodium salts
1675-54-3	Diglycidyl Ether of Bisphenol A	7440-43-9	Cadmium
2426-08-6	Butyl(n-) glycidyl ether	7440-47-3	Chromium
2451-62-9	Triglycidyl isocyanurate	7440-48-4	Cobalt and inorganic compounds as Co
2499-95-8	Hexyl acrylate	7440-50-8	Copper (Fume, Dusts and Mists) as Cu
2807-30-9	Propoxyethanol(2-)	7440-56-4	Germanium
2921-88-2	Chlorpyrifos (Dursban)	7440-57-5	Gold
4098-71-9	Isophorone diisocyanate (IPDI)	7440-62-2	Vanadium
4994-16-5	Phenylcyclohexene(4-)	7440-65-5	Yttrium and compounds, as Y
5124-30-1	Methylene bis(4-cyclohexylisocyanate)	7440-66-6	Zinc
7085-85-0	Ethyl 2-cyanoacrylate	7440-67-7	Zirconium
7429-90-5	Aluminum Metal and insoluble compounds	7440-69-9	Bismuth
7439-91-0	Lanthanum	7440-70-2	Calcium
7439-92-1	Lead and inorganic compounds as Pb	7440-74-6	Indium
7439-93-1	Lithium	7440-09-5	Sulfur dioxide
7439-95-4	Magnesium	7553-56-2	Iodine
7439-96-5	Manganese, elemental and Inorganic compounds as Mn	7646-85-7	Zinc chloride fume
7439-97-6	Mercury as Hg (Elemental and inorganic forms)	7647-01-0	Hydrogen chloride
7439-98-7	Molybdenum	7664-38-2	Phosphoric acid
7440-02-0	Nickel and inorganic compounds as Ni	7664-39-3	Hydrogen fluoride or fluorides as F
7440-05-3	Palladium	7664-41-7	Ammonia
		7664-93-9	Sulfuric acid

7697-37-2	Nitric Acid	10035-10-6	Hydrogen bromide
7722-84-1	Hydrogen peroxide	10049-04-4	Chlorine dioxide
7723-14-0	Phosphorus	10102-43-9	Nitric oxide
7726-95-6	Bromine	10102-44-0	Nitrogen dioxide
7758-97-6	Lead chromate as Cr	11097-69-1	Chlorodiphenyl (54% chlorine)
7778-18-9	Calcium sulfate	12001-26-2	Mica
7782-42-5	Graphite	12125-02-9	Ammonium chloride
7782-49-2	Selenium and compounds as Se	13463-67-7	Titanium dioxide
7782-50-5	Chlorine	13494-80-9	Tellurium
7782-65-2	Germanium tetrahydride	13765-19-0	Calcium chromate as Cr
7783-06-4	Hydrogen sulfide	13838-16-9	Enflurane (Ethrane)
7784-42-1	Arsine	14464-46-1	Silica cristobalite
7789-06-2	Strontium chromate as Cr	14807-96-6	Talc
7803-51-2	Phosphine	14808-60-7	Silica quartz
8002-74-2	Paraffin wax fume	22204-53-1	Naproxen sodium
8003-34-7	Pyrethrum	26675-46-7	Forane (Isoflurane)
8006-61-9	Gasoline	28182-81-2	Hexamethylene diisocyanate(1,6-) Homopolymer (HDI Polymer)
8008-20-6	Kerosene	28523-86-6	Sevoflurane (Sevofrane)
8012-95-1	Mineral Oil	34590-94-8	Dipropylene glycol methyl ether
8032-32-4	VM&P Naphtha	53469-21-9	Chlorodiphenyl (42% chlorine)
8050-09-7	Resin Acids	57041-67-5	Desflurane (Suprene)
8052-41-3	Stoddard solvent	64742-95-6	Aromatic 100
8052-42-4	Asphalt fume	65996-93-2	Coal tar pitch volatiles
9002-86-2	Polyvinyl chloride (PVC)	65997-15-1	Portland cement
9005-25-8	Starch	88917-22-0	Dipropylene glycol methyl ether acetate (DPGMEA)
10024-97-2	Nitrous oxide		
10028-15-6	Ozone		



The illustrations, instructions, and principles contained in the material are general in scope and, to the best of our knowledge, current at the time of publication. Our risk control services are advisory only. We assume no responsibility for: managing or controlling customer safety activities, implementing any recommended corrective measures, or identifying all potential hazards.

No attempt has been made to interpret any referenced codes, standards, or regulations. Please refer to the appropriate government authority for interpretation or clarification.

Insurance underwritten by Liberty Mutual Insurance Co. or its affiliates or subsidiaries.

© 2020 Liberty Mutual Insurance, 175 Berkeley Street, Boston, MA 02116.

[libertymutualgroup.com/riskcontrolservices](https://libertymutualgroup.com/riskcontrolservices)

2020 Edition Lab Sampling Guide